

Product datasheet for **RG231795**

UQCRB (NM_001254752) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: UQCRB (NM_001254752) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: UQCRB
Synonyms: MC3DN3; QCR7; QP-C; QPC; UQBC; UQBP; UQCR6; UQPC
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG231795 representing NM_001254752
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGGTAAGCAGGCCGTTTCAGCATCAGGCAAGTGGCTGGATGGTATTCGAAAATGGTATTACAATG
CTGCAGGATCAATAAACTGGGGTTAATGCGAGATGATAAATACGAGGATGAAGATGTAAAAGAAGC
CATAAGAAGACTTCCTGAGAACCTTTAATGACAGGATGTTTCGCATTAAGAGGGCACTGGACCTGAAC
TTGAAGCATCAGATCTTGCTAAAGAGCAGTGGACCAAATATGAAGAGGTCTTTGCTGTTCCAGCTCTGC
ACTCTGCTTCTACTTAGATGAAAAGATCAGCCATTGAGTGTCCCTCCAGATCCCAAGAAGAGCTTCTG
TGAAAGCAAACCTCACCCCTTGAAGTGCATTGAAACTAGGCAAAGGAAAATTTCTACCTGAACCGTATC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG231795 representing NM_001254752
Red=Cloning site **Green**=Tags(s)
MAGKQAVSASGKWLGDGIRKYYNAAGFNKLGLMRDDTIYEDEDVKEAIRRLPENLYNDRMFRIKRALDLN
LKHQILPKEQWTKYEEVFAVPALHSASYLDEKISPLSVPPDPKKSFCANSHPLNCIETRQRKISTLNRI

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001254752

ORF Size: 420 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001254752.2
RefSeq Size:	4976 bp
RefSeq ORF:	423 bp
Locus ID:	7381
UniProt ID:	P14927
Cytogenetics:	8q22.1
Protein Pathways:	Alzheimer's disease, Cardiac muscle contraction, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease
Gene Summary:	This gene encodes a subunit of the ubiquinol-cytochrome c oxidoreductase complex, which consists of one mitochondrial-encoded and 10 nuclear-encoded subunits. The protein encoded by this gene binds ubiquinone and participates in the transfer of electrons when ubiquinone is bound. This protein plays an important role in hypoxia-induced angiogenesis through mitochondrial reactive oxygen species-mediated signaling. Mutations in this gene are associated with mitochondrial complex III deficiency. Alternatively spliced transcript variants have been found for this gene. Related pseudogenes have been identified on chromosomes 1, 5 and X. [provided by RefSeq, Dec 2011]