

Product datasheet for **RG231957**

ZNF695 (NM_001204221) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF695 (NM_001204221) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZNF695
Synonyms:	SBZF3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG231957 representing NM_001204221 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGACTATTGGCATTAGGGATGTGGCTCTAGAATTCTCTCCAGAGGAGTGGGAATGCCTGGACCCAG
CTCAGCGGAGTTTGTATAGGGATGTGATGTTAGAGAACTACAGAAACCTGATCTCCCTTGGTGAGGATAG
CTTCAATATGCAATTCCTATTTACAGTCTTGCTATGTCTAAGCCAGAACTGATCATCTGTCTGGAGGCA
AGGAAAGAGCCCTGGAACGTGAACACAGAGAAGACAGCCAAACACTCAGTTTGTCTTCTTATCTTACTG
AAGACATTTTGCAGAGCAGGGCTGCAAGTTTCATTCCAAAAAGTGATGCTGAGAAGATATGAAAGATG
TTGTCTTGAGAAATTACGCTTAAGGAATGACTGGGAAATTCATGTGAAGATGTGCTTGCTTCCCTTTG
CCTTCTGCCATGATTCTAAGTTTCCTGAGGCCTCCCCAGAAGCAGAAGCATGTAAAGCCACAGAACCG
TCCAGTCGTCGACCATCGAGCTGCTT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:	>RG231957 representing NM_001204221 Red=Cloning site Green=Tags(s)
-------------------	---


MGLLAFRDVALEFSPEWECLDPAQRSLYRDVMLENYRNLI SLGEDSFNMQFLFHSLAMSKPELIICLEA
RKEPWNVNTAKHSLSSYLTEDILPEQGLQVSFQKVMRLRRYERCCLEKLRLRNDWEIPCEDVLASPL
PSAMILSFLRPPQKQKHVKPTEPVQSSTIELL

TRTRPLE - GFP Tag - V

Restriction Sites:	SgfI-MluI
--------------------	-----------



[View online »](#)

Cloning Scheme:


ACCN: NM_001204221

ORF Size: 516 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204221.1](#), [NP_001191150.1](#)

RefSeq Size: 954 bp


RefSeq ORF: 519 bp

Locus ID: 57116

UniProt ID: [Q8IW36](#)

Cytogenetics: 1q44
Protein Families: Transcription Factors
Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG231957