

## Product datasheet for **RG232057**

### CLECL1 (NM\_001253750) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** CLECL1 (NM\_001253750) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** CLECL1  
**Synonyms:** DCAL-1; DCAL1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG232057 representing NM\_001253750  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGTTAGTAATTTCTTCCATGTCATACAAGTATTCGAGAAATCTGCTACCTTGATTAGTAAGACTGAAC  
ACATTGGTTTTGTCATTTATTTCATGGAGGAAGTCCACCACCCACTTGGGGAGCAGAAGGAAATTTGCCAT  
CTCAATTTACTTATCAGAAGTTTCTTTCAGAAATATGATTGTCCTTCAGTGGGACATCATTTGTGGTC  
TTCTCTCTTTTTGATCTGTGCAATGGCTGGAGATGTAGTCTACGCTGACATCAAACTGTTCCGACTT  
CCCCGTTAGAAGCTCGCGTTTCCACTTCAGAGATCTGTTTCTTTCACTTTTCTACTGTCCATAAATCATG  
TCCTGCCAAAGACTGGAAGGTGCATAAGGGAAAATGTTACTGGATTGCTGAACTAAGAAATCTTGGAAC  
AAAAGTCAAAATGACTGTGCCATAAACAATTCATATCTCATGGTGATTCAAGACATTACTGCTATGCTGA  
TCACCTTTTTGCTGTACCTTCAAATGGCATTCTTCAATGTGTGTACATCTGTGATGCTCTGTGTATCC  
AAATTTCTTCTTG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG232057 representing NM\_001253750  
Red=Cloning site Green=Tags(s)


MVSNFFHVIQVFEKSATLISKTEHIGFVIYSWRKSTTHLGSRRKFAISIYLSEVSLQKYDCPFSGTSFV  
FSLFLICAMAGDVVYADIKTVRSTPLELAFPLQRSVSNFSTVHKSCPAKDWKVKHKCYWIAETKKS  
KSQNDCAINNSYLMVIQDITAMLITFLLYLQMAFLQCVYISDVSVPNFLL

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



**Cloning Scheme:**

**Plasmid Map:**


**ACCN:** NM\_001253750

**ORF Size:** 573 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001253750.1](#), [NP\\_001240679.1](#)

**RefSeq Size:** 1586 bp

**RefSeq ORF:** 576 bp

**Locus ID:** 160365

**UniProt ID:** [Q8IZS7](#)

**Cytogenetics:** 12p13.31

**Protein Families:** Druggable Genome

**Gene Summary:** This gene encodes a type II transmembrane, C-type lectin-like protein that is highly expressed on dendritic and B cells. This protein may act as a T-cell costimulatory molecule that enhances interleukin-4 production, and maybe involved in the regulation of the immune response. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011]