

## Product datasheet for **RG232081**

### **ZA20D3 (ZFAND6) (NM\_001242917) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ZA20D3 (ZFAND6) (NM\_001242917) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** ZFAND6  
**Synonyms:** AWP1; ZA20D3; ZFAND5B  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG232081 representing NM\_001242917  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCTCAAGAACTAATCACAGCCAAGTGCCTATGCTTTGTTCCACTGGCTGTGGATTTTATGAAACC  
CTCGTACAAATGGCATGTGTTTCAGTATGCTATAAAGAACATCTTCAAAGACAGAATAGTAGTAATGGTAG  
AATAAGCCACCTGTTCAATGCACAGATGGCAGTGTGCCAGAAGCCAGTCAGCATTAGACTCTACATCT  
TCATCTATGCAGCCAGCCCTGTATCAAATCAGTCACTTTATCAGAATCTGTAGCATCTTCTCAATTGG  
ACAGTACATCTGTGGACAAAGCAGTACCTGAAACAGAAGATGTGCAGGCTTCAGTATCAGACACAGCACA  
GCAGCCATCTGAAGAGCAAAGCAAGTCTCTTGA AAAACCGAAACAAAAAAGAATCGCTGTTTCATGTGC  
AGGAAGAAAGTGGGACTTACTGGGTTTGAATGCCGGTGTGAAATGTTTACTGTGGTGTACACCGTACT  
CAGATGTACAAATTGCTCTTACAATTACAAAGCCGATGCTGCTGAGAAAATCAGAAAAGAAAATCCAGT  
AGTTGTTGGTAAAAGATCCAAAAGATT

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >RG232081 representing NM\_001242917  
**Red=Cloning site Green=Tags(s)**



MAQETNHSQVPMCLCSTGCGFYGNPRTNGMCSVCYKEHLQRQNSSNGRISPPVQCTDGSVPEAQSALDSTS  
SSMQPSPVSNQSLLESVASSQLDSTSVDKAVPETEDVQASVSDTAQQPSEEQSKSLEKPKQKKNRCFMC  
RKKVGLTGFECCGNVYCGVHRYSDVHNCYNYKADAAEKIRKENPVVVGEKIQKI

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**

**Plasmid Map:**

**ACCN:** NM\_001242917

**ORF Size:** 588 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)
**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001242917.1](#), [NP\\_001229846.1](#)

**RefSeq Size:** 1672 bp

**RefSeq ORF:** 591 bp

**Locus ID:** 54469

**UniProt ID:** [Q6FIF0](#)

**Cytogenetics:** 15q25.1

**Gene Summary:** Involved in regulation of TNF-alpha induced NF-kappa-B activation and apoptosis. Involved in modulation of 'Lys-48'-linked polyubiquitination status of TRAF2 and decreases association of TRAF2 with RIPK1. Required for PTS1 target sequence-dependent protein import into peroxisomes and PEX5 stability; may cooperate with PEX6. In vitro involved in PEX5 export from the cytosol to peroxisomes (By similarity).[UniProtKB/Swiss-Prot Function]