

Product datasheet for **RG232112**

RPL15 (NM_001253380) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RPL15 (NM_001253380) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: RPL15
Synonyms: DBA12; EC45; L15; RPL10; RPLY10; RPYL10
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG232112 representing NM_001253380
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGTGCATACAAGTACATCCAGGAGCTATGGAGAAAGAAGCAGTCTGATGTCATGCGCTTTCTTCTGA
 GGGTCCGCTGCTGGCAGTACCGCCAGCTCTGCTCTCCACAGGGCTCCCCGCCCCACCCGGCCTGATAA
 AGCGCGCCGACTGGGCTACAAGGCCAAGCAAGGTTACGTTATATATAGGATTCGTGTTCCCGTGGTGGC
 CGAAAACGCCAGTTCCTAAGGGTGCAACTACGGCAAGCCTGTCCATCATGGTGTAAACCAGCTAAAGT
 TTGCTCGAAGCCTTCAGTCCGTTGCAGAGGAGCGAGCTGGACGCCACTGTGGGGCTCTGAGAGTCTGAA
 TTCTTACTGGGTTGGTGAAGATTCCACATACAAATTTTTGAGGTTATCCTCATTGATCCATTCCATAAA
 GCTATCAGAAGAAATCCTGACACCCAGTGGATCACAAACCAGTCCACAAGCACAGGGAGATGCGTGGGC
 TGACATCTGCAGGCCGAAAGAGCCGTGGCCTTGAAAGGGCCACAAGTCCACCACACTATTGGTGGCTC
 TCGCCGGGCGACTTGAGAAAGGCCAATACTCTCCAGCTCCACCGTTACCGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG232112 representing NM_001253380
 Red=Cloning site Green=Tags(s)

MGAYKYIQELWRKKQSDVMRFLLRVRCWQYRQLSALHRAPRTRPDKARRLGYKAKQGYVIYRIRVRRGG
 RKRPVPGATYGKPVHHGVNQLKFARSLQSVAEERAGRHCALRVLNSYWVGEDSTYKFFEVIIDPFHK
 AIRRNPDTQWITKPVHKHREMRGLTSAGRKSRGLGKHKFHHTIGGSRAAWRRRNTLQLHRYR

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001253380.2](#)

RefSeq Size: 2125 bp

RefSeq ORF: 615 bp

Locus ID: 6138

UniProt ID: [P61313](#)

Cytogenetics: 3p24.2

Protein Pathways: Ribosome

Gene Summary: Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of four RNA species and approximately 80 structurally distinct proteins. This gene encodes a member of the L15E family of ribosomal proteins and a component of the 60S subunit. This gene shares sequence similarity with the yeast ribosomal protein YL10 gene. Elevated expression of this gene has been observed in esophageal tumors and gastric cancer tissues, and deletion of this gene has been observed in a Diamond-Blackfan anemia (DBA) patient. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Mar 2017]