

Product datasheet for **RG232148**

FAM156A (NM_001242491) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM156A (NM_001242491) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	FAM156A
Synonyms:	PRO0659; TMEM29
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG232148 representing NM_001242491 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGATCCACTCCAGAAACGGAATCCAGCATCGCCTTCCAAATCTTCCCGATGACAGCTGCAGAGACTT
CCCAGGAAGGTCCAGCGCCCTCTCAGCCTTCGTAAGCAGGCGATGATGGCCTCAGTAACCTGAG
CCCCGGTCTGGCCCCAGCCAGGCGTGCCTCTCCAGAGGGGCTGCTCCGCCAGCGGTACAGAGAGGAG
AAGACCCTGGAAGAGCGGCGGTGGGAGAGGCTGGAGTTCTTCAGAGGAAGAAAGCATTCTCGGCATG
TGAGGAGGAGACACCGCATCACATGGCCCCATGCTGTGGGAGGGAAGCCAGAATCTCCCCATTAGG
TGACAGAAGTCAGAAATCGATTCCGATGTGAATGTGATACTGCCAGAGCCACAGGCCGAATCTTTCTGGG
ATCCCTGGGGAGAGTAACAGGGCCCCACATCCCTCCTCCTGGGAGACGCTGGTGCAGGGCCTCAGTGGCT
TGACTCTCAGCCTAGGCACCAACAGCCCGGCCTCTGCCTGAAGCGGCACTCCAGCCACAGGAGACAGA
GGAGAAGCGCCAGCGAGAGAGGCAGCAGGAGAGCAAAATAATGTTTCAGAGGCTGCTCAAGCAGTGGTTA
GAGGAAAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >RG232148 representing NM_001242491
Red=Cloning site Green=Tags(s)

MDPLQKRNPASPSKSSPMTAAETSQEGPAPSQPSYSEQPMGLSNLSPGPGPSQAVPLPEGLLRQRYREE
 KTLLEERRWERLEFLQRKKAFLRHVRRRRHRDHMAPYAVGREARISPLGDRSQNRFRCCECRYCQSHRPNL
 SGIPGESNRAPHSPSSWETLVQGLSGLTSLGTNPQGPLPEAALQPQETEEKRQERQEQESKIMFQRLKQWL
 EEN


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001242491

ORF Size: 639 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001242491.1 , NP_001229420.1
RefSeq Size:	2060 bp
RefSeq ORF:	642 bp
Locus ID:	29057
UniProt ID:	Q8NDB6
Cytogenetics:	Xp11.22