

# **Product datasheet for RG232250**

## CD63 (NM\_001257389) Human Tagged ORF Clone

### **Product data:**

**Product Type:** Expression Plasmids

**Product Name:** CD63 (NM\_001257389) Human Tagged ORF Clone

Tag: TurboGFP

Symbol: CD63

Synonyms: LAMP-3; ME491; MLA1; OMA81H; TSPAN30

Mammalian Cell

Selection:

Neomycin

Vector: pCMV6-AC-GFP (PS100010)

E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide >RG232250 representing NM\_001257389
Sequence: Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC

GCCGCGATCGCC

ATGGCGGTGGAAGGAATGAAATGTGTGAAGTTCTTGCTCTACGTCCTCGCTGGCCTTTTGCGCCT
GTGCAGTGGAACGACCATGATTGCCGTGGGTGTCGGGGCACAGCTTGTCCTGAGTCAGACCATAATCCAGGGGGC
TACCCCTGGCTCTCTGTTGCCAGTGGTCATCATCGCAGTGGGTGTCTTCCTCTTCCTGGTGGCTTTTGTG
GGCTGCTGCGGGGCCTGCAAGGAGAACTATTGTCTTATGATCACGTTTGCCATCTTTCTGTCTCTTATCA
TGTTGGTGGAGGTGGCCGCAGCCATTGCTGGCTATGTGTTTAGAGATAAGGTGATGTCAGAGTTTAATAA
CAACTTCCGGCAGCAGATGGAGAAATTACCCGAAAAACAACCACACTGCTTCGATCCTGGACAGGATGCAG
GCAGATTTTAAGTGCTGTGGGGCTGCTAACTACACAGATTGGGAGAAAATCCCTTCCATGTCGAAGAACC
GAGTCCCCGACTCCTGCTGCATTAATGTTACTGTGGGCTTGTGGGATTAATTTCAACGAGAAGGCGATCCA
TAAGGAGGGCTGTTGGGAGAAGATTTGGGGGCTGGCTGAGGAAAAATGTGCTGGTGGTAGCTGCAGCAGCC
CTTGGAATTGCTTTTTGTCGAGGTTTTTGGGAATTGTCTTTTGCCTGCTCCTCGTGAAGAAGTTCAGAAGTG
GCTACGAGGTGATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Protein Sequence: >RG232250 representing NM\_001257389

Red=Cloning site Green=Tags(s)


MAVEGGMKCVKFLLYVLLLAFCACAVGLIAVGVGAQLVLSQTIIQGATPGSLLPVVIIAVGVFLFLVAFV GCCGACKENYCLMITFAIFLSLIMLVEVAAAIAGYVFRDKVMSEFNNNFRQQMENYPKNNHTASILDRMQ ADFKCCGAANYTDWEKIPSMSKNRVPDSCCINVTVGCGINFNEKAIHKEGCVEKIGGWLRKNVLVVAAAA LGIAFVEVLGIVFACCLVKSIRSGYEVM

**Restriction Sites:** 


Sgfl-Mlul

**Cloning Scheme:** 





### Plasmid Map:



**ACCN:** NM\_001257389

ORF Size: 714 bp

#### CD63 (NM\_001257389) Human Tagged ORF Clone - RG232250

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001257389.1</u>, <u>NP 001244318.1</u>

RefSeq Size:1238 bpRefSeq ORF:717 bpLocus ID:967

 UniProt ID:
 P08962

 Cytogenetics:
 12q13.2

**Protein Families:** Druggable Genome, Transmembrane

**Protein Pathways:** Lysosome

**Gene Summary:** The protein encoded by this gene is a member of the transmembrane 4 superfamily, also

known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. The encoded protein is a cell surface glycoprotein that is known to complex with integrins. It may function as a blood platelet activation marker. Deficiency of this protein is associated with Hermansky-Pudlak syndrome. Also this gene has been associated with tumor progression. Alternative splicing results in multiple transcript variants encoding different

protein isoforms. [provided by RefSeq, Apr 2012]