

## Product datasheet for **RG232262**

### MFF (NM\_001277063) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	MFF (NM_001277063) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MFF
Synonyms:	C2orf33; EMPF2; GL004
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG232262 representing NM_001277063 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCAGAAATTAGTCGAATTCAGTACGAAATGGAATATACTGAAGGCATTAGTCAGCGAATGAGGGTCC  
CAGAAAAGTTAAAAGTAGCACCGCCAACGCTGACCTGGAACAAGGATTCCAAGAAGGAGTTCCAAATGC  
TAGTGTGATAATGCAAGTCCGGAGAGGATTGTTGTAGCAGGAAATAATGAAGATGTTTCATTTCAAGA  
CCAGCAGATCTTGACCTTATTCAGTCAACTCCCTTTAAACCCCTGGCACTGAAAACACCACCTCGTGAC  
TTACGCTGAGTGAAGACCACTAGATTTTCTGGATTTAGAAAGACCTCCTACAACCCCTCAAATGAAGA  
AATCCGAGCAGTTGGCAGACTAAAAGAGAGCGGTCTATGAGTGAATAATGCTGTTTCGCAAAATGGACAG  
CTGGTCAGAAATGATTCTCTGTGGCACAGATCAGATTCTGCCCAAGAAATAAAATTTCAAGGTTCCAGG  
CACCGATTTCTGCACCGGAGTACACGTATGGCATTCAAATATAGATAACAACATTGAAGGAACGTCAGA  
TGACCTGACTGTTGTAGATGCAGCTTCACTAAGACGACAGATAATCAAATAAATAGACGCTACAACCT  
CTGGAAGAGGAGAAACAAGAACGTGCTAAAAGAGAAATGGTCATGTATTCAATTACTGTAGCTTTCTGGC  
TGCTTAATAGCTGGCTCTGTTTCGCCGC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG232262 representing NM\_001277063  
Red=Cloning site Green=Tags(s)

MAEISRIQYEMEYTEGISQRMRVPEKLVAPPNADLEQGFQEGVFNASVIMQVPERIVVAGNNEDVSFSR  
 PADLDLIQSTPFKPLALKTPPRVLTLSERPLDFDLERPPTTPQNEEIRAVGRLKRERSMSENAVRQNGQ  
 LVRNSLWHRSDSAPRNKISRFAQPISAPEYTYGINSIDTTIEGTSDDLTVVDAASLRRIKLNRRLLQL  
 LEEENKERAKREVMYSITVAFWLLNSLWFRR

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** NM\_001277063

**ORF Size:** 729 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001277063.2</a>
<b>RefSeq Size:</b>	1880 bp
<b>RefSeq ORF:</b>	732 bp
<b>Locus ID:</b>	56947
<b>UniProt ID:</b>	<a href="#">Q9GZY8</a>
<b>Cytogenetics:</b>	2q36.3
<b>Protein Families:</b>	Transmembrane
<b>Gene Summary:</b>	This is a nuclear gene encoding a protein that functions in mitochondrial and peroxisomal fission. The encoded protein recruits dynamin-1-like protein (DNM1L) to mitochondria. There are multiple pseudogenes for this gene on chromosomes 1, 5, and X. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2013]