

## Product datasheet for **RG232276**

### **C4ORF28 (PACRGL) (NM\_001258345) Human Tagged ORF Clone**

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Expression Plasmids  |
| Product Name:             | C4ORF28 (PACRGL) (NM_001258345) Human Tagged ORF Clone                         |
| Tag:                      | TurboGFP   |
| Symbol:                   | PACRGL   |
| Synonyms:                 | C4orf28  |
| Mammalian Cell Selection: | Neomycin   |
| Vector:                   | pCMV6-AC-GFP (PS100010)  |
| E. coli Selection:        | Ampicillin (100 ug/mL)   |
| ORF Nucleotide Sequence:  | >RG232276 representing NM_001258345<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCAGAAATCAGAGGGCTCTGGAGGTACACAGTTGAAAAACAGAGCAACAGGTAACATGATCAAAGGA  
CATCATCAAGCACACAGTTAAAACACAGGAATGCAGTTCAGGGAAGCAAATCCTCATTGTCAACCAGTTC  
TCCAGAGTCTGCAAGAAAACCTTCATCCTAGACCAAGTGATAAACTGAACCCTAAAACAATTAATCCGTTT  
GGTGAACAGTACGAGTGCCTTCTGCATTTGCAGCTATTTACTCTAAAGGAGGTATTCCTTGAGATTGG  
TACATGGTTCAGTAAAACACAGATTACAGTGGGAATGTCTCCTGAAAGTCTTTCATTTGATCCACTTCT  
TATTACTTTAGCTGAGGGTCTGAGAGAGACTAAGCATCCATACACTTTTGTGTCAAAGGAGGGTTTTAGA  
GAATTACTTTTGGTCAAAGGTGCTCCTGAAAAAGCTATTCCTTTGCTACCTAGACTGATTCTGTGCTAA  
AGGCAGCTCTGGTCCATTCGGATGATGAAGTGTGTTGAAAGAGGATTGAATGCTCTAGTTCAGCTAAGTGT  
CGTTGTTGGTCTTCTCTAAACGACCATCTGAAGCATCTGCTTACAAGCCTTCCAAAGAGATTAATGGAC  
AAGAAATTCAAAGAGCCAATCACCAGCGCATTACAAAAGCTAGAGCAACATGGTGGAAAGTGGGAGCCTTA  
GCATCATCAAATCTAAAATTC AACATACTGCTCCATATGCTGT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG232276 representing NM\_001258345  
Red=Cloning site Green=Tags(s)


MQKSEGGTQLKNRATGNVDQRTSSSTQLKHRNAVQGSKSSLSTSSPESARKLHPRPSDKLNPKTINPF  
 GEQSRVPSAFAAIYSKGGIPCRLVHGSVKHRLQWECPPESLSFDPLLITLAEGLRETKHPYTFVSKEGFR  
 ELLLKGAPEKAIPLLPRLIPVLKAALVHSDDEVFERGLNALVQLSVVVGPSLNDHLKHLTSLSKRLMD  
 KKFKEPITSALQKLEQHGGSGLSIIKSKIPTYCSICC

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** NM\_001258345

**ORF Size:** 744 bp

|                               |   |
|-------------------------------|---|
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>  |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <a href="#">NM_001258345.3</a>  |
| <b>RefSeq Size:</b>           | 2171 bp   |
| <b>RefSeq ORF:</b>            | 747 bp  |
| <b>Locus ID:</b>              | 133015  |
| <b>UniProt ID:</b>            | <a href="#">Q8N7B6</a>  |
| <b>Cytogenetics:</b>          | 4p15.31   |