

Product datasheet for **RG232412**

FSBP (NM_001256141) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FSBP (NM_001256141) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	FSBP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG232412 representing NM_001256141 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTAGGAAAGGCTAGATCTTCCAATTTACCTTATCCGAAAAGCTTGATTTGCTAAAGCTTGTGAAGC
CATATGTGAAAATTCTCGAAGAACACACTAATAAACATTCACTAATAGTAGAAAAGAATAGATGTTGGGA
TATCATAGCAGTTAACTATAATGCAATTGGAGTAGACCGCCCTCCTCGAACAGCACAGGGCTACGCACC
CTTTATAAAAGGCTCAAAGAATATGCCAAACAGGAGCTATTGCAGCAAAAAGAGACCCAATCAGATTTTA
AAAGCAATATTTCTGAGCCAACCAAGAAAGTTATGGAGATGATTTCCAGATTTCCAGTTTTGCTGGT
AAGAGACAGGAACACATACAAAGTGCAAACTTGATGAGGAGGCACAGGCTGGTACCAAGTTCACCTACAG
GTAATGTTGGATCACCATCCTGTTGCTATTACAGTGGAGGTGAAGCAAGAAGAAGACATTAACACCTC
CTCCACTGGTTTTAAATTCTCAACAGAGTGATACTTTAGAGCAAAGAGAAGAATGAATTAGTACATGT
TATGGAAGATCCTTGTGCGCTTCACTTTCCTCTGTTGATATGAGAATGACATCGTCTCCATCTTCTATT
CCAAGGAGAGATGATTTTTTTCGGCATGAGAGTGGTGAACACTTTAGGTCACTATTAGGGTATGATCCTC
AGATCCTGCAAAATGTTGAAAGAGGAGCATCAGATAATTTAGAAAATCAAAAAATTTGGATTGTATGT
TCAGGAGAAGAGGATGGATTGAAAAGAAGGCAGCAGCTAGAGGAAGAGCTACTAAGAGCAAAAATTGAA
GTGGAGAAGCTGAAAGCAATTCGCTTACGGCATGATCTACCTGAATATAACAGTCTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG232412 representing NM_001256141
Red=Cloning site Green=Tags(s)



MVGKARSSNFTLSEKDLLKLKPYVKILEEHTNKHSHVIVEKNRCWDIIAVNYNAIGVDRPPRTAQGLRT
LYKRLKEYAKQELLQKETQSDFKSNISEPTKKVMEMIPQISSFCLVRDRNHIQSANLDEEAQAGTSSLQ
VMLDHHHPVAITVEVKQEEDIKPPPLVLNSQQSDTLEQREEHEL VHVMERSLSPSLSSVDMRMTSSPSSI
PRRDDFFRHESGEHFRSLLGYDPQILQMLKEEHQIILENQKNFGLYVQEKRDGLKRRQQLLEEELLRAKIE
VEKLKAIRLRHDLPEYNSL

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



ACCN: NM_001256141

ORF Size: 897 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001256141.2](#)

RefSeq Size: 5418 bp


RefSeq ORF: 900 bp

Locus ID: 100861412

UniProt ID: [O95073](#)

Cytogenetics: 8q22.1

Product images:



Circular map for RG232412