

## Product datasheet for **RG232450**

### STOML1 (NM\_001256676) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	STOML1 (NM_001256676) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	STOML1
Synonyms:	hUNC-24; SLP-1; STORP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG232450 representing NM_001256676 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGATTGTGTTCCGCTGGGCCGGATCCGCACCCCCAGGGACCTGGCATGGTTCTGCTCTTGCCCTTCA  
TTGACTCCTTTCAGAGGGTGGATCTGAGGACACGAGCCTTCAACGTCCTCCCTGCAAGCTGGCCTCTAA  
GGACGGGGCTGTGCTGTCCGTGGGAGCCGATGTCCAGTTTCGCATCTGGGACCCGGTGTGTCGGTGATG  
ACTGTGAAAGACCTGAACACAGCCACACGCATGACAGCCAGAAGCCATGACCAAGGCCCTGCTCAAGA  
GGCCGCTGCGGGAGATCCAGATGGAGAAGCTCAAGATCAGCGACCAGCTTCTGCTGGAGATCAACGATGT  
GACCAGGGCCTGGGGCTGGAGGTAGACCGCTGGAGCTGGCAGTGGAGGCCGTGCTCCAGCCGCCCCAG  
GACAGCCAGCTGGGCCAACCTGGACAGCACCTCCAGCAGCTGGCCCTGCACTTCTGGGAGGAAGCA  
TGAACCAATGGCAGGAGGTGCCCGTCCCCGGGGCCAGACACCGTGGAGATGGTGAAGTTGAGCC  
ACCTGCCCTCAAGTTGGTGCCAGGTCCAGTCCGAAGCAGCCTCTGGCGGAGGGGCTACTGACTGCTCTA  
CAGCCCTTCTGTCTGAGGCCCTGGTCAAGCAAGTCGGGGCTGTACCAGTTCAATGTCGTCCTGCCCA  
GCGCACCCAAAGCGCCTACTTCTGGACCTCACTACAGGACGAGGAAGAGTGGGACACGGGGTGCCTGA  
TGGCATCCCTGATGTGGTGGTGGAGATGGCCGAGGCAGACCTGCGGGCCCTGCTATGCAGAGAGCTGCGG  
CCCCTGGGGCCTACATGAGTGGACGGCTGAAGGTGAAGGGCGACCTGGCTATGGCCATGAAGCTGGAGG  
CTGTCTCAGGGCCTTGAAG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG232450 representing NM\_001256676  
 Red=Cloning site Green=Tags(s)

MIVFRLGRIRTPQGGMVLLLPFIDSFQRVDLRTRAFNVPPCKLASKDGAVL SVGADVQFRIWDPVLSVM  
 TVKDLNATRMATAQNAMTKALLKRPLREIQMEKLIKISDQLLEINDVTRAWGLEVDRELAVEAVLQPPQ  
 DSPAGPNLDSTLQQLALHFLGGSMNSMAGGAPSPGPDVTVMVSEVEPPAPQVGARSSPKQPLAEGLLTAL  
 QPFLSEALVSQVGACYQFNVLPSGTQSAYFLDLTTGRGRVGHGVPDGI PDVVVEMAEADLRALLCRELR  
 PLGAYMSGRLKVKGD LAMAMKLEAVLRALK

TRTRPLE - GFP Tag - V

**Restriction Sites:**


SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**Plasmid Map:**



**ACCN:** NM\_001256676

**ORF Size:** 930 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001256676.1</a> , <a href="#">NP_001243605.1</a>
<b>RefSeq Size:</b>	1852 bp
<b>RefSeq ORF:</b>	933 bp
<b>Locus ID:</b>	9399
<b>UniProt ID:</b>	<a href="#">Q9UBI4</a>
<b>Cytogenetics:</b>	15q24.1
<b>Protein Families:</b>	Transmembrane
<b>Gene Summary:</b>	May play a role in cholesterol transfer to late endosomes (PubMed:19696025). May play a role in modulating membrane acid-sensing ion channels. Can specifically inhibit proton-gated current of ASIC1 isoform 1. Can increase inactivation speed of ASIC3. May be involved in regulation of proton sensing in dorsal root ganglions (By similarity). May play a role in protecting FBXW7 isoform 3 from degradation (PubMed:23082202).[UniProtKB/Swiss-Prot Function]