

## Product datasheet for **RG232566**

### PRKACB (NM\_001242859) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PRKACB (NM_001242859) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PRKACB
Synonyms:	CAFD2; PKA C-beta; PKACB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG232566 representing NM_001242859 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGATTGTCACGCAAATCATCAGATGCATCTGCTTCTCTCAGAAATATCTGTGAAAGAGTTTC  
TAGCCAAAGCCAAAGAAGACTTTTTGAAAAATGGGAGAATCCAACCTCAGAATAATGCCGGACTGAAGA  
TTTTGAAAGGAAAAAACCCCTTGAACAGGTTTCATTTGGAAGAGTCATGTTGGTAAACACAAAGCCACT  
GAACAGTATTATGCCATGAAGATCTTAGATAAGCAGAAGTTGTTAACTGAAGCAAATAGAGCATACTT  
TGAAATGAGAAAAGAATATTACAGGCAGTGAATTTCTTTCTTGTTCGACTGGAGTATGCTTTTAAGGA  
TAATTCTAATTTATACATGGTTATGGAATATGTCCCTGGGGTGAAATGTTTTACATCTAAGAAGAATT  
GGAAGGTTCAAGTACCCATGCACGGTCTATGCAGCTCAGATAGTGCTAACATTTCGAGTACCTCCATT  
CACTAGACCTCATCTACAGAGATCTAAAACCTGAAAATCTCTTAATTGACCATCAAGGCTATATCCAGGT  
CACAGACTTTGGGTTTGCCAAAAGAGTTAAAGGCAGAACTTGGACATTATGTGGAACCTCAGAGTATTTG  
GCTCCAGAAAATTTCTCAGCAAGGGCTACAATAAGGCAGTGGATTGGTGGGCATTAGGAGTGCTAATCT  
ATGAAATGGCAGCTGGCTATCCCCATTCTTGCAGACCAACCAATTGAGATTTATGAAAAGATTGTTTC  
TGGAAAGTCCGATCCCATCCCACTTCAGTTCAGATCTCAAGGACCTTCTACGGAACCTGCTGCAGGTG  
GATTTGACCAAGAGATTTGAAAATCTAAGAATGGTGTGATATAAAAACTCACAAGTGTTTGCCA  
CGACAGATTGGATTGCTATTTACAGAGGAAGGTTGAAGCTCCATTCATACCAAAGTTTAGAGGCTCTGG  
AGATACCAGCAACTTTGATGACTATGAAGAAGAAGATATCCGTGTCTCTATAACGAAAAATGTGAAAA  
GAATTTGGTGAATTT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG232566 representing NM\_001242859  
Red=Cloning site Green=Tags(s)

MGLSRKSSDASACSSSEISVKEFLAKAKEDFLKKWENPTQNNAGLEDFERKKTLLGTGSFGRVMLVKHKAT  
 EQYYAMKILDKQKVVVLLKQIEHTLNEKRILQAVNFPFLVRLLEYAFKDNSNLYMVMEEVPPGGEMFSLRRI  
 GRFSEPHARFYAAQIVLTFEYLHSLDLIYRDLKPENLLIDHQGYIQVTDGFGAKRVKGRWTWLCGTPEYL  
 APEIILSKGYNKAVDWWALGVL IYEMAAGYPPFFADQPIQIYEKIVSGKVRFP SHFSSDLKDLLRNLLQV  
 DLTKRFGNLKNGVSDIKTHKWFATTDWIAIYQRKVEAPFIPKFRGSGDTSNFDDYEEEDIRVSITEKCAK  
 EFGF

TRTRPLE - GFP Tag - V


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:**

NM\_001242859

<b>ORF Size:</b>	1065 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001242859.2</a>
<b>RefSeq Size:</b>	4308 bp
<b>RefSeq ORF:</b>	1068 bp
<b>Locus ID:</b>	5567
<b>UniProt ID:</b>	<a href="#">P22694</a>
<b>Cytogenetics:</b>	1p31.1
<b>Protein Families:</b>	Druggable Genome, Protein Kinase
<b>Protein Pathways:</b>	Apoptosis, Calcium signaling pathway, Chemokine signaling pathway, Dilated cardiomyopathy, Gap junction, GnRH signaling pathway, Hedgehog signaling pathway, Insulin signaling pathway, Long-term potentiation, MAPK signaling pathway, Melanogenesis, Olfactory transduction, Oocyte meiosis, Prion diseases, Progesterone-mediated oocyte maturation, Taste transduction, Vascular smooth muscle contraction, Vibrio cholerae infection, Wnt signaling pathway
<b>Gene Summary:</b>	The protein encoded by this gene is a member of the serine/threonine protein kinase family. The encoded protein is a catalytic subunit of cAMP (cyclic AMP)-dependent protein kinase, which mediates signalling through cAMP. cAMP signaling is important to a number of processes, including cell proliferation and differentiation. Multiple alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2014]