

Product datasheet for **RG232686**

AKT2 (NM_001243028) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AKT2 (NM_001243028) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	AKT2
Synonyms:	HIHGH; PKBB; PKBBETA; PRKBB; RAC-BETA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG232686 representing NM_001243028 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGACCGAGAGGCCGCGACCCAACACCTTTGTACATCGCTGCCTGCAGTGGACCACAGTCATCGAGA
GGACCTTCCACGTGGATTCTCCAGACGAGAGGGAGGAGTGGATGCGGGCCATCCAGATGGTCGCAACAG
CCTCAAGCAGCGGGCCCCAGGCGAGGCCCATGGACTACAAGTGTGGCTCCCCAGTGACTCCTCCACG
ACTGAGGAGATGGAAGTGGCGTCAGCAAGGCACGGGCTAAAGTGACCATGAATGACTTCGACTATCTCA
AACTCCTTGGCAAGGAACTTTGGCAAAGTCATCCTGGTGCAGGAGAAGGCCACTGGCCGCTACTACGC
CATGAAGATCCTGCGGAAGGAAGTCATCATTGCCAAGGATGAAGTCGCTCACACAGTCACCGAGAGCCGG
GTCTCCAGAACCAGGCACCCGTTCTCACTGCGCTGAAGTATGCCTCCAGACCCACGACCCGCTGT
GCTTTGTGATGGAGTATGCCAACGGGGGTGAGCTGTTCTTCCACCTGTCCCGGGAGCGTGTCTTACAGA
GGAGCGGGCCCGTTTTATGGTGCAGAGATTGTCTCGGCTCTTGAGTACTTGCACTCGCGGGACGTGGTA
TACCGGCATCAAGCTGGAAAACCTCATGCTGGACAAAGATGGCCACATCAAGTCACTGACTTTGGCC
TCTGAAAAGAGGGCATCAGTGACGGGGCCACCATGAAAACCTTCTGTGGGACCCCGAGTACCTGGCGCC
TGAGGTGCTGGAGACAATGACTATGGCCGGCCGTGGACTGGTGGGGCTGGGTGTGGTCATGTACGAG
ATGATGTGCGGGCCCTGCCCTTCTACAACCAGGACCACGAGCGCCTCTTCGAGTCATCCTCATGGAAG
AGATCCGCTTCCCGGCACGCTCAGCCCCGAGGCCAAGTCCCTGCTTGGTGGGCTGCTTAAGAAGGACCC
CAAGCAGAGGCTTGGTGGGGGGCCAGCGATGCCAAGGAGGTGATGGAGCACAGGTTCTTCTCAGCATC
AACTGGCAGGACGTGGTCCAGAAGAAGCTCCTGCCACCCTTCAAACCTCAGGTACGTCAGGTCGACA
CAAGGTAATTCGATGATGAATTTACCGCCAGTCCATCACAATCACACCCCTGACCGCTATGACAGCCT
GGGCTTACTGGAGCTGGACCAGCGGACCCACTTCCCCAGTTCTCCTACTCGGCCAGCATCCGGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG232686 representing NM_001243028
 Red=Cloning site Green=Tags(s)


MKTERPRPNTFVIRCLQWTTVIERTFHVDSPDEREEWMRAIQMVANSLKQRAPGEDPMDYKCGSPSDSST
 TEEMEVAVSKARAKVTMNDFDYLKLLGKGTFGKVLVREKATGRYYAMKILRKEVIAKDEVAHTVTESR
 VLQNRHPFLTALKYAFQTHDRLCFVMEYANGGELFFHLSRERVFTEERARFYGAEIVSALEYLHSRDVV
 YRDIKLENMLDKDGHIKITDFGLCKEGISDGATMKTFCGTPEYLAPEVLEDNDYGRAVDWVWGLGVVME
 MMCGRLPFYNQDHERLFELILMEEIRFPRTLSPKAKSLLAGLLKKDPKQRLGGGSPDAKEVMEHRFFLSI
 NWQDVVQKLLPPFKPQVTSEVDTRYFDDEFTAQSITITPPDRYDSLGLLELDQRTHFPQFSYASIRE

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_001243028

ORF Size:	1257 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001243028.3
RefSeq Size:	5187 bp
RefSeq ORF:	1260 bp
Locus ID:	208
Cytogenetics:	19q13.2
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase
Protein Pathways:	Acute myeloid leukemia, Adipocytokine signaling pathway, Apoptosis, B cell receptor signaling pathway, Chemokine signaling pathway, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Glioma, Insulin signaling pathway, Jak-STAT signaling pathway, MAPK signaling pathway, Melanoma, mTOR signaling pathway, Neurotrophin signaling pathway, Non-small cell lung cancer, Pancreatic cancer, Pathways in cancer, Progesterone-mediated oocyte maturation, Prostate cancer, Renal cell carcinoma, Small cell lung cancer, T cell receptor signaling pathway, Tight junction, Toll-like receptor signaling pathway, VEGF signaling pathway

Gene Summary:

This gene is a putative oncogene encoding a protein belonging to a subfamily of serine/threonine kinases containing SH2-like (Src homology 2-like) domains, which is involved in signaling pathways. The gene serves as an oncogene in the tumorigenesis of cancer cells. For example, its overexpression contributes to the malignant phenotype of a subset of human ductal pancreatic cancers. The encoded protein is a general protein kinase capable of phosphorylating several known proteins, and has also been implicated in insulin signaling. [provided by RefSeq, Nov 2019]