

## Product datasheet for **RG232692**

### **IKZF3 (NM\_001257410) Human Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	IKZF3 (NM_001257410) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	IKZF3
Synonyms:	AIO; AIOLOS; ZNFN1A3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG232692 representing NM_001257410 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAAGATATACAAACAAATGCGGAAGTAAAAGCACTCAGGAGCAGTCTGTGCCCGCAGAAAAGTGCAG  
CGGTTTTGAATGACTACAGTTTAAACAAATCTCATGAAATGGAAAATGTGGACAGTGGAGAAGGCCCAGC  
CAATGAAGATGAAGACATAGGAGGTGAACGCCATTCCAGTGTAAATCAGTGTGGGCATCTTTACTCAG  
AAAGGTAACCTCCTCCGCCACATTAACATGCACACAGGGGAAAAACCTTTAAGTGTACCTCTGCAACT  
ATGCATGCCAAAGAAGAGATGCGCTCACGGGCATCTTAGGACACATTCTGTGGAGAAACCTACAATG  
TGAGTTTTGTGGAAGGAGTTACAAGCAGAGAAGTCCCTTGAGGAGCACAAAGGAGCGCTGCCGTACATTT  
CTTCAGAGCACTGACCCAGGGGACACTGCAAGTGCAGGCAAGACACATCAAAGCAGAGATGGGAAGT  
AAAGAGCTCTCGTACTGGACAGATTAGCAAGCAATGTGGCAAAACGAAAAAGCTCAATGCCTCAGAAATT  
CATTGGTGAAGCGCCACTGCTTTGATGTCAACTATAATTCAAGTTACATGTATGAGAAAGAGAGTGAG  
CTCATACAGACCCGCATGATGGACCAAGCCATCAATAACGCCATCAGCTATCTTGCGCCGAAAGCCCTGC  
GCCCTTGGTCCAGACACCGCTGCTCCACCTCGGAGATGGTCCAGTTATCAGCAGCATGTATCCCAT  
AGCCCTCACCCGGGCTGAGATGTCAAACGGTCCCCCTCAAGAGCTGAAAAAGAAAAGCATCCACCTTCCA  
GAGAAGAGCGTGCCTTCTGAGAGAGGCCCTCTCTCCCAACAATAGTGCCACGACTCCACGGACACTGACA  
GCAACCATGAAGAACGCCAGAATCACATCTATCAGCAAAATCACATGGTCTGTCTCGGGCCCGCAATGG  
GATGCCACTTCTGAAGGAGGTTCCCGCTCTTACGAACTCCTCAAGCCCCGCCATCTGCCAAGAGAC  
TCCGTCAAAGTGTCAACAAGGAAGGGGAGGTGATGGATGTGTATCGGTGTGACCACTGCCGCGTCTCT  
TCCTGGACTATGTGATGTTACGATTACATGGGCTGCCACGGCTCCGTGACCCTTTCGAGTGTAAACAT  
GTGTGGATATCGAAGCCATGATCGGTATGAGTTCTCGTCTCACATAGCCAGAGGAGAACACAGAGCCCTG  
CTGAAG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG232692 representing NM\_001257410  
 Red=Cloning site Green=Tags(s)

MEDIQTNALKSTQEQSVAESA AVLNDYSLTKSHEMENVDSGEGPANEDDIGGERPFQCNQCGASFTQ  
 KGNLLRHIKLTGEKPFKCHLCNYACQRRDALTGHLRTHSVEKPYKCEFCGRSYQRSSLEEHEKRCRTF  
 LQSTDPGDTASAEARHIKAEMGSERALVLDRLASNVAKRKSMPQKF IGEKRHCFD VNYNSSYMEKESE  
 LIQTRMMDQAINNAISYLGAEALRPLVQTPAPTSEMPVPISSMYPIALTRAEMSGAPQELEKKS IHLP  
 EKSVPSE RGLSPNNSGHDSTDTDSNHEERQNHIIYQQNHMVL SRARNGMPLLKEVPRSYELLKPPPI CPRD  
 SVKVINKEGEVMDVYRCDHCRVLF LDYVMFTIHMGC HGRDPFECNMCGYRSHDRYEFSSHIARGEHRAL  
 LK

TRTRPLE - GFP Tag - V


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:**

NM\_001257410

<b>ORF Size:</b>	1266 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001257410.2</a>
<b>RefSeq Size:</b>	9425 bp
<b>RefSeq ORF:</b>	1269 bp
<b>Locus ID:</b>	22806
<b>UniProt ID:</b>	<a href="#">Q9UKT9</a>
<b>Cytogenetics:</b>	17q12-q21.1
<b>Gene Summary:</b>	This gene encodes a member of the Ikaros family of zinc-finger proteins. Three members of this protein family (Ikaros, Aiolos and Helios) are hematopoietic-specific transcription factors involved in the regulation of lymphocyte development. This gene product is a transcription factor that is important in the regulation of B lymphocyte proliferation and differentiation. Both Ikaros and Aiolos can participate in chromatin remodeling. Regulation of gene expression in B lymphocytes by Aiolos is complex as it appears to require the sequential formation of Ikaros homodimers, Ikaros/Aiolos heterodimers, and Aiolos homodimers. Several alternative transcripts encoding different isoforms have been described, as well as some non-protein coding variants. [provided by RefSeq, Apr 2012]