

## Product datasheet for **RG232694**

### **PAX6 (NM\_001258465) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** PAX6 (NM\_001258465) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** PAX6  
**Synonyms:** AN; AN1; AN2; ASGD5; D11S812E; FVH1; MGDA; WAGR  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG232694 representing NM\_001258465  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGCAGAACAGTCACAGCGGAGTGAATCAGCTCGGTGGTGTCTTTGTCAACGGCGGCCACTGCCGGACT  
 CCACCCGGCAGAAGATTGTAGAGCTAGCTCACAGCGGGGCCCGGCGTGCACATTTCCCGAATTCTGCA  
 GGTGTCCAACGGATGTGTGAGTAAAATTCTGGGCAGGTATTACGAGACTGGCTCCATCAGACCCAGGGCA  
 ATCGGTGGTAGTAAACCGAGAGTAGCGACTCCAGAAGTTGTAAGCAAAAAGCCAGTATAAGCGGGAGT  
 GCCCGTCCATCTTTGCTTGGGAAATCCGAGACAGATTACTGTCCGAGGGGTCTGTACCAACGATAACAT  
 ACCAAGCGTGTCAATAAACAGAGTTCTTCGCAACCTGGCTAGCGAAAAGCAACAGATGGGCGCAGAC  
 GGATGTATGATAAACTAAGGATGTTGAACGGGACAGCCGGAAGCTGGGGCACCCGCCCTGGTTGGTATC  
 CGGGGACTTCGGTGCCAGGGCAACCTACGCAAGATGGCTGCCAGCAACAGGAAGGAGGGGAGAGAATAC  
 CACTCCATCAGTTCCAACGGAGAAGATTCAGATGAGGCTCAAATGCGACTTCAGCTGAAGCGGAAGCTG  
 CAAAGAAATAGAATCCTTTACCAAGAGCAAATGAGGCCCTGGAGAAAGAGTTTGAGAGAACCATT  
 ATCCAGATGTGTTGCCGAGAAAGACTAGCAGCAAAAATAGATCTACCTGAAGCAAGAATACAGGTATG  
 GTTTTCTAATCGAAGGGCCAAATGGAGAAGAGAAGAAAACTGAGGAATCAGAGAAGACAGGCCAGCAAC  
 ACACCTAGTCATATTCTATCAGCAGTAGTTTCAGCACCAGTGTCTACCAACCAATCCACAACCCACCA  
 CACCGGTTTCCTCCTCACATCTGGCTCCATGTTGGCCGAACAGACACAGCCCTCACAAACACCTACAG  
 CGCTCTGCCGCCTATGCCAGCTTCACCATGGCAATAACCTGCCTATGCAACCCCCAGTCCCCAGCCAG  
 ACCTCCTCATACTCCTGCATGCTGCCACCAGCCCTTCGGTGAATGGGCGGAGTTATGATACCTACACCC  
 CCCCACATATGCAGACACACATGAACAGTCAGCCAATGGGCACCTCGGGCACCCTTCAACAGGACTCAT  
 TCCCCTGGTGTTCAGTTCCAGTTCAAGTTCCCGGAAGTGAACCTGATATGTCTCAATACTGGCCAAGA  
 TTACAG

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG232694 representing NM\_001258465  
 Red=Cloning site Green=Tags(s)

MQNSHSGVNLGGVFVNGRPLPDSTRQKIVELAHSGARPCDISRILQVSNCGVSKILGRYYETGSIRPRA  
 IGGSKPRVATPEVYSKIAQYKRECPISIFAWIIRDRLLSEGVCNTNDNIPSVSSINRVLRLNLASEKQQMGAD  
 GMYDKLRMLNGQTGSWGTTRPGWYPTGTPVGPQPTQDGCQQQEGGGENTNSISSNGEDSDEAQMRLQLKRKL  
 QRNRTSFTQEIEALEKEFERHTHYPDVFAERERLAAKIDLPEARIQVWFSNRRAKWRREEKLRNQRQASN  
 TPSHIPISSSFSTSVYQPIPPQTPVSSFTSGSMLGRDTDALTNTYSALPPMPSTMANLPMQPPVPSQ  
 TSSYSCLPTSPSVNGRSYDTPPHMQTHMSQPMGTSGTTSTGLISPGVSVPVQVPGSEPDMSQYWPRLQ


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001258465

<b>ORF Size:</b>	1266 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001258465.3</a>
<b>RefSeq Size:</b>	6857 bp
<b>RefSeq ORF:</b>	1269 bp
<b>Locus ID:</b>	5080
<b>UniProt ID:</b>	<a href="#">P26367</a>
<b>Cytogenetics:</b>	11p13
<b>Protein Families:</b>	Adult stem cells, Druggable Genome, Embryonic stem cells, Transcription Factors
<b>Protein Pathways:</b>	Maturity onset diabetes of the young
<b>Gene Summary:</b>	<p>This gene encodes paired box protein Pax-6, one of many human homologs of the Drosophila melanogaster gene <i>prd</i>. In addition to a conserved paired box domain, a hallmark feature of this gene family, the encoded protein also contains a homeobox domain. Both domains are known to bind DNA and function as regulators of gene transcription. Activity of this protein is key in the development of neural tissues, particularly the eye. This gene is regulated by multiple enhancers located up to hundreds of kilobases distant from this locus. Mutations in this gene or in the enhancer regions can cause ocular disorders such as aniridia and Peter's anomaly. Use of alternate promoters and alternative splicing results in multiple transcript variants encoding different isoforms. Interestingly, inclusion of a particular alternate coding exon has been shown to increase the length of the paired box domain and alter its DNA binding specificity. Consequently, isoforms that carry the shorter paired box domain regulate a different set of genes compared to the isoforms carrying the longer paired box domain. [provided by RefSeq, Mar 2019]</p>