

Product datasheet for **RG232732**

5HT3E (HTR3E) (NM_198313) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	5HT3E (HTR3E) (NM_198313) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	HTR3E
Synonyms:	5-HT3-E; 5-HT3c1; 5-HT3E
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG232732 representing NM_198313 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAGGAAGCTGGTTCACAGGAAAAGATTTTCCTTCTACCTCCTTCTCGGTTTTCTGCTTCAAGGAA
GGGGCGTTACTTTACCATCAATTGCTCAGGGTTTGGCCAGCACGGGGCGGATCCCCTGCTCTGAATTC
AGTGTAAATAGAAAGCCCTTCGTCGGTCCCAACATCAGCGTCCCCACCAAGTCAACATCTCCTTC
GCGATGCTGCCATCCTAGATGTGGTTGGGATAACCCATTTATCAGCTGGAACCCAGAGGAATGTGAGG
GCATCACGAAGATGAGTATGGCAGCCAAGAACCCTGTGGCTCCAGACATTTTCATCATTGAACCTCATGGA
TGTGGATAAGACCCAAAAGGCCTCACAGCATATGTAAGTAATGAAGGTCGCATCAGGTATAAGAAACCC
ATGAAGGTGGACAGTATCTGTAACCTGGACATCTTCTACTTCCCCTTCGACCAGCAGAAGTGCACACTCA
CCTTCAGCTCATTCTCTACACAGTGGACAGCATGTTGCTGGACATGGAGAAAGAAGTGTGGGAAATAAC
AGACGCATCCCGGAACATCCTTCAGACCATGGAGAATGGGAGCTCCTGGGCCTCAGCAAGGCCACCCGA
AAGTTGTCAGGGGAGGCAACCTGTATGATCAGATCGTGTCTATGTGGCCATCAGGCGCAGGCCAGTC
TCTATGCATAAACCTTCTCGTGCCAGTGGCTTTCTGGTTGCCATCGATGCCCTCAGCTTCTACCTGCC
AGTGAAAAGTGGGAATCGTGTCCATTCAAGATAACGCTCCTGCTGGGCTACAACGTCTTCTGCTCATG
ATGAGTGACTTGTCCCCACAGTGGCACCCCTCATCGGTGCTACTTCGCCCTGTGCCTGTCCCTGATG
TGGTGGGCAGCTGCTGGAGACCATTTTCATCACCCACCTGCTGCACGTGGCCACCACCCAGCCCCACC
CCTGCCTCGGTGGCTCCACTCCCTGCTGCTCCACTGCAACAGCCCGGGGAGATGCTGTCCACTGCGCCC
CAGAAGGAAAATAAGGGCCCGGTCTACCCCCACCCACCTGCCCGGTGTGAAGGAGCCAGAGGTATCAG
CAGGGCAGATGCCGGCCCTGCGGAGGCAGAGCTGACAGGGGCTCAGAAATGGACAAGGGCCAGCGGGA
ACACGAGGCCAGAAAGCAGCACTCAGTGGAGCTGTGGTTGCAGTTCAGCCACGGATGGACGCCATGCTC
TCCGCTCTACCTGCTTTCATGGCTCCTCTATCATCACCGTCATATGCCTCTGGAACACC

ACCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG232732 representing NM_198313
Red=Cloning site Green=Tags(s)


MEGSWFHRKRFSFYLLLGFLQGRGVFTTINCSGFGQHGDPTALNSVFNRPFRPVNTISVPTQVNI
 AMSAILDVVWDNPFISWNPEECEGITKMSMAAKNLWLPDIFI IELMDVDKTPKGLTAYVSN
 EGRIRYKPKMKVDSICNLDIFYFPDQQNCTLTFSSFLYTVDSMLLDMEKEVWEITDASRNILQTHGEWELLGLSKATA
 KLSRGGNLYDQIVFYVAIRRRPSLYVINLLVPSGFLVAIDALSFYLPVKSNGRVPFKITLLLGYNVFLLM
 MSDLLPTSGTPLIGVYFALCLSLMVGSLLETIFITHLLHVATTQPPPLRWLHLLLLHCNSPGRCCPTAP
 QKENKGPLTPTHLPGVKEPEVSAGQMPGPAEAEELTGGSEWTRAQREHEAQKQHSVELWLQF
 SHAMDAMLFRLLYLLFMASIIITVICLWNT

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_198313

ORF Size:	1323 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_198313.2 , NP_938055.1
RefSeq Size:	2109 bp
RefSeq ORF:	1326 bp
Locus ID:	285242
UniProt ID:	A5X5Y0
Cytogenetics:	3q27.1
Protein Families:	Druggable Genome, Ion Channels: Cys-loop Receptors, Transmembrane
Gene Summary:	This locus encodes a 5-hydroxytryptamine (serotonin) receptor subunit. The encoded protein, subunit E, may play a role in neurotransmission in myenteric neurons. Genes encoding subunits C, D and E form a cluster on chromosome 3. Alternatively spliced transcript variants encoding distinct isoforms have been described. [provided by RefSeq, Feb 2012]