

Product datasheet for **RG232735**

LRRC31 (NM_001277127) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRC31 (NM_001277127) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	LRRC31
Synonyms:	HEL-S-293
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG232735 representing NM_001277127 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGAGTCAAACAAGGAAGAAAACCTTCCTCAGAAGGAGAACTAAGCCCCAGACTTCAACTGTCAACAAAT
TTCTCAGGGGCTCCAATGCTGAAAGCAGAAAAGAGGACAATGACCTTAAACAAGTGATCCCAACCCAG
CGACTGGATACAGAAGACAGCCACCTCAGAGACTGCTAAGCCTCTCAGTTCAGAAATGGAATGGAGATCC
AGTATGGAGAAAAATGAGCATTTCTGCAGAAGCTGGGCAAAAAGGCTGTCAACAAGTGCTAGATTGGA
ATAACTGTGGATTAACAACAGCGGACATGAAAGAAATGGTTGCCTTGCTGCCTTTTCTCCAGACTTGGA
AGAACTGGATATCTCCTGGAATGGTTTTGTAGGTGGAACCTCCTTTCCATCACTCAGCAATGCATCTG
GTCAGCAAGTTAAAAATCTTGAGGCTGGGTAGCTGCAGACTCACCACTGACGATGTTCAAGCACTGGGAG
AAGCATTTGAGATGATTCCTGAACCTTGAAGAGCTAAATTTGTCTTGAACAGTAAAGTGGGAGGAAATTT
GCCTCTGATCCTTCAGAAGTTCCAAAAGGGAGCAAGATACAAATGATTGAGCTTGTGGATTGCTCCCTC
ACGTCAGAAGATGGGACATTTCTGGGTCAACTGCTACCTATGCTGCAAGTCTCGAAGTACTTGATCTTT
CCATTACAGAGACATTGTTGGCAGTCTGAACAGTATTGCTCAGGGATTAAGAAGCACCTCAAATCTGAA
AGTACTGAAGTTACATTCATGTGGATTATCAGAAAGAGTGTCAAATATTGGATGCTGCTTTTAGGTAT
TTGGGTGAGCTGAGGAAATTAGATCTTTCTGCAATAAGGATCTAGGTGGAGGTTTGAAGACTCGCCGG
CTCAGTTGGTCATGCTAAAGCATCTACAAGTCTAGATCTTACCAGTGCTCACTAACAGCAGATGACGT
GATGTCACTGACCCAGGTCATTCTTTACTTTCAATCTTCAAGAATTGATTTATCAGCCAACAAAAAG
ATGGGCAGTTCTTCTGAAAACCTTACTCAGCAGGCTCCGATTTTTACCAGCATTGAAGTCATTAGTTATCA
ACAACGTGCTTTGGAGAGTGAGACTTTTACAGCTCTTGTGAAGCCTCTGTTACCTCTCTGCTCTGGA
AGTATTCAACCTTTCTTGAACAAGTGTGTTGGTGGCAACTTGAAGCTGCTTCTGGAAACACTAAAGCTT
TCCATGTCTCTTCAAGTCTGAGGCTGAGCAGCTGTTCCCTGGTGACAGAGGATGTGGCTCTCTGGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG232735 representing NM_001277127
 Red=Cloning site Green=Tags(s)


MSQTRKKTSSGETKPQTSTVNKFLRGSNAESRKEDNDLKTSDSQPSDWIQKTATSETAKPLSSEMEWRS
 SMEKNEHFLQKLGGKAVNKCLDLNCCGLTTADMKEMVALLPFLPDLEELDISWNGFVGGTLLSITQQMHL
 VSKLKIIRLGLSCLRTDDVQALGEAFEMIELEELNLWNKSVGGNPLILQKFQKGSKIQMIELVDCSL
 TSEDGTFGLQLPMLQSLEVLDSLNRDIVGSLNSIAQGLKSTSNLKVLLHSCGLSQKSVKILDAAFRY
 LGELRKLDLSCNKDLGGGFEDSPAQLVMLKHLQVLDLHQCSLTADDVMSLTQVPIPLLSNLQELDLANKK
 MGSSSENLLSRLRFLPALKSLVINNCALESETFTALAEASVHLSALEVFNLWNKCVGGNKLKLLLETCLK
 SMSLQVLRLLSSCLVTEVALLE

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001277127

ORF Size: 1329 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001277127.2](#)

RefSeq Size: 2561 bp


RefSeq ORF: 1332 bp

Locus ID: 79782

UniProt ID: [Q6UY01](#)

Cytogenetics: 3q26.2

Product images:



Circular map for RG232735