

## Product datasheet for **RG232762**

### TEA domain family member 2 (TEAD2) (NM\_001256661) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	TEA domain family member 2 (TEAD2) (NM_001256661) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TEAD2
Synonyms:	ETF; TEAD-2; TEF-4; TEF4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG232762 representing NM\_001256661  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGGGGGAACCCCGGGCTGGGGCCGCCCTGGACGATGGCAGCGGCTGGACGGCAGTGAGGAAGGCAGTG  
 AGGAGGGTACCGCGGCAGTGAGGGGGCTGGGGGTGACGGGGGCCCGGATGCAGAGGGGGTGTGGAGCCC  
 AGACATTGAGCAGAGCTTCCAGGAGGCCCTGGCCATCTATCCACCCTGCGGCCGCCGAAAAATAATTTTG  
 TCTGATGAAGGCAAGATGTATGGTGGAAATGAACTGATCGCCCGTACATCAAGCTGAGAACGGGGAAGA  
 CCCGAACCGAAAAACAGGTTTCTAGTCACATCCAGTTTTGGCCCGAAGGAAATCAAGGAAATCCAGTC  
 CAAGTTGAAGGCTCTGAACGTGGACCAGGTTTCCAAGGACAAGGCTTCCAGACAATGGCAACCATGTCC  
 TCTGCCAGCTCATCTCCGCGCTTCTCTGCAGGCCAACTGGGTCCTGCTCAGGCCTCTGAGC  
 TTTTCCAGTTTTGGTCTGGAGGATCTGGCCCCCTGGAATGTTCCAGATGTGAAGCCATTCTCACAGAC  
 ACCGTTACCTTGTCACTGACTCCCCATCTACTGACCTCCAGGGTACGAGCCCCCAAGCCCTCTCA  
 CCCTGCCCCACCTACCCATCGCCCCAGCCTGGCAGGCTCGGGGCTGGGCACCGCCCGTTGCAGC  
 TGGTAGAGTTCTCAGCCTTCGTGGAACCGCCAGATGCAGTTGATTCTTACCAGAGGCACCTGTTCTGTGA  
 CATCAGCCAGCACTGCCCCAGCCCCGGAGCGCCCGCTCGAGAGTGTGGACGTCCGGCAGATCTACGAC  
 AAATTCCTGAGAAAAAGGGTGGCTCCGAGAGCTATATGATCGTGGCCCCCCCATGCCTTCTTCTCTGG  
 TCAAGTTCTGGGCGACCTGAACGGGGCCAAGTGGTGGAGGAGCGGGGCCGGTGGCAGCATCAGCAG  
 TGGTGGCTTCTACGGAGTGAGCAGCCAGTATGAGAGCCTGGAACACATGACCCTCACCTGTTCTCTCAAG  
 GTCTGCTCTTTGGCAAGCAGGTGGTGGAGAAGGTGGAGACGGAACGGGCCAGCTGGAGGACGGCAGAT  
 TTGTGTACCGCTGCTGCGCTCGCCATGTGCGAGTACCTGGTGAATTTCTTGACAAGTTGCGGCAGCT  
 GCCTGAGCGATACATGATGAACAGCGCTCTGGAAAACCTCACCATCCTCCAGGTGGTGACAAACAGAGAC  
 ACCCAGGAACTGCTGCTCTGCACCGCTATGTCTTCGAGGTCTCCACCAGCGAGCGTGGGGCCAGCATC  
 ACATTTACCGCTGGTCAAGGAC

**ACGCGT**ACGCGGCCGCTCGAG – GFP Tag – GTTTAA

**Protein Sequence:**

>RG232762 representing NM\_001256661  
 Red=Cloning site Green=Tags(s)

MGEPRAGAALDDGSGWTGSEEGSEEGTGGSEAGGDDGGPDAEGVWSPDIEQSFQEALAIYPPCGRRKIIL  
 SDEGKMYGRNELIARYIKLRTGKTRTRKQVSSHQVLARRKSREIQSKLKALNVDQVSKDAFQTMATMS  
 SAQLISAPSLQAKLGPTGPOASELFQFWSGGSGPPWNPVDPKPFSTPFTLSLTPPSTDLPGYEPQALS  
 PLPPPTSPPAWQARGLGTARLQLVEFSAFVEPPDAVDSYQRHLFVHISQHCPSGAPPLESVDVRIYD  
 KFPEKKGGLREL YDRGPPHAFVLKFWADLNWGPSGEEAGAGGSISSGGFYGVSSQYESLEHMTLTCSSK  
 VCSFGKQVVEKVEVERAQLEDGRFVYRLLRSPMCEYLVNFLHKLRLQLPERYMMNSVLENFTILQVVTNRD  
 TQELLLCTAYVFEVSTSERGAQHIIYRLVRD

**TRTRPLE** – GFP Tag – V


**Restriction Sites:**

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001256661

ORF Size: 1353 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001256661.2</a>
<b>RefSeq Size:</b>	2211 bp
<b>RefSeq ORF:</b>	1356 bp
<b>Locus ID:</b>	8463
<b>UniProt ID:</b>	<a href="#">Q15562</a>
<b>Cytogenetics:</b>	19q13.33
<b>Protein Families:</b>	Transcription Factors
<b>Gene Summary:</b>	<p>Transcription factor which plays a key role in the Hippo signaling pathway, a pathway involved in organ size control and tumor suppression by restricting proliferation and promoting apoptosis. The core of this pathway is composed of a kinase cascade wherein MST1/MST2, in complex with its regulatory protein SAV1, phosphorylates and activates LATS1/2 in complex with its regulatory protein MOB1, which in turn phosphorylates and inactivates YAP1 oncoprotein and WWTR1/TAZ. Acts by mediating gene expression of YAP1 and WWTR1/TAZ, thereby regulating cell proliferation, migration and epithelial mesenchymal transition (EMT) induction. Binds to the SPH and GT-IIC 'enhansons' (5'-GTGGAATGT-3'). May be involved in the gene regulation of neural development. Binds to the M-CAT motif.</p> <p>[UniProtKB/Swiss-Prot Function]</p>