

## Product datasheet for **RG232805**

### Estrogen Receptor beta (ESR2) (NM\_001214903) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Estrogen Receptor beta (ESR2) (NM_001214903) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ESR2
Synonyms:	ER-BETA; Erb; ESR-BETA; ESRB; ESTRB; NR3A2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG232805 representing NM\_001214903  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGATATAAAAACTCACCATCTAGCCTTAATTCCTCTCCTACAACCTGCAGTCAATCCATCTTAC  
 CCCTGGAGCACGGCTCCATATACATACCTCCTCTATGTAGACAGCCACCATGAATATCCAGCCATGAC  
 ATTCTATAGCCCTGCTGTGATGAATTACAGCATTCCCAGCAATGCACTAACTTGGAAAGTGGGCCTGGT  
 CGGCAGACCACAAGCCAAATGTGTTGTGGCCAACACCTGGGCACCTTTCTCCTTTAGTGGTCCATCGCC  
 AGTTATCACATCTGTATGCGGAACCTCAAAGAGTCCCTGGTGTGAAGCAAGATCGCTAGAACACACCTT  
 ACCTGTAAACAGAGAGACTGAAAAGGAAGGTTAGTGGGAACCGTTGCGCCAGCCCTGTTACTGGTCCA  
 GGTTCAAAGAGGGATGCTCACTTCTGCGCTGTCTGCAGCATTACGCATCGGGATACACTATGGAGTCT  
 GGTCTGTGAAGGATGTAAGGCCTTTTTAAAAGAAGCATTCAAGGACATAATGATTATATTTGTCCAGC  
 TACAAATCAGTGTACAATCGATAAAAACCGGCGCAAGAGCTGCCAGGCCCTGCCGACTTCGGAAGTGTAC  
 GAAGTGGGAATGGTGAAGTGTGGCTCCCGAGAGAGAGATGTGGGTACCGCCTTGTGCGGAGACAGAGAA  
 GTGCCGACGAGCAGCTGCACTGTGCCGGCAAGGCCAAGAGAAGTGGCGGCCACGCGCCCCGAGTGGGGGA  
 GCTGCTGCTGGACGCCCTGAGCCCCGAGCAGCTAGTGCTCACCCCTCCTGGAGGCTGAGCCGCCCCATGTG  
 CTGATCAGCCGCCCAGTGCGCCCTTACCAGGCCTCCATGATGATGTCCTGACCAAGTTGGCCGACA  
 AGGAGTTGGTACACATGATCAGCTGGGCCAAGAAGATCCCGGCTTTGTGGAGCTCAGCCTGTTCCGACCA  
 AGTGGCGCTCTTGAGAGCTGTTGGATGGAGGTGTTAATGATGGGGCTGATGTGGCGTCAATTGACCAC  
 CCCGGCAAGCTCATCTTTGCTCCAGATCTGTTCTGGACAGGGATGAGGGGAAATGCGTAGAAGGAATTC  
 TGGAAATCTTTGACATGCTCCTGGCAACTACTCAAGGTTTCGAGAGTTAAAACCAACACAAAGAATA  
 TCTCTGTGTCAAGGCCATGATCCTGCTCAATTCCAGTATGTACCCTCTGGTCACAGCGACCCAGGATGCT  
 GACAGCAGCCGGAAGCTGGCTCACTTGCTGAACGCCGTGACCGATGCTTTGGTTTGGGTGATTGCCAAGA  
 GCGGCATCTCTCCAGCAGCAATCCATGCGCCTGGCTAACCTCCTGATGCTCCTGTCCACGTCAGGCA  
 TCGGAGGTACGCGCCC

**ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA**

**Protein Sequence:**

>RG232805 representing NM\_001214903  
 Red=Cloning site Green=Tags(s)


MDIKNSPSSLNSPSSYNCSQSILPLEHGSIIYIPSSYVDSHHEYPAMTFYSPAVMNYIIPSNVTNLEGGPG  
 RQTTSPNVLWPTPGHLSPLVVHRQLSHLYAEPQKSPWCEARSLEHTLPVNRETLKRKVSNGRCASPVTGP  
 GSKRDAHFCAVCSDYASGYHYGVWSCEGCKAFFKRSIQGHNDYICPATNQCTIDKNRRKSCQACRLRKCY  
 EVGMVKCGSRRERCGYRLVRRQRSADQLHCAGKAKRSGGHAPRVRELLLDALSPEQLVLTLLAEPPHV  
 LISRPSAPFTEASMMMSLTKLADKELVHMISWAKKIPGFVELSLFDQVRLLESCWMEVLMGLMWRSIDH  
 PGKLIFAPDLVLDREDEGKCEGILEIFDMLLATTSRFRELKQLQKEYLVCAMILLNSSMYPLVTATQDA  
 DSSRKLALLNAVTDALVWVIAKSGISSQQSMRLANLLMLLSHVRHARYAP

**TRTRPLE - GFP Tag - V**

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**

**Plasmid Map:**


**ACCN:** NM\_001214903

**ORF Size:** 1416 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001214903.1, NP_001201832.1</u>
<b>RefSeq Size:</b>	1579 bp
<b>RefSeq ORF:</b>	1418 bp
<b>Locus ID:</b>	2100
<b>Cytogenetics:</b>	14q23.2-q23.3
<b>Protein Families:</b>	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors
<b>Gene Summary:</b>	<p>This gene encodes a member of the family of estrogen receptors and superfamily of nuclear receptor transcription factors. The gene product contains an N-terminal DNA binding domain and C-terminal ligand binding domain and is localized to the nucleus, cytoplasm, and mitochondria. Upon binding to 17beta-estradiol or related ligands, the encoded protein forms homo- or hetero-dimers that interact with specific DNA sequences to activate transcription. Some isoforms dominantly inhibit the activity of other estrogen receptor family members. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been fully characterized. [provided by RefSeq, Jul 2008]</p>