

## Product datasheet for **RG233060**

### **IL31RA (NM\_001242638) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	IL31RA (NM_001242638) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	IL31RA
Synonyms:	CRL; CRL3; GLM-R; GLMR; GPL; hGLM-R; IL-31RA; PLCA2; PRO21384; zcytoR17
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:


>RG233060 representing NM\_001242638  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCCGCATCGCC

ATGAAGCTCTCTCCCAGCCTTCATGTGTAACTGGGGATGATGTGGACCTGGGCACTGTGGATGCTCC  
CCTCACTCTGCAAATTCAGCCTGGCAGCTCTGCCAGCTAAGCCTGAGAACATTCCTGTGTCTACTACTA  
TAGGAAAAATTTAACCTGCACCTGGAGTCCAGGAAAGGAAACCAGTTATACCCAGTACACAGTTAAGAGA  
ACTTACGCTTTTGGAGAAAAACATGATAATTGTACAACCAATAGTTCTACAAGTAAAAATCGTGCTTCGT  
GCTCTTTTTCTTCCAAGAATAACGATCCAGATAATTATACCATTGAGGTGGAAGCTGAAAAATGGAGA  
TGGTGTAAATTAATCTCATATGACATACTGGAGATTAGAGAACATAGCGAAAACTGAACCACTAAGATT  
TTCCGTGTGAAACCGTTTTGGGCATCAAACGAATGATTCAAATTGAATGGATAAAGCCTGAGTTGGCGC  
CTGTTTCATCTGATTTAAAATACACACTTCGATTCAAGACAGTCAACAGTACCAGCTGGATGGAAGTCAA  
CTTCGCTAAGAACCGTAAGGATAAAAACCAACGTACAACCTCACGGGGCTGCAGCCTTTTACAGAATAT  
GTCATAGCTCTGCGATGTGCGGTCAAGGAGTCAAAGTTCTGGAGTGACTGGAGCCAAAGAAAAATGGGAA  
TGACTGAGGAAGAAGCTCCATGTGGCCTGGAACCTGTGGAGAGTCTGAAACCAGCTGAGGCGGATGGAAG  
AAGGCCAGTGCAGTTGTTATGGAAGAAGGCAAGAGGAGCCCCAGTCTAGAGAAAACACTTGGCTACAAC  
ATATGGTACTATCCAGAAAGCAACACTAACCTCACAGAAACAATGAACACTACTAACCAGCAGCTTGAAC  
TGCATCTGGGAGGCGAGAGCTTTTGGGTGTCTATGATTTCTTATAATTCTCTTGGGAAGTCTCCAGTGCC  
CACCTGAGGATTCAGCTATTCAGAAAAATCATTTCAGTGCATTGAGGTGATGCAGGCCTGCGTTGCT  
GAGGACCAGCTAGTGGTGAAGTGGCAAAGCTCTGCTCTAGACGTGAACACTTGGATGATTGAATGGTTTC  
CGGATGTGGACTCAGAGCCCACCACCTTTCTGGGAATCTGTGTCTCAGGCCACGAACTGGACGATCCA  
GCAAGATAAAATTAACCTTTCTGGTGTATAACATCTCTGTGTATCCAATGTTGCATGACAAAGTTGGC  
GAGCCATATTCCATCCAGGCTTATGCCAAAGAAGGCGTTCCATCAGAAGTCTGAGACCAAGGTGGAGA  
ACATTGGCGTGAAGACGGTCACGATCACATGGAAAAGAGATTCCAAGAGTGAGAGAAAGGTATCATCTG  
CAACTACACCATCTTTTACCAAGCTGAAGGTGGAAAAGGATTCTCCAAGACAGTCAATTCAGCATCTTG  
CAGTACGGCCTGGAGTCCCTGAAACGAAAGACCTTTACATTGTTTCAGGTGATGGCCAGCACCAGTGCTG  
GGGAACCAACGGGACCAGCATAAATTTCAAGACATTGTCATTGAGTGTCTTTGAGATTATCCTCATAAC  
TTCTCTGATTGGTGGAGGCTTTCTATTCTCATTATCCTGACAGTGGCATATGGTCTCAAAAAACCAAC  
AAATTGACTCATCTGTGTTGGCCACCCTTCCCAACCCTGCTGAAAGTAGTATAGCCACATGGCATGGAG  
ATGATTTCAAGGATAAGCTAAACCTGAAGGAGTCTGATGACTCTGTGAACACAGAAGACAGGATCTTAAA  
ACCATGTTCCACCCAGTGACAAGTTGGTGAATGACAAGTTGGTGGTGAACCTTTGGGAATGTTCTGCAA  
GAAATTTTACAGATGAAGCCAGAACGGGTGAGAAAACAATTTAGGAGGGGAAAAGATGGGACTAGAA  
TTCTGTCTTCTGCCCACTTCAATA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



**Plasmid Map:**


**ACCN:** NM\_001242638

**ORF Size:** 1986 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001242638.2](#)

**RefSeq Size:** 3036 bp

**RefSeq ORF:** 1989 bp

**Locus ID:** 133396

**UniProt ID:** [Q8NI17](#)

**Cytogenetics:** 5q11.2

**Protein Families:** Druggable Genome, Transmembrane

**Gene Summary:** The protein encoded by this gene belongs to the type I cytokine receptor family. This receptor, with homology to gp130, is expressed on monocytes, and is involved in IL-31 signaling via activation of STAT-3 and STAT-5. It functions either as a monomer, or as part of a receptor complex with oncostatin M receptor (OSMR). Several alternatively spliced transcript variants encoding different isoforms have been noted for this gene.[provided by RefSeq, Jun 2011]