

Product datasheet for **RG233281**

SRP19 (NM_001204194) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SRP19 (NM_001204194) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: SRP19
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG233281 representing NM_001204194
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTTGCCTGCCGCGCGGTCCCCGGCCGACCAGGACAGGTTATTTGTATCTATCCTGCTTATTTAA
ATAATAAGAAGACCATCGCAGAGGGAAGGCGAATCCCCATAAGTAAGAAAAATAAAATGTA CTCTAGAGA
ATGGAATCGTGATGTCCAATACAGAGGCAGAGTCCGGGTCCAGCTCAAACAGGAAGATGGGAGCCTCTGC
CTTGTACAGTTCCCATCACGTAAGTCAGTAATGTTGTATGCAGCAGAAATGATACCTAACTAAAAACAA
GGACACAAAAACAGGAGGTGCTGACCAAAGTCTTCAACAAGGAGAGGGAAGTAAAAAGGGAAAGGAAA
GAAAAAGAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG233281 representing NM_001204194
Red=Cloning site Green=Tags(s)
MACAAARSPADQDRFIC IYPAYLNNKKTIAEGRRIPI SKKNKMSREWNRDVQYRGRVVRVQLKQEDGSLC
LVQFSPSRKSVMLYAAEMIPK LKTRTQKTGGADQSLQQEGEGSKKGKGGKKK

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001204194

ORF Size: 360 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204194.1](#), [NP_001191123.1](#)

RefSeq Size: 1430 bp

RefSeq ORF:	363 bp
Locus ID:	6728
UniProt ID:	P09132
Cytogenetics:	5q22.2
Protein Pathways:	Protein export
Gene Summary:	Signal-recognition-particle assembly, binds directly to 7S RNA and mediates binding of the 54 kDa subunit of the SRP.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG233281