

## Product datasheet for **RG233290**

### **SPINK2 (NM\_001271718) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** SPINK2 (NM\_001271718) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** SPINK2  
**Synonyms:** HUSI-II; SPGF29  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG233290 representing NM\_001271718  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

**ATGGCGTGTGCGTGTGCGCTTGGCGTGTGCTCCTGGCAGTTACCTTCGCAGGTAGCGCTCGGAGCG  
GTCCTGGCGAGCGGGGACCTCCGGAGAAAAGCGGGTTGGGAGTCAGACCGGCGCGGACCTGCCCTGC  
TCCGGCGGCTCGGCGACGGTACCCGCGCGCCGTTACTGGCGGTTCCCAGAGGACCTGCCAGCCTCT  
CTGATCCCTCAATTTGGTCTGTTTTCAAATATAGAACGCCAACTGCTCTCAGTATAGATTACCAGGAT  
GTCCCAGACACTTAACCCTGTGTGTGGCAGTGACATGTCCACTTATGCCAATGAATGTACTCTGTGCAT  
GAAATCAGGGAAGGTGGTCATAATATTAATATTCGAAATGGACCCTGC**

**ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA**

**Protein Sequence:** >RG233290 representing NM\_001271718  
**Red**=Cloning site **Green**=Tags(s)  
MALSVLRLALLLLAVTFAGSARSGPGERPPEKSGFGSQTGGGPCPAPGGLGDGTRAPVTGGSPEDLPAS  
LIPQFGLFSKYRTPNCSQYRLPGCPRHFNPCVCGSDMSTYANECTLCMKIREGHNKIIRNGPC

**TRTRPLE - GFP Tag - V**

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001271718

**ORF Size:** 402 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


**RefSeq:** [NM\\_001271718.2](#)

**RefSeq Size:** 786 bp

**RefSeq ORF:** 405 bp

**Locus ID:** 6691  
**Cytogenetics:** 4q12  
**Protein Families:** Secreted Protein, Transmembrane  
**Gene Summary:** This gene encodes a member of the family of serine protease inhibitors of the Kazal type (SPINK). The encoded protein acts as a trypsin and acrosin inhibitor in the genital tract and is localized in the spermatozoa. The protein has been associated with the progression of lymphomas. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2012]

**Product images:**



Circular map for RG233290