

Product datasheet for **RG233295**

PSMC3IP (NM_001256015) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PSMC3IP (NM_001256015) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: PSMC3IP
Synonyms: GT198; HOP2; HUMGT198A; ODG3; TBPIP
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG233295 representing NM_001256015
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTGAGTGATGCTGACCTCAAGTCCTAGATGGCAAATCGTGGCCCTCACTGCTAAGGTGCAGAGCT
 TGCAGCAGAGCTGCCGTACATGGAGGCTGAGCTCAAGGAATTATCTAGTGCCCTGACCACACCAGAGAT
 GCAGAAAGAAATCCAGGAGTTAAAGAAGGAATGCGCTGGCTACAGAGAGAGATTGAAGAACATTAAGCA
 GCTACCAATCATGTGACTCCAGAAGAGAAAGAGCAGGTGTACAGAGAGAGGCAGAAGTACTGTAAGGAGT
 GGAGGAAGAGGAAGAGGATGGCTACAGAGCTGTCTGATGCAATACTTGAAGGATACCCCAAGAGCAAGAA
 GCAGTTCTTTGAGGAAGTTGGGATAGAGACGGATGAAGATTACAACGTCACACTCCAGACCCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG233295 representing NM_001256015
 Red=Cloning site Green=Tags(s)
 MVSADLQVLDGKIVALAKVQSLQQSCRYMEAEKELSSALTPEMQKEIQELKKECAGYRERLKNIKA
 ATNHVTPEEKEQVYRERQKYCKEWRKRKRMATELSDAILEGYPKSKKQFFEEVGIETDEDYNTLPDP

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-Mlul



[View online »](#)

Cloning Scheme:



ACCN: NM_001256015

ORF Size: 414 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001256015.1](#), [NP_001242944.1](#)

RefSeq Size: 1609 bp

RefSeq ORF: 417 bp

Locus ID: 29893


UniProt ID: [Q9P2W1](#)

Cytogenetics: 17q21.2

Protein Families: Druggable Genome

Gene Summary: This gene encodes a protein that functions in meiotic recombination. It is a subunit of the PSMC3IP/MND1 complex, which interacts with PSMC3/TBP1 to stimulate DMC1- and RAD51-mediated strand exchange during meiosis. The protein encoded by this gene can also co-activate ligand-driven transcription mediated by estrogen, androgen, glucocorticoid, progesterone, and thyroid nuclear receptors. Mutations in this gene cause XX female gonadal dysgenesis. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Dec 2011]

Product images:



Circular map for RG233295