

Product datasheet for **RG233390**

B7H4 (VTCN1) (NM_001253849) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: B7H4 (VTCN1) (NM_001253849) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: B7H4
Synonyms: B7-H4; B7h.5; B7H4; B7S1; B7X; PRO1291; VCTN1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG233390 representing NM_001253849
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTT**CAGAGCCGGACAGCAGT**GTTT**GCTGATCAAGT**GATAGTTGGCAATGCCTCTTTGCGGCTGAAAA
 ACGTGCAACTCACAGATGCTGGCACCTACAAATGTTATATCATCACTTCTAAAGGCAAGGGGAATGCTAA
 CCTTGAGTATAAACTGGAGCCTCAGCATGCCGGAAGTGAATGTGGACTATAATGCCAGCTCAGAGACC
 TTGCGGTGTGAGGCTCCCGATGGTCCCGCAGCCACAGTGGTCTGGGCATCCCAAGTTGACCAGGGAG
 CCAACTTCTCGGAAGTCTCCAATACCAGCTTTGAGCTGAACTCTGAGAATGTGACCATGAAGTTGTGTC
 TGTGCTCTACAATGTTACGATCAACAACATACTCCTGTATGATTGAAAATGACATTGCCAAAGCAACA
 GGGGATATCAAAGTGACAGAATCGGAGATCAAAAGCGGAGTACCTACAGCTGCTAAACTCAAAGGCTT
 CTCTGTGTGCTCTTCTTTGCCATCAGCTGGGCACCTTCTGCCTCTCAGCCCTACCTGATGCTAAA
 A

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG233390 representing NM_001253849
 Red=Cloning site Green=Tags(s)

MFRGRTAVFADQVIVGNASLRLKNVQLTDAGTYKCYIITSKGKGNANLEYKTGAFSMPEVNVNDYNASSET
 LRCEAPRWFPPQTVVWASQVDQGANFSEVSNTSFELNSENVTMKVSVLYNVTINNTYSCEIENDIAKAT
 GDKVTESEIKRRSHLQLLNKASLKVSSFFAISWALLPLSPYMLK

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001253849

ORF Size: 561 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001253849.1](#), [NP_001240778.1](#)

RefSeq Size: 2734 bp

RefSeq ORF: 564 bp

Locus ID: 79679


UniProt ID: [Q7Z7D3](#)

Cytogenetics: 1p13.1-p12

Protein Families: Transmembrane

Gene Summary: This gene encodes a protein belonging to the B7 costimulatory protein family. Proteins in this family are present on the surface of antigen-presenting cells and interact with ligand bound to receptors on the surface of T cells. Studies have shown that high levels of the encoded protein has been correlated with tumor progression. A pseudogene of this gene is located on chromosome 20. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011]

Product images:



Circular map for RG233390