

## Product datasheet for **RG233437**

### MLC1SA (MYL6B) (NM\_001199629) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** MLC1SA (MYL6B) (NM\_001199629) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** MYL6B  
**Synonyms:** MLC1SA  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG233437 representing NM\_001199629  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCTCCCAAGAAGGATGTTCCCGTGAAGAAACCAGCAGGGCCCTCCATCTCCAAACCTGCTGCTAAGC  
CAGCAGCAGCAGGGGCTCCTCCAGCCAAGACAAAGCTGAGCCAGCTGTCCCCAGGCCCTCAGAAAAC  
CCAGGAGCCTCCAGTCGATCTCTCAAAGTGGTGATCGAGTTTAAACAAGGACCAGCTGGAGGAGTTCAAG  
GAGGCCTTCGAGCTGTTTGACCGAGTGGGGATGGCAAGATCCTGTACAGCCAGTGTGGGACGTGATGA  
GGCCCTGGGCCAGAACCCACCAACGCCGAGGTGCTCAAGGTCCTGGGGAACCCAAGAGTGATGAGCT  
GAAGTCGCGGCGTGTGGACTTTGAGACTTTCCTGCCATGCTCCAGGCAGTGGCCAAGAACCAGGCCAA  
GGCACATATGAGGACTACTTGGAGGGGTTTCGTGTGTTTGACAAGGAGGGGAACGGCAAAGTCATGGGAG  
CAGAGCTCAGACATGTTCTCACCACCCTTGAGAGAAGATGACTGAGGAGGAGGTGGAGACCGTTCTGGC  
AGGACACGAGGACAGCAACGGCTGCATCAACTACGAGGCCTTCTTGAACACATCCTAAGCGTC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG233437 representing NM\_001199629  
Red=Cloning site Green=Tags(s)

MPPKDVVVKKPAAGPSISKPAAKPAAAGAPPAKTKAEPVAVPQAPQKTQEPVVDLSKVVIEFNKDQLEEFK  
EAFELFDRVGDGKILYSQCGDVMRALGQNPTNAEVLKVLGNPKSDELKSRVDFETFLPMLQAVAKNRGQ  
GTYEDYLEGFRVFDKEGNGKVMGAELRHVLTTLGKMTTEEEVETVLAGHEDSNGCINYEAFLLKHLVSV


**TRTRPLE** - GFP Tag - V

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**

**Plasmid Map:**


**ACCN:** NM\_001199629

**ORF Size:** 624 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001199629.2](#)

**RefSeq Size:** 1072 bp

**RefSeq ORF:** 627 bp

**Locus ID:** 140465

**UniProt ID:** [P14649](#)

**Cytogenetics:** 12q13.2

**Protein Pathways:** Vascular smooth muscle contraction

**Gene Summary:** Myosin is a hexameric ATPase cellular motor protein. It is composed of two heavy chains, two nonphosphorylatable alkali light chains, and two phosphorylatable regulatory light chains. This gene encodes a myosin alkali light chain expressed in both slow-twitch skeletal muscle and in nonmuscle tissue. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2010]