


Product datasheet for RG233455

EMA (MUC1) (NM_001204296) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EMA (MUC1) (NM_001204296) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	EMA
Synonyms:	ADMCKD; ADMCKD1; ADTKD2; CA 15-3; CD227; EMA; H23AG; KL-6; MAM6; MCD; MCKD; MCKD1; MUC-1; MUC-1/SEC; MUC-1/X; MUC1/ZD; PEM; PEMT; PUM
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	<pre>>RG233455 representing NM_001204296 Red=Cloning site Blue=ORF Green=Tags(s)</pre>
	TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCC <mark>GCGATCGC</mark> C
	ATGACACCGGGCACCCAGTCTCCTTTCTTCCTGCTGCTGCTCCTCACAGTGCTTACAGCTACCACAGCCC CTAAACCCGCAACAGTTGTTACGGGTTCTGGTCATGCAAGCTCTACCCCAGGTGGAGAAAAGGAGACTTC
	GGCTACCCAGAGAAGTTCAGTGCCCAGCTCTACTGAGAAGAATGCTATCCCAGCACCGACTACTACCAAG
	AGCTGCAGAGAGACATTTCTGAAATGTTTTTGCAGATTTATAAACAAGGGGGTTTTCTGGGCCTCTCCAA
	TATTAAGTTCAGTGAGTGATGTGCCATTTCCTTTCTCTGCCCAGTCTGGGGCTGGGGGTGCCAGGCTGGGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US Phone: +1-888-267-4436

https://www.origene.com techsupport@origene.com

	EMA (MUC1) (NM_001204296) Human Tagged ORF Clone – RG233455
Protein Sequence:	: >RG233455 representing NM_001204296 Red=Cloning site Green=Tags(s)
	MTPGTQSPFFLLLLLTVLTATTAPKPATVVTGSGHASSTPGGEKETSATQRSSVPSSTEKNAIPAPTTTK SCRETFLKCFCRFINKGVFWASPILSSVSDVPFPFSAQSGAGVPGWGIALLVLVCVLVALAIVYLIALAV CQCRRKNYGQLDIFPARDTYHPMSEYPTYHTHGRYVPPSSTDRSPYEKVSAGNGGSSLSYTNPAVAATSA NL
	TRTRPLE - GFP Tag - V
Restriction Sites:	Sgfl-Mlul
Cloning Scheme:	Cloning sites used for ORF Shuttling:
	Sgfi ORF Miui GCGATOGC C ATG SXX ACG CGT
	EcoR I BamH I Kpn I RBS Sgf I Asc I CTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGGAGGAGATCTGCCGCCGCGGCGCGCGC
	Hind III Nhe I Rsr II MIU I Not I Xho I GFP Tag CAAGCTTAACTAGCTAGCGGACCG ACG CGT ACG CGG CCG CTC GAG ATG GAG AGC GAC T R T R P L E M E S D
	GAA GAA AGA GTT TAA ACGGCCGGCCGCGGAGCT E E R V Stop

ACCN:	NM_001204296
ORF Size:	636 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Service Servic

Reconstitution Method:	 Centrifuge at 5,000xg for 5min. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. Close the tube and incubate for 10 minutes at room temperature. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM 001204296.2</u>
RefSeq Size:	1037 bp
RefSeq ORF:	639 bp
Locus ID:	4582
UniProt ID:	<u>P15941</u>
Cytogenetics:	1q22
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
Gene Summary:	This gene encodes a membrane-bound protein that is a member of the mucin family. Mucins are O-glycosylated proteins that play an essential role in forming protective mucous barriers on epithelial surfaces. These proteins also play a role in intracellular signaling. This protein is expressed on the apical surface of epithelial cells that line the mucosal surfaces of many different tissues including lung, breast stomach and pancreas. This protein is proteolytically cleaved into alpha and beta subunits that form a heterodimeric complex. The N-terminal alpha subunit functions in cell-adhesion and the C-terminal beta subunit is involved in cell signaling. Overexpression, aberrant intracellular localization, and changes in glycosylation of this protein have been associated with carcinomas. This gene is known to contain a highly polymorphic variable number tandem repeats (VNTR) domain. Alternate splicing results in

multiple transcript variants.[provided by RefSeq, Feb 2011]

Product images:



Circular map for RG233455

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US