

## Product datasheet for **RG233469**

### IST1 (NM\_001270978) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	IST1 (NM_001270978) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	IST1
Synonyms:	CHMP8; OLC1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG233469 representing NM_001270978 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**

ATGCACAAGCTGAGTGTGGAAGCCCCACCCAAAATCCTGGTGGAGAGATACCTGATTGAAATTGCAAAGA  
ATTACAACGTACCCTATGAACCTGACTCTGTGGTCATGGCAGAAGCTCCTCCTGGGGTAGAGACAGATCT  
TATTGATGTTGGATTCACAGATGATGTGAAGAAAGGAGGCCCTGGAAGAGGAGGGAGTGGTGGCTTCACA  
GCACCAGTTGGTGGACCTGATGGAACGGTGCCAAATGCCCATGCCATGCCATGCCTATGCCATCTGCAA  
ATACGCCCTTCTCATATCCACTGCCAAAGGACCATCAGATTTCAATGGACTGCCAATGGGGACTTATCA  
GGCCTTCCCAAATATTCATCCACCTCAGATACCAGCAACTCCCCATCGTATGAATCTGTAGATGACATT  
AATGCTGATAAGAATATCTCTTCTGCACAGATTGTTGGTCTGGACCCAAGCCAGAAGCCTCTGCAAAGC  
TTCTTCCAGACCTGCAGATAACTATGACAACCTTGTCTACCAGAGTTGCCATCTGTGCCAGACACACT  
ACCAACTGCATCTGCTGGTGCCAGCACCTCAGCATCTGAAGACATTGACTTTGATGATCTTCCCGGAGG  
TTTGAAGAGCTGAAAAAGAAAACA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG233469 representing NM\_001270978  
 Red=Cloning site Green=Tags(s)

MHKL SVEAPPKILVERYLIEIAKNYNVPYEPDSVVM AEAPPGVETDLIDVGF TDDVKKGGPGRGSGGFT  
 APVGGPDGTVPMPPMPMPMP SANTSPPFSYPLPKGPSDFNGLPMGTYQAFPNIHPPQIPATPPSYESVDDI  
 NADKNISSAQIVGPGKPEASAKLPSRPADNYDNFVLP ELPSPVDTLPTASAGASTASEIDFDL SRR  
 FEELKKKT

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001270978

**ORF Size:** 654 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001270978.2](#)

**RefSeq Size:** 3667 bp

**RefSeq ORF:** 657 bp


**Locus ID:** 9798

**UniProt ID:** [P53990](#)

**Cytogenetics:** 16q22.2

**Gene Summary:** This gene encodes a protein with MIT-interacting motifs that interacts with components of endosomal sorting complexes required for transport (ESCRT). ESCRT functions in vesicle budding, such as that which occurs during membrane abscission in cytokinesis. There is a pseudogene for this gene on chromosome 19. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Aug 2012]

**Product images:**



Circular map for RG233469