

## Product datasheet for **RG233471**

### ANTKMT (NM\_001271285) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ANTKMT (NM\_001271285) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** ANTKMT  
**Synonyms:** ANT-KMT; C16orf24; FAM173A  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG233471 representing NM\_001271285  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAGCAGGACGACCCGGTCGAGGCGCTGACGGAGCTGCGGAGCGCGCGGCTGGGCGCGCTGGAGCTGC  
TGCAGGCGGGCGCCGGCTCGGGCTTGGCAGCCTACGCGGTGTGGGCGCTGCTGCTCCAGCCCGGCTCCG  
GCGCGTCCGCTGCGGCTGCAGGTGCCCTACGTCGGCGGAGCGCGCGGAGGTGGAGCACGTGTTGTCG  
CTGCTGCGAGGACGCCCCGAAAAACGGTGGATCTGGGCTCTGGCGACGGCAGGATCGTGTGGCGGCC  
ACAGGTGCGGCCCTCCGCCGGCCGTGGGCTACGAGCTGAACCCCTGGCTGGTGGCGCTGGCGCGGCTGCA  
CGCTGGAGGGCCGGCTGTGCCGCGAGCTCTGCTATCGCCGCAAGGATCTCTGGAAGCTCCCGCTGCTG  
GAGGACAAGCTGCGGACAGAGCTGCCTGCTGGGCCCCGCGTGGTGTCTGGGCGCTTCCCACTCCCCACCT  
GGCAGCCTGTGACCGCGGTTGGCGAGGGCCTGGACCGAGTATGGGCTTATGATGTTCTGAGGGTGGGCA  
GGCTGGGAGGCCGCTCCTCGCGGATACCCATCCAGGCTGCCCGGACCTAGTTCTGCCCCATCCCG  
GGGGCCCTATTCTCAGGCCAGC

**ACGCGTACGCGGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG233471 representing NM\_001271285  
 Red=Cloning site Green=Tags(s)

MEQDDPVEALTELRERRLGALELLQAAAGSGLAAYAVWALLLQPGFRRVPLRLQVPYVGASARQVEHVLS  
 LLRGRPGKTVDLGSGDGRIVLAAHRCGLRPAVGYELNPWLVALARLHAWRAGCAGSVCYRRKDLWKLPLL  
 EDKLRTELPAGARVVSGRFPLPTWQPVTAVGEGLDREVWAYDVPEGGQAGEAASSRIPIQAAPGPSSAPIP  
 GGLISQAS

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001271285

**ORF Size:** 654 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001271285.2](#)

**RefSeq Size:** 830 bp

**RefSeq ORF:** 657 bp


**Locus ID:** 65990

**Cytogenetics:** 16p13.3

**Protein Families:** Transmembrane

**Gene Summary:** S-adenosyl-L-methionin-dependent protein-lysine N-methyltransferase.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for RG233471