

## Product datasheet for **RG233543**

### COLEC11 (NM\_001255982) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** COLEC11 (NM\_001255982) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** COLEC11  
**Synonyms:** 3MC2; CL-11; CL-K1-I; CL-K1-II; CL-K1-IIa; CL-K1-IIb; CLK1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG233543 representing NM\_001255982  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGGGGAATCTGGCCCTGGTGGCGTTCTAATCAGCCTGGCCTTCCTGTCAGTCTGCCATCTGGAC  
ATCCTCAGCCGGCTGGCGATGACGCCTGCTGTGCAGATCCTCGTCCCTGGCCTCAAAGGAGACATGGG  
GGACAAAGGACAGAAAGGCAGTGTGGGTCGTATGGAAAATTGGTCCCATTGGCTCTAAAGGTGAGAAA  
GGAGATCCGGTGACATAGGACCCCTGGTCTAATGGAGAACCAGGCCTCCCATGTGAGTGCAGCCAGC  
TGGCAAGGCCATCGGGGAGATGGACAACCAGGTCTCTCAGCTGACCAGCGAGCTCAAGTTCATCAAGAA  
TGATGTCGCCGGTGTGCGCGAGACGGAGCAAGATCTACCTGCTGGTGAAGGAGGAGAAGCGCTACGCG  
GACGCCAGCTGTCTGCCAGGGCCGCGGGGACGCTGAGCATGCCAAGGACGAGGCTGCCAATGGCC  
TGATGGCCGATACCTGGCGCAAGCCGGCTGGCCCGTCTTTCATCGGCATCAACGACCTGGAGAAGGA  
GGGCGCCTTCGTGTAATCTGACCACTCCCCATGCGGACCTTCAACAAGTGGCGCAGCGGTGAGCCCAAC  
AATGCCTACGACGAGGAGGACTGCGTGGAGATGGTGGCCTCGGGCGGCTGGAACGACGTGGCCTGCCACA  
CCACCATGTAATCATGTGTGAGTTTGACAAGGAGAACATG

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG233543 representing NM\_001255982  
Red=Cloning site Green=Tags(s)

M R G N L A L V G V L I S L A F L S L L P S G H P Q P A G D D A C S V Q I L V P G L K G D M G D K G Q K G S V G R H G K I G P I G S K G E K  
 G D S G D I G P P G P N G E P L P C E C S Q L R K A I G E M D N Q V S Q L T S E L K F I K N A V A G V R E T E S K I Y L L V K E E K R Y A  
 D A Q L S C Q G R G G T L S M P K D E A A N G L M A A Y L A Q A G L A R V F I G I N D L E K E G A F V Y S D H S P M R T F N K W R S G E P N  
 N A Y D E E D C V E M V A S G G W N D V A C H T T M Y F M C E F D K E N M

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001255982

**ORF Size:** 741 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001255982.2](#)

**RefSeq Size:** 1633 bp

**RefSeq ORF:** 744 bp


**Locus ID:** 78989

**UniProt ID:** [Q9BWP8](#)

**Cytogenetics:** 2p25.3

**Gene Summary:** This gene encodes a member of the collectin family of C-type lectins that possess collagen-like sequences and carbohydrate recognition domains. Collectins are secreted proteins that play important roles in the innate immune system by binding to carbohydrate antigens on microorganisms, facilitating their recognition and removal. The encoded protein binds to multiple sugars with a preference for fucose and mannose. Mutations in this gene are a cause of 3MC syndrome-2. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2011]

**Product images:**



Circular map for RG233543