

## Product datasheet for **RG233637**

### CREM (NM\_001267564) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CREM (NM_001267564) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CREM
Synonyms:	CREM-2; hCREM-2; ICER
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG233637 representing NM_001267564 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTGGTGGCATCAGCATAATCTATGTTTCAGGCGTCCTATAGAAGAGGATTATTCTTCAGGGGATGTGG  
AAGAAAAGGTAGCAGCAATTGCAGAGACAGATGAATCTGCAGAATCAGAAGGTGAATTGATTCTCATAA  
ACGTAGAGAAATCCTTTACGAAGACCCTTTATAGGAAAATACTGAATGAAGTGCCTCTGATGTGCCT  
GGTGTTCCTCAAGATTGAAGAAGAGAGATCAGAGGAAGAAGGAACACCACCTAGTATTGCTACCATGGCAG  
TACCAACTAGCATATATCAGACTAGCACGGGCAATACATTGCTATAGCCCAAGGTGGAACAATCCAGAT  
TTCTAACCAGGATCTGATGGTGTTCAGGGACTGCAGGCATTAACAATGACAAATTCAGGAGCTCCCTCCA  
CCAGGTGCTACAATTGTACAGTACGCAGCACAATCAGCTGATGGCACACAGCAGTTCTTTGTCCCAGGCA  
GCCAGGTGTGTGTTCAAGATGAGGAAACTGAACTTGCCCAAGTCACATGGCTGCTGCCACTGGTGACAT  
GCCAECTTACCAGATCCGAGCTCCTACTGCTGCTTTGCCACAGGGAGTGGTGTGCTGCTGCATCGCCCGGA  
AGTTTGCACAGTCCCAGCAGCTGGCAGAAGAAGCAACACGCAACGAGAGCTGAGGCTAATGAAAAACA  
GGGAAGCTGCCAAGAATGTCGACGTCGAAAGAAAAGAAATGTAAAATGTCTGGAGAGCCGAGTTGCAGT  
GCTGGAAGTCCAGAACAAGAAGCTTATAGAGGAACTTGAACCTTGAAAGACATTTGTTCTCCAAAAACA  
GATTAC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG233637 representing NM\_001267564  
Red=Cloning site Green=Tags(s)

MWWHQHNLCFRRPIEEDYSSGDVEEKVAAIAETDESAESEGVIDSHKRREILSRRPSYRKILNELSSDVP  
 GVPKIEEERSEEEGTPPSIATMAVPTSIYQTSTGQYIAIAQGGTIQISNPGSDGVQGLQALMTNSGAPP  
 PGATIVQYAAQSADGTQQFFVPGSQVVVQDEETELAPSHMAAATGDMPTYQIRAPTAALPQGVVMAASPG  
 SLHSPQQLAEEATRRELRLMKNREAAKECRRRKKEYVKCLESRVAVLEVQNKKLEEELETLKDICSPKT  
 DY

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001267564

**ORF Size:** 846 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001267564.2](#)

**RefSeq Size:** 2196 bp

**RefSeq ORF:** 849 bp

**Locus ID:** 1390


**UniProt ID:** [Q03060](#)

**Cytogenetics:** 10p11.21

**Protein Families:** Druggable Genome, Transcription Factors

**Gene Summary:** This gene encodes a bZIP transcription factor that binds to the cAMP responsive element found in many viral and cellular promoters. It is an important component of cAMP-mediated signal transduction during the spermatogenic cycle, as well as other complex processes. Alternative promoter and translation initiation site usage allows this gene to exert spatial and temporal specificity to cAMP responsiveness. Multiple alternatively spliced transcript variants encoding several different isoforms have been found for this gene, with some of them functioning as activators and some as repressors of transcription. [provided by RefSeq, Jul 2008]

**Product images:**



Circular map for RG233637