

Product datasheet for **RG233879**

C20orf96 (NM_080571) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C20orf96 (NM_080571) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	C20orf96
Synonyms:	dj1103G7.2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG233879 representing NM_080571 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGCTTACAAAAGTAGACCAAGCACTCTGGGACTCACTCCATAGTCCAGGAGTTCAGGTTCCGGATT
ATGTTCCATGGCAGCAGTCCAAGCAGGAAACCAAGCCATCTACTCTGCCTCCAGTCCAACAAGCCAACAG
CCTTCATACAAGCAAAATGAAGACTTTGACTAGGGTCCAACCAAGTGTTCCTCAAGCCCACTACGGTG
GTGACAAGCTGCCAGCCGAAGAATCCAAGAGAACTACATAGAAGGCGGAAGTTGGACCCTGGGAAGATGC
ATGCCAAATCTGGTTAATGAAGACCTCGCTCAGGAGCGGGAGGCCGCTCTGCGAGAGCTCCGAAGCCG
TGAGAACTTCCTCAGCAAGCTCAACCGGAGCTGATCGAGACCATCCAGGAGATGGAGAACAGCAGCAGC
CTGCACGTGCGGGCCCTGCTGCAGCAGCAGGACACCCCTGGCGACCATCATCGACATCTTGGAGTACTCAA
ACAAGAAGAGGCTGCAGCAATTGAAATCTGAGCTTCAGGAGTGGGAAGAAAAGAAATGCAAGATGAG
CTATCTTGAGCAGCAGGCAGAGCAGCTGAATGCCAAGATTGAGAAGACCCAGGAGGAAGTGAAGTTCCTG
AGCACTTACATGGACCATGAGTATTCATCAAGTCTGTCCAGATCTCACTCTTATGCGCCAGCTGCAGC
AGGTTAAGGACAGCCAGCAGGATGAGCTGGATGACCTCGGTGAGATGCGCAGAAAGGTCTTGAATCCTT
GTCTGACAAGATTCAGAAGAAGAAGAAAAAATCTGAGTCTGTGGTGGCGGAAACCCAGCGTCCCTAT
GAAGAGGCTCTCTACAGAAGATGTGGAAAGCCAGGACTTCTGAAATGCATGCAAGGTTTCAGAGAAA
TTATTGACCAGTTTGAGGAGAACATGCCTGTATTAAGGGCCGAGGTGGAAGAGCTCCAAGCCAGACCCG
GGAACCCCGAGAGGTATATTTGAGGATGTTCTGCTTCGGAGACCAAGTGCACCCAGACATGGATGTC
ATCCTCAACATTCTGTGAAGAGCCACTACCTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG233879 representing NM_080571
Red=Cloning site Green=Tags(s)


MSYKSRPKHSGTHSIVQEFQVPDYVPWQQSKQETKPSLPPVQQANSLHTSKMKTLTRVQPVFHFKPTTV
VTSCQPKNPRELHRRRKLDPGKMHAKIWLTKSLRSGRAALRELSRENFLSKLNRELIETIQEMENSTT
LHVRALLQQQDTLATIIDILEYSNKKRLQQLKSELQEWEEKKKCKMSYLEQQAEQLNAKIEKTQEEVNFL
STYMDHEYSIKSVQISTLMRQLQQVKDSQQDELDDL GEMRRKVLESLSDKI QKKKKKILSSVVAETQRPY
EEALLQKMWESQDFLKCMQRFREIIDQFEENMPVLR AEVEELQAQTREPREVIFEDVLLRRPKCTPDMDV
ILNIPVEEPLPF

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_080571

ORF Size: 1086 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_080571.1](#), [NP_542138.1](#)

RefSeq Size: 1615 bp


RefSeq ORF: 1089 bp

Locus ID: 140680

UniProt ID: [Q9NUD7](#)

Cytogenetics: 20p13

Product images:



Circular map for RG233879