

Product datasheet for **RG233919**

SERPINB6 (NM_001271824) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	SERPINB6
Synonyms:	CAP; DFNB91; MSTP057; PI-6; PI6; PTI; SPI3
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG233919 representing NM_001271824
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGATGTTCTCGCAGAAGCAAATGGCACCTTTGCCTTAAACCTTTTGAAAACGCTGGGTAAAGACAAC
CGAAGAATGTGTTTTCTCACCCATGAGCATGTCCTGTGCCCTGGCCATGGTCTACATGGGGGCAAGGG
AAACACCGCTGCACAGATGGCCAGATACTTTCTTCAATAAAAGTGCCGGTGGTGGAGACATCCACCAG
GGCTTCCAGTCTCTCTCACCGAAGTGAACAAGACTGGCACGCAGTACTTGCTTAGGATGGCCAACAGGC
TCTTTGGGGAAAAGTCTTGATTTCTCTCATCTTTTAGAGATTCTGCCAAAAATTCTACCAAGCAGA
GATGGAGGAGCTTGACTTTATCAGCGCCGTAGAGAAGTCCAGAAAACACATAAACACCTGGGTAGCTGAA
AAGACAGAAGGTAAAATTGCGGAGTTGCTCTCTCCGGGCTCAGTGGATCCATTGACAAGGCTGGTCTGG
TGAATGCTGTCTATTTAGAGGAACTGGGATGAACAGTTTGACAAGGAGAACACCGAGGAGAGACTGTT
TAAAGTCAGCAAGAATGAGGAGAACTGTGCAATGATGTTTAAGCAATCTACTTTTAAGAAGACCTAT
ATAGGAGAAATATTTACCCAAATCTTGGTGCTTCCATATGTTGGCAAGGAACGAATATGATCATCATGC
TTCCGGACGAGACCACTGACTTGAGAACGGTGGAGAAAGAACTCACTTACGAGAAGTTCGTAGAATGGAC
GAGGCTGGACATGATGGATGAAGAGGAGGTGGAAGTGCCCTCCCGCGGTTTAACTAGAGGAAAGCTAC
GACATGGAGAGTGTCTGCGCAACCTGGGCATGACTGATGCCTTCGAGCTGGGCAAGGCAGACTTCTCTG
GAATGTCCCAGACAGACCTGTCTGTCCAAGGTCGTGCACAAGTCTTTTGTGGAGGTCAATGAGGAAGG
CACGGAGGCTGCAGCCGCCACAGCTGCCATCATGATGCGGTGTGCCAGATTCTCCCCCGCTTCTGC
GCCGACCCCTTCTTTTCTTCATCCAGCACAGCAAGACCAACGGGATTCTTCTGCGGCCGCTTTT
CCTCTCCG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence:

>RG233919 representing NM_001271824

Red=Cloning site Green=Tags(s)


MDVLAENGTFALNLLKTLGKDNSKNVFFSPMSMSCALAMVYMGAKGNTAAQMAQILSFNKGSGGGDIHQ
GFQSLLTEVNKTGTQYLLRMANRLFGEKSCDFLSSFRDSCQKFYQAEMEELDFISAVEKSRKHINTWVAE
KTEGKIAELLSPGSVDPLTRLVLVNAVYFRGNWDEQFDKENTEERLFKVSKEEKPQMMFKQSTFKKTY
IGEIFTQILVLPYVGKELNMIIMLPDETDLRTVEKELTYEKFVWTRLDMMDEEEVSLPRFKLEESY
DMESVLRNLGMTDAFELGKADFSGMSQTDLSLKVHKSFEVNEEGTEAAAATAAIMMRCARFVPRFC
ADHPFLFFIQHKTNGILFCGRFSSP

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001271824

ORF Size:

1128 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:


This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001271824.2</u>
RefSeq Size:	1972 bp
RefSeq ORF:	1131 bp
Locus ID:	5269
UniProt ID:	<u>P35237</u>
Cytogenetics:	6p25.2
Protein Families:	Druggable Genome
Gene Summary:	The protein encoded by this gene is a member of the serpin (serine proteinase inhibitor) superfamily, and ovalbumin(ov)-serpin subfamily. It was originally discovered as a placental thrombin inhibitor. The mouse homolog was found to be expressed in the hair cells of the inner ear. Mutations in this gene are associated with nonsyndromic progressive hearing loss, suggesting that this serpin plays an important role in the inner ear in the protection against leakage of lysosomal content during stress, and that loss of this protection results in cell death and sensorineural hearing loss. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Sep 2010]

Product images:



Circular map for RG233919