

Product datasheet for **RG233948**

Corticotropin Releasing Factor Receptor 2 (CRHR2) (NM_001202483) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Corticotropin Releasing Factor Receptor 2 (CRHR2) (NM_001202483) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Corticotropin Releasing Factor Receptor 2
Synonyms:	CRF-RB; CRF2; CRFR2; HM-CRF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG233948 representing NM_001202483
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGACGCGGCACTGCTCCACAGCCTGCTGGAGGCCAACTGCAGCCTGGCGCTGGCTGAAGAGCTGCTCT
 TGGACGGCTGGGGCCACCCCTGGACCCGAGGGTCCCTACTCCTACTGCAACACGACCTTGGACCAGAT
 CGGAACGTGCTGGCCCCGAGCGCTGCCGAGCCCTCGTGGAGAGGCCGTGCCCGAGTACTTCAACGGC
 GTCAAGTACAACACGACCCGGAATGCCTATCGAGAATGCTTGGAGAATGGGACGTGGGCCTCAAAGATCA
 ACTACTCACAGTGTGAGCCATTTGGATGACAAGCAGAGGAAGTATGACCTGCACTACCGCATCGCCCT
 TGTCTGCAACTACCTGGGCCACTGCGTATCTGTGGCAGCCCTGGTGGCCGCTTCTGCTTTTCTGGCC
 CTGCGGAGCATTGCTGTCTGCGGAATGTATTACTGGAACCTCATCACCACCTTTATCCTGCGAAATG
 TCATGTGGTTCTGCTGCAGCTCGTTGACCATGAAGTGCACGAGAGCAATGAGGTCTGGTCCGCTGCAT
 CACCACCATCTCAACTACTTCTGTTGACCAACTTCTTCTGGATGTTTGTGGAAGGCTGCTACCTGCAC
 ACGGCCATTGTCATGACCTACTCCACTGAGCGCCTGCGCAAGTGCCTTCTCTTTCATCGGATGGTGCA
 TCCCTTCCCCATCATCGTCGCCTGGGCCATCGGCAAGCTCTACTATGAGAATGAACAGTCTGGTTTGG
 CAAGGAGCCTGGCGACCTGGTGGACTACATCTACCAAGGCCCATCATTCTCGTGCTCCTGATCAATTTT
 GTATTTCTGTTAACATCGTCAGGATCCTAATGACAAAGTTACGCGCGTCCACCACATCCGAGACAATCC
 AGTACAGGAAGGCAGTGAAGGCCACCCTGGTGTCTCTGCCCTCCTGGGCATCACCTACATGCTCTTCTT
 CGTCAATCCGGGGAGGACGACCTGTACAGATCATGTTTCATCTATTTCAACTCCTTCTGCAGTCGTTT
 CAGGGTTTCTTCGTGTCTGCTTCTACTGCTTCTTCAATGGAGAGAGCTGGGTGAGCAAGGAAGCCAGG
 CCGCAGGGCCTCATGGAAGAGAGAAGCCAGAGCAGAGGTGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG233948 representing NM_001202483
 Red=Cloning site Green=Tags(s)

MDAALLHSLLEANCSLALAEELLLDGGWPPLDPEGPYSYCNTTLDQIGTCWPRSAAGALVERPCPEYFNG
 VKYNTTRNAYRELENGTWASKINYSQCEPILDDKQRKYDLHYRIALVVNYLGHCVSVAALVAAFLLFLA
 LRSIRCLRNVIHWNLIITFILRNVMWFLQLVDHEVHESNEVWCRCITTFNYFVVTNFFWFMFVEGCVLH
 TAIVMTYSTERLRKCLFLFIGWCIPFPIIIVAWAIGKLYENEQCWFKEPGDLVDYIYQGP IILVLLINF
 VFLFNIVRILMTKLRASSTSETIQYRKAVKATLVLLPLLGITIYMLFFVNPGEDDLSQIMFIYFNSFLQSF
 QGFFVSVFYCFNNGESWVSKEAQAAGPHGREKPEQRW

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:



ACCN: NM_001202483

ORF Size: 1161 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001202483.1](#), [NP_001189412.1](#)

RefSeq Size: 3109 bp

RefSeq ORF: 1164 bp

Locus ID: 1395

UniProt ID: [Q13324](#)


Cytogenetics: 7p14.3

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Neuroactive ligand-receptor interaction

Gene Summary: The protein encoded by this gene belongs to the G-protein coupled receptor 2 family, and the subfamily of corticotropin releasing hormone receptor. This receptor shows high affinity for corticotropin releasing hormone (CRH), and also binds CRH-related peptides such as urocortin. CRH is synthesized in the hypothalamus, and plays an important role in coordinating the endocrine, autonomic, and behavioral responses to stress and immune challenge. Studies in mice suggest that this receptor maybe involved in mediating cardiovascular homeostasis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.[provided by RefSeq, Jan 2011]

Product images:



Circular map for RG233948