

Product datasheet for **RG233967**

B4GALT3 (NM_001199873) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	B4GALT3 (NM_001199873) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	B4GALT3
Synonyms:	beta4Gal-T3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG233967 representing NM_001199873 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTGCGGAGGCTGCTGGAGCGGCCTTGACAGCTGGCCCTGCTTGTGGGCTCCCAGCTGGCTGTCATGA
TGTACCTGCTACTGGGGGCTTCCGAAGTCTCAGTGCCCTATTTGGCCGAGATCAGGGACCGACATTTGA
CTATTCTACCCCTCGTATGTCTACAGTAACCTCAGTCACCTGCCTGGGGCCCCAGGGGGTCTCCAGCT
CCTCAAGGTCTGCCCTACTGTCCAGAACGATCTCCTCTCTTAGTGGTCTCTGTGCGGTGTCCTTTAGCC
CAGTGCCATCACTGGCAGAGATTGTGGAGCGGAATCCCCGGGTAGAACAGGGGGCCGGTACCGCCCTGC
AGGTTGTGAGCCCGCTCCCGAACAGCCATCATTGTGCCTCATCGTGGCCGGGAGACCACCTGCGCCTG
CTGCTCTACCACCTGCACCCCTTCTTGACGCGCCAGCAGCTTGCTTATGGCATCTATGTCATCCACCAGG
CTGAAATGGAACATTTAACAGGGCAAACTGTTGAACGTTGGGGTGCAGAGGCCCTGCGTGATGAAGA
GTGGGACTGCCTGTTCTTGACGATGTGGACCTTTGCCAGAAAATGACCACAATCTGTATGTGTGTGAC
CCCCGGGACCCCGCATGTTGCCGTTGCTATGAACAAGTTGGATACAGCCTCCCGTACCCCACTACT
TCGGAGGAGTCTCAGCACTTACTCCTGACCAGTACCTGAAGATGAATGGCTTCCCAATGAATACTGGGG
CTGGGGTGGTGAAGATGACGACATTGCTACCAGGTGCGCCTGGCTGGGATGAAGATCTCTCGGCCCCCC
ACATCTGTAGGACACTATAAGATGGTGAAGCAGGAGGATAAGGGCAATGAGGAAAATCCCCACAGAT
TTGACCTCTGGTCCGTACCCAGAATTCCTGGACGCAAGATGGGATGAACCTCACTGACATACCAAGTTGCT
GGCTCGAGAGCTGGGCCTCTTTATACCAACATCACAGCAGACATTGGGACTGACCTCGGGGTCCTCGG
GCTCCTCTGGGCCACGTTACCCACCTGGTTCTCCCAAGCCTTCCGTCAAGAGATGCTGCAACGCCGGC
CCCCAGCCAGGCCTGGCCTCTACTACTGCCAACACACAGCCCTCCGAGGTTACAC

ACCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG233967 representing NM_001199873
 Red=Cloning site Green=Tags(s)

MLRRLLEPCTLALLVGSQAVMMYLSLGGFRSLALFGRDQGPTFDYSHPRDVYNSLSHLPGAPGGPPA
 PQGLPYCPERSPLL VGPVSVSFSVPVSLAEIVERNPRVEPGGRYRPAGCEPRRTAIIVPHRAREHHLRL
 LLYHLHPFLQRQQLAYGIYVIHQAGNGTFNRAKLLNVGVREALRDEEWDCLFLHDVDLLPENDINGLYVCD
 PRGPRHVAVAMNKFGYSLPYPQYFGGVSALTPDQYLMNGFPNEYWGWGGEDDDIATRVLAMKISRPP
 TSVGHYKMKVHRGDKGNEENPHRFDLLVRTQNSWTQDGMNSLT YQLLARELGPLYTNITADIGTDPRGPR
 APSGPRYPPGSSQAFRQEMLQRRPPARPGPLSTANHTALRGSH


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001199873

ORF Size:	1179 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001199873.1 , NP_001186802.1
RefSeq Size:	2438 bp
RefSeq ORF:	1182 bp
Locus ID:	8703
UniProt ID:	O60512
Cytogenetics:	1q23.3
Protein Families:	Transmembrane
Protein Pathways:	Glycosphingolipid biosynthesis - lacto and neolacto series, Keratan sulfate biosynthesis, Metabolic pathways, N-Glycan biosynthesis
Gene Summary:	This gene is one of seven beta-1,4-galactosyltransferase (beta4GalT) genes. They encode type II membrane-bound glycoproteins that appear to have exclusive specificity for the donor substrate UDP-galactose; all transfer galactose in a beta1,4 linkage to similar acceptor sugars: GlcNAc, Glc, and Xyl. Each beta4GalT has a distinct function in the biosynthesis of different glycoconjugates and saccharide structures. As type II membrane proteins, they have an N-terminal hydrophobic signal sequence that directs the protein to the Golgi apparatus and which then remains uncleaved to function as a transmembrane anchor. By sequence similarity, the beta4GalTs form four groups: beta4GalT1 and beta4GalT2, beta4GalT3 and beta4GalT4, beta4GalT5 and beta4GalT6, and beta4GalT7. This gene encodes an enzyme that may be mainly involved in the synthesis of the first N-acetyllactosamine unit of poly-N-acetyllactosamine chains. Multiple alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Dec 2010]