

Product datasheet for **RG233983**

Clorf141 (NM_001276351) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Clorf141
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PSI00010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG233983 representing NM_001276351
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCAGAAAAATCCTAGAGAAGTTGGATGTCCTTGATAAGCAAGCAGAGATAATCTTGGCCAGAAGAA
CAAAGATAAACAGGCTTCAGAGTGAAGGAAGAAAAACAATATGGCTATACCCCTGACATTTGATTTTCA
GTTGGAATTTGAAGAAGCTCTTGCTACATCCGCTCTAAGGCAATATCAAAGATCAAAGAAGACAAAGTCA
TGCAGCATTACAAAATCAAAAATGCATGTCTCTTTCAAATGTGAGCCTGAACCTAGAAAGAGTAATTTTG
AAAAGTCAAATTTAAGACCATCTTTATTCAAACAAATGTAAAAATAAGAAAGTGAGTCAACAGCTCA
AATTGAAAAAACCTAGGAAACCATTTGGATTCTGTTGGTCTCTTAGAAGGTGATAGAAATAAAGAAAA
AAATCTCCACAGATGAACGATTTTAAATATAAAGAAAAACAATCGGTGAGAAATTATCAATTAAGTAAGT
ATAGGTGAGTAAGAAAGAAAGCTTGCTCCGTTGTGCTTTGAGGATGAATTGAAAAATCCACATGCCAA
GATAGTCAACGTTAGTCCAACAAAGACAGTAACCTCTCACATGGAACAAAGGACACAAATCCCATAATT
TTCCATGACACAGAATATGTACGAATGTTACTTTGACAAAAATAGATTTTCTCTCATCTTTGGAAA
ATGAAACATTTACCCACATAAAGAACAAATTTCAATTTAGAAAGAAATGTGAAATCCTCAATCTAT
AATTGGCAATCAATCTATTTCTCTTTCAAACCCCAAAAACTATGCCTACAGTACAGAGAAAGATATA
CAGATCCCTATGTCTTTAAAGCGGGCCACAACTGTAGATGATAAACTAAAGAAGAAACTAATAAGC
AGAGACTAGAAAAAGATCTTGAATACACTCTATAATTTCTACAGAATTTTCTAGCCTAACAAAACA
ATTTGTGGGTTACCTTGATAAAGCTGTTATTCATGAAATGAGTGCCCAAACTGGAAAAATTTGAAAGAATG
TTTTCTGCAGGAAACCAACGAGCATACCCACATCCAGTGCCTTACCTGTCAAATGTTACTCAAAGCCTT
TTAAATATATATGAATAAATAATGTAACGCCACTGGATAATTTGTTAACTTATCAAATGAAATTTT
AAATGCCTCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >RG233983 representing NM_001276351
 Red=Cloning site Green=Tags(s)

MAEKILEKLDVLDKQAEIILARRTKINRLQSEGRKTTMAIPLTFDFQLEFEEALATSASKAISKIKEDKS
 CSITKSKMHVSFKCEPEPRKSNFEKSNLRPFFIQTNVKNKESESTAQIEKKPRKPLDSVGLLEGDRNKRK
 KSPQMNDNFNIKENKSVRNYQLSKYRSVRKKSLLPLCFEDELKNPHAKIVNVSPKTKTVTSHMEQKDTNP
 FHDTEYVRMLLLTKNRFSSHPLENENIYPHKRTNFILERNCEILKSIIGNQSISLFPQKTMPTVQRKDI
 QIPMSFKAGHTTVDDKLKKKTNKQTLNRSWNTLYNFSQNFSSLTKQFVGYLDKAVIHEMSAQTGKFERM
 FSAGKPTS IPTSSALPVKCYSKPFKYIYELNNVTPLDNLLNLSNEILNAS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001276351

ORF Size: 1200 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001276351.2](#)

RefSeq Size: 2177 bp


RefSeq ORF: 1203 bp

Locus ID: 400757

UniProt ID: [Q5J VX7](#)

Cytogenetics: 1p31.3

Product images:



Circular map for RG233983