

Product datasheet for **RG233991**

ZNF322A (ZNF322) (NM_001242798) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZNF322A
Synonyms:	HCG12; ZNF322A; ZNF388; ZNF489
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG233991 representing NM_001242798
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTACACTTCAGAAGAGAAATGTAATCAGAGAACTCAAAAAGGAAAATATATAATGTATGCCCTCGGA
AGGGTAAAAAGATTTTATTCATATGCATGAGATTATTCAGATAGATGGTCATATATACCACTGCCTTGA
ATGCAAGCAAACTTCTGTGAAAACCTTAGCTCTTATTATGTGTGAGAGAACCATACTGGGGAGAAACCT
TATAATGTGATATGTGTGAGAAAACCTTTGTCCAAAGCTCAGATCTTACTTCACACCAGAGGATCCACA
ATTACGAGAAAACCTTATAAATGTAGCAAAATGTGAGAAGAGCTTTTGGCATCACTTAGCGCTTTCAGGACA
TCAGAGAACACATGCAGGTAAAAATTCTATACATGTGACATTTGTGGCAAGAATTTTGGTCAGAGTTCT
GATCTGCTTGTCCACCAGCGAAGCCATACTGGCGAGAAACCATATCTATGTAGTGAGTGTGACAAATGCT
TCAGTAGAAGTACAAACCTCATAAGGCATCGAAGAACTCACACAGGTGAGAAACCATTTAAGTGTCTCGA
GTGTGAAAAAGCTTTTAGTGGGAAATCAGATCTTATTAGCCACCAGAGAACTCACACTGGGGAAAGGCC
TACAAATGTAATAAGTGTGAGAAAAGTTACCGACACCGTTTCAGCCTTCATTGTACATAAAAGAGTTCATA
CTGGGGAGAAAGCCCTATAAGTGTGGTGCCTGTGAAAAATGCTTTGGCCAGAAATCAGACCTTATCGTGCA
CCAGAGAGTCCACACAGGTGAGAAGCCGTATAAATGCCTGGAATGTATGAGAAGTTTACTCGGAGTGCC
AACCTAATTAGGCACCGCAACTCACACTCACACTTTTAAATGCCTTGAATATGAAAAAAGCTTTAACT
GTAGCTCAGATCTTATTGTACATCAGAGAATTCACATGGAAGAGAAACCACATCAGTGGTCTGCGTGTGA
GAGTGGCTTCTCCTAGGAATGGACTTTGTTGCCAACAGAAAATGAGAACTCAAACAGAGGAGCTACAC
TATAAATCACTGTATGTGATAAAAGCTTCCACCAGAGTTCAGCCCTTCTCAACATCAGACAGTACACA
TTGGTGAAAAACCGTTTGTCTGTAATGTGAGTGAAAAAGGTCTTGAGCTTAGCCCTCCCATGCGTCAGA
AGCCTCACAGATGTCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG233991 representing NM_001242798
 Red=Cloning site Green=Tags(s)

MYTSEEKCNQRTQKRKIYNVCPRKGGKIFIMHEIIQIDGHIYQCLECKQNFCEALIMCERHTHTGEKP
 YKCDMCEKTFVQSSDLTSHQRIHNYEKPYKCSKCEKSFVHHLALSGHQRTHAGKKFYTCDICGKNFGQSS
 DLLVHQRSHTEKPYLCSECDKCFSRSTNLIRHRTHTGEKPKCLECEKAFSGKSDLIHQRTHTGERP
 YKCNKCEKSYRHRSAFIVHKRVHTGEKPYKCGACEKCFGQKSDLIHQRVHTGEKPYKCLECMRSFTRSA
 NLIRHQATHHTFKCLEYEKSFNCSSDLIVHQRIHMEEKPHQWSACESGFLLGMDFVAQQKMRTQTTEELH
 YKYTECDKSFHQSSALLQHQTIVHIGEKPFVCNVSEKGLELSPPHASEASQMS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001242798

ORF Size: 1206 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001242798.1](#), [NP_001229727.1](#)

RefSeq Size: 4747 bp

RefSeq ORF: 1209 bp


Locus ID: 79692

UniProt ID: [Q6U7Q0](#)

Cytogenetics: 6p22.2

Gene Summary: ZNF322A is a member of the zinc-finger transcription factor family and may regulate transcriptional activation in MAPK (see MAPK1; MIM 176948) signaling pathways (Li et al., 2004 [PubMed 15555580]).[supplied by OMIM, Mar 2008]

Product images:



Circular map for RG233991