

Product datasheet for **RG234010**

Aminoacylase 1 (ACY1) (NM_001198895) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Aminoacylase 1 (ACY1) (NM_001198895) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ACY1
Synonyms:	ACY-1; ACY1D; HEL-S-5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG234010 representing NM_001198895 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACCAGCAAGGGTCCCGAGGAGGAGCACCCATCGGTGACGCTCTTCCGCCAGTACCTGCGTATCCGCA
CTGTCCAGCCCAAGCCTGACTATGGAGCTGCTGTGGCTTTCTTTGAGGAGACAGCCCGCCAGCTGGGCT
GGGCTGTCAGAAAGTAGAGGTGGCACCTGGCTATGTGGTGACCGTGTGACCTGGCCAGGCCAACCCCT
ACACTCTCCTCCATCTTGCTCACTCCCACACGGATGTGGTGCCTGTCTCAAGGAACATTGGAGTCAGC
ACCCCTTTGAGGCCTTCAAGGATTCTGAGGGCTACATCTATGCCAGGGGTGCCAGGACATGAAGTGCCT
CAGCATCCAGTACCTGGAAGCTGTGAGGAGGCTGAAGGTGGAGGGCCACCGGTTCCCCAGAACCATCCAC
ATGACCTTTGTGCCTGATGAGGAGGTTGGGGTACCAAGGCATGGAGCTGTTCTGTCAGCGGCCTGAGT
TCCACGCCCTGAGGGCAGGCTTTGCCCTGGATGAGGGCATAGCCAATCCCACTGATGCCTTCACTGTCTT
TTATAGTGAGCGGAGTCCCTGGTGGGTGCGGGTACCAGCACTGGGAGGCCAGGCCATGCCTCACGCTTC
ATGGAGGACACAGCAGCAGAGAAGCTGCACAAGTTGTAACCTCCATCCTGGCATTCCGGGAGAAGGAAT
GGCAGAGGCTGCAGTCAAACCCCACTGAAAGAGGGGTCGCTGACCTCCGTGAACCTGACTAAGCTAGA
GGGTGGCGTGGCTATAACGTGATACCTGCCACCATGAGCGCCAGCTTTGACTTCCGTGTGCACCCGAT
GTGGACTTCAAGGCTTTTGGAGGACAGCTGCAGAGCTGGTCCAGGACAGTGGCGAGGGGTCACCCTAG
AGTTTGCTCAGAAGTGGATGCACCCCAAGTGACACCTACTGATGACTCAAACCTTGGTGGGCGAGCTTT
TAGCCGGGTCTGCAAGGATATGAACCTCACTCTGGAGCCTGAGATCATGCCTGCTGCCACTGACAACCGC
TATATCCGCGCGGTGGGGTCCCAGCTCTAGGCTTCTACCCATGAACCGCACACCTGTGCTGCTGCACG
ACCACGATGAACGGCTGCATGAGGCTGTGTTCTCCGTGGGGTGGACATATATACACGCCTGCTGCCTGC
CCTTGCCAGTGTGCTGCCCTGCCAGTGACAGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online >](#)

Protein Sequence: >RG234010 representing NM_001198895
 Red=Cloning site Green=Tags(s)

MTSKGPEEEHPSVTLFRQYLRIQVQPKPDYGAAVAFFEETARQLGLGCQKVEVAPGYVVTVLTPWGTNP
 TLSSILLNSHTDVVPVFKHWSHDPFEAFKDSEGYIYARGAQMCKVSIQYLEAVRRLKVEGHRFPRTIH
 MTFVDPDEEVGGHQGMELFVQRPEFHALRAGFALDEGIANPTDAFTVFYSERSPPWVVRTSTGRPGHASRF
 MEDTAAEKLHKVNSILAFREKEWQRLQSNPHLKEGVSVNLTKLEGGVAYNVIPATMSASFDFRVAPD
 VDFKAFEEQLQSWCQAAGEGVTLLEFAQKWMHPQVTPDSDSNPWWAAF SRVCKDMNLTLEPEIMPAATDNR
 YIRAVGV PALGFSPMNRTPVLLHDHDERLHEAVFLRGVDIYTRLLPALASVPALPSDS

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_001198895

ORF Size:	1224 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001198895.2
RefSeq Size:	1673 bp
RefSeq ORF:	1227 bp
Locus ID:	95
UniProt ID:	Q03154
Cytogenetics:	3p21.2
Protein Families:	Protease
Protein Pathways:	Arginine and proline metabolism, Metabolic pathways
Gene Summary:	This gene encodes a cytosolic, homodimeric, zinc-binding enzyme that catalyzes the hydrolysis of acylated L-amino acids to L-amino acids and an acyl group, and has been postulated to function in the catabolism and salvage of acylated amino acids. This gene is located on chromosome 3p21.1, a region reduced to homozygosity in small-cell lung cancer (SCLC), and its expression has been reported to be reduced or undetectable in SCLC cell lines and tumors. The amino acid sequence of human aminoacylase-1 is highly homologous to the porcine counterpart, and this enzyme is the first member of a new family of zinc-binding enzymes. Mutations in this gene cause aminoacylase-1 deficiency, a metabolic disorder characterized by central nervous system defects and increased urinary excretion of N-acetylated amino acids. Alternative splicing of this gene results in multiple transcript variants. Read-through transcription also exists between this gene and the upstream ABHD14A (abhydrolase domain containing 14A) gene, as represented in GeneID:100526760. A related pseudogene has been identified on chromosome 18. [provided by RefSeq, Nov 2010]