

Product datasheet for **RG234034**

KIAA0895 (NM_001199708) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIAA0895 (NM_001199708) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	KIAA0895
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG234034 representing NM_001199708 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAAAGCCCTGGTGCCTGTGACAAGCAGTGGTACTTGGTACTGCCTGGAGAGGCGGCCTGCTGTTTTTG
TCACTAGTTTCAGTGTCAAGTCCTGTAAAGTTTACACATGATATCTCTGTTACAGGGAATGGCATAGTACT
GCCACCTAAACCCAAAAGCAAGGTCAAGTGGTGCCATTTCTCCACTCTTCCAAAGCCAAAGCCTCAGCTG
TCTAGAAGCTTTGAAAAGGGAGATGACTTTTCTGGGAAGAAATTTGTATATTGACTGCTATAAAACCCA
CCAACCTAGAGAAAGAAAACTGAGATTCTTCAAATCTGACTATACCTACAATCCTCAGTTTGAGTATGC
CAATCCTGCTCTGCCAAGCGTATTAGCTAAGCATAGCCACGCATCTGACCGATTTCTTAAGCAGTCTATC
AATATTATGGAAGTACTTTACAGAAATATGGAAGTTATGAAAAATTTGAACAAGCCACTGGTGGTAGCT
TGCTCTCTAAACTCGAATCTGGAGTCACGTTAGGAAGTACATGATGAAGGAAGGCTGCCTAGGGGAGAT
TGTAAGTTCATCTCACTGAGGACCTGCTTTCCCGAGCGTCAATGACAGTAGTAAATGGATGTCGACTCTG
ACTATCAATGTGTCCACTGCACGTGAGCATTGGCTGGAGGGAATGCTGAGGCATGAAATAGGTACACATT
ATTTTCGAGGTATTAACAACCTCCAGCAGCCATGGAACAGTTGGACTGGACGTAAAAACATGAGCTAAA
GCCAAATAATCCACAGAGGAAGGACTAGCAAGCATTACAGTGTCTGTTTAGAAAAGACCCCTTTTTTA
TGGAGGGCTGCCCTCCTCTACTACACTGTTTATCAAGCCAGCCAAATGTCTTTTGTGAAGTCTTTAAAG
ATATTGGCAGTTTGTCAAGGACCCCAATACAAGATGGGATTATTGTGTACGGGCCAAGAGGGGATGGAC
TGATACTTCCCAACCAGGTGTTTTAGTAAAGACCAGGTATACTTGGATGGAATCTTCAAATCCTCCGA
TACAGAGATACCATAGACTTCCATCTGCTGACTGCTCTCGGGAAGGTTTCTTATGAAGATGTGGATCGCT
TAAAGGATTGGCAGTTACCGAAAACATGAGGGTCCCTCATTTCTGCAGGACCATGGCCGATATATGGA
ACACTTAGAGAAGATCATGGAAGTGAATGAACTGACTGACAGGGAAGTGAAGATCTTATA

ACGCGTACGCGGCGGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >RG234034 representing NM_001199708
 Red=Cloning site Green=Tags(s)

MKALVPVTSSGTWYCLERRPAVFVTSSVSPVKFTHDISVTGNGIVLPPKPKSKVKWCHFSTLPKPKPQL
 SRSFEKGDDFSGKKFCILTAIKPTNLEKEKLRFFKSDYTYPQFEYANPALPSVLAKHSHASDRFLKQSI
 NIMELTLQKYGSYEKFEQATGGSLLSKTRIWSHVRKYYMMKEGCLGEIVVHLTEDLLSRASMTVVNGCPTL
 TINVSTAREHWLEGLMRHEIGTHYFRGINNLQQPWNSWTGRKKHELKPNNPTEGLASIHSVLFRKDPFL
 WRAALLYTVYQASQMSFCELFKDIGRFVKDPNTRWDYCVRAKRGWTDTSQPGCFSKDQVYLDGILQILR
 YRDTIDFHLLTALGKVSIEDVDRLKGLAVTENMRVPFLQDHGRYMEHLEKIMEVNELTDRELKDLI

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001199708

ORF Size: 1251 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001199708.2

RefSeq Size: 4304 bp


RefSeq ORF: 1254 bp

Locus ID: 23366

UniProt ID: Q8NCT3

Cytogenetics: 7p14.2

Product images:



Circular map for RG234034