

## Product datasheet for **RG234059**

### **CSNK1G3 (NM\_001270572) Human Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CSNK1G3 (NM_001270572) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CSNK1G3
Synonyms:	CKI-gamma 3; CSNK1G3L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG234059 representing NM_001270572 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAAAATAAAAAGAAAGACAAGGACAAATCAGATGATAGAATGGCACGACCTAGTGGTCGATCGGGAC  
ACAACACTCGAGGAAGCTGGTCTTCATCGTCTGGAGTTTTAATGGTTGGACCTAACTTTAGAGTTGGAAA  
AAAAATGGATGTGGCAATTTGGAGAATTACGATTAGGGAAAAATTTATACACAAATGAATATGTGGCA  
ATTAAGTTGGAGCCCATGAAATCAAGAGCACACAGCTACATTTGGAATACAGATTCTATAAGCAGTTAG  
GATCTGGAGATGGTATACCTCAAGTTTACTATTTTCGGCCCTTGTGGTAAATACAATGCTATGGTGTGGA  
ACTGCTGGGACCTAGTTTGAAGACTTGTGGTGTGACAGAACATTTTCTCTAAAACAGTTCTC  
ATGATAGCTATACTAACTGATTTCTCGCATGGAATATGTCCATTCAAGAAGCTTGATATACAGAGATGTAA  
AACCTGAGAACTTCTAATAGGACGACCAGGAAACAAAACCCAGCAAGTTATTCACATTATAGATTTTGG  
TTTGGCAAAGGAATATATTGATCCGGAGACAAAGAAACACATACCATACAGAGAACACAAGAGCCTTACA  
GGAACAGCTAGATATATGAGCATAAACACACATTTAGGAAAAGAACAAAGTAGAAGAGACGATTTAGAAG  
CTTTAGGTCATATGTTTCATGATTTTCTGAGAGGCAGTCTTCTTGGCAAGGCTTAAAGGCTGACACATT  
AAAGGAGAGGTATCAGAAAATTTGGAGATACAAAACGGGCTACACCAATAGAAGTGTATGTGAAAATTTT  
CCAGAAGAAATGGCAACATATCTTCGTTATGTAAGAAGGCTAGATTTTTTTGAAAACAGACTATGACT  
ACTTAAGAAAGCTTTTTACTGACTTGTGGTATGTAAGAAGGATATATGTTGATTATGAATATGACTGGAT  
TGGTAAACAGTTGCCTACTCCAGTGGTGCAGTTCAGCAAGATCTGCTCTGTATCAAACAGAGAAGCA  
CATCAACACAGAGATAAGATGCAACAATCCAAAACAGGTTGTAAGTTCTACAAATGGAGAGTTAAACA  
CAGATGACCCACCGCAGGACGTTCAAATGCACCCATCACAGCCCTACTGAAGTAGAAGTATGATGGATGA  
AACCAACTGCCAGAAAGTGTGAACATGTGGTGTCTGTTTTTTCAAACGAAGGAAAAGGAAAACCAT  
CAGCGCCACAAA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG234059 representing NM\_001270572  
Red=Cloning site Green=Tags(s)

MENKKKDKDKSDDRMARPSGRSGHNTRGTGSSSSGVLVMPNFRVGGKIGCGNFGELRLGKNLYTNEYVA  
 IKLEPMKSRAPQLHLEYRFYKQLGSGDGIPOVYFPGCGKYNAMVLELLGPSLEDLFDLCDRTFSLKTVL  
 MIAIQLISRMEYVHSKNLIYRDVKPENFLIGRPGNKTQQVIHIIDFGLAKEYIDPETKKHIPPYREHKS  
 LTGTARYMSINTHLGKEQSRDDLEALGHMFMYFLRGLPWQGLKADTLKERYQKIGDTRATPIEVLNENF  
 PEEMATYLRVYRRLDFFEKPDYDLRKLFDLDRKGYMFDYEYDWIGQLPTPVGAVQQDPALSSNREA  
 HQHRDKMQQSKNQVVSSTNGELNTDDPTAGRSNAPITAPTEVEVMDETNCQKVLNMWCCCFKRRKRKTI  
 QRHK

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001270572

**ORF Size:** 1272 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001270572.1](#), [NP\\_001257501.1](#)

**RefSeq Size:** 4182 bp

**RefSeq ORF:** 1275 bp

**Locus ID:** 1456

**UniProt ID:** [Q9Y6M4](#)


**Cytogenetics:** 5q23.2

**Protein Families:** Druggable Genome, Protein Kinase

**Protein Pathways:** Hedgehog signaling pathway

**Gene Summary:** This gene encodes a member of a family of serine/threonine protein kinases that phosphorylate caseins and other acidic proteins. A related protein in the African clawed frog participates in the transmission of Wnt/beta-catenin signaling. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012]

**Product images:**



Circular map for RG234059