

Product datasheet for **RG234197**

ALDH2 (NM_001204889) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ALDH2 (NM_001204889) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ALDH2
Synonyms:	ALDH-E2; ALDHI; ALDM
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG234197 representing NM_001204889
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGTTGCGCGCTGCCGCCGCTTCGGGCCCGCCCTGGGCCCGCCCTCTTGTGACGCCGCCACCCAGG
 CCGTGCCTGCCCAACAGCAGCCGAGGTCTTCTGCAACCAGATTTTCATAAAACAATGAATGGCACGA
 TGCCGTACAGCAGGAAAACATTCCCCACCGTCAATCCGTCCACTGGAGAGGTCATCTGTGACGTAGCTGAA
 GGGGACAAGGCCCTGGAGACCCTGGACAATGGCAAGCCCTATGTCATCTCTACCTGGTGGATTTGGACA
 TGGTCTCAAATGTCTCCGGTATTATGCCGGCTGGGCTGATAAGTACCACGGGAAAACCTCCCCATTGA
 CGGAGACTTCTCAGCTACACACGCCATGAACCTGTGGGGGTGTGCGGGCAGATCATTCCGTGGAATTTCC
 CCGTCTGATGCAAGCATGGAAGCTGGGCCAGCCTTGCAACTGGAACCTGGTTGTGATGAAGGTAG
 CTGAGCAGACACCCCTCACCGCCCTATGTGGCCAACCTGATCAAGGAGGCTGGCTTTCCCCCTGGTGT
 GGTCAACATTGTGCTGGATTTGGCCCCACGGCTGGGGCCGCATTGCCTCCCATGAGGATGTGGACAAA
 GTGGCATTACAGGCTCCACTGAGATTGGCCCGTAATCCAGGTTGCTGCTGGGAGCAGCAACCTCAAGA
 GAGTGACCTTGGAGCTGGGGGGGAAGAGCCCAACATCATCATGTGAGATGCCGATATGGATTGGGCCGT
 GGAACAGGCCACTTCGCCCTGTTCTTCAACCAGGGCCAGTGTGCTGTGCCGGCTCCCGGACCTTCGTG
 CAGGAGGACATCTATGATGAGTTTGTGGAGCGGAGCGTTGCCCGGCCAAGTCTCGGGTGGTCCGGAAACC
 CCTTTGATAGCAAGACCGAGCAGGGGCCGAGGTGGATGAAACTCAGTTAAGAAGATCCTCGGCTACAT
 CAACACGGGAAGCAAGAGGGGGCGAAGCTGCTGTGTGGTGGGGCATTGCTGCTGACCGTGGTTACTTC
 ATCCAGCCCACTGTGTTGGAGATGTGCAGGATGGCATGACCATCGCAAGGAGGAGATCCTCGGGCCAG
 TGATGCAGATCCTGAAGTTCAAGACCATAGAGGAGGTTGTTGGGAGAGCCAACAATCCAGTACGGGCT
 GGCCGCAGCTGTCTTACAAAGGATTTGGACAAGGCAATTACCTGTCCAGGCCCTCCAGGCGGGCACT
 GTGTGGTCAACTGCTATGATGTGTTTGGAGCCAGTCACCCCTTGGTGGCTACAAGATGTCCGGGAGTG
 GCCGGGAGTTGGGCGAGTACGGGCTGCAGGCATACACTGAAGTAAAACCTGCACAGTCAAAGTGCCTCA
 GAAGAACTCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG234197 representing NM_001204889
 Red=Cloning site Green=Tags(s)

MLRAAARFGPRLGRLLSAAATQAVPAPNQPEVFCNQIFINNEWHDAVSRKTFPTVNPSTGEVICQVAE
 GDKALETLDNKPYVISYLVLDMLKCLRYAGWADKYHGKTIPIIDGDFSYTRHEPVGVCGQIIPWNF
 PLLMQAWKLGALATGNVVMKVAEQTPLTALYVANLIKEAGFPVGVVNI VPGFGPTAGAAIASHEDVVK
 VAFGTSTEIGRVIQVAAGSSNLKRVTLLELGGKSPNIIMSDADMDWAVEQAHFALFFNQGCCAGSRTFV
 QEDIYDEFVRSVARAKSRVGNPFDSKTEQGPVDETQFKKILGYINTGKQEGAKLLCGGGIAADRGYF
 IQPTVFGDVQDGMTIAKEEIFGPVMQILKFKTIEEVVGRANNSTYGLAAAFTKDLDKANYLSQALQAGT
 VVWNCYDVFQAQSPFGGYKMSGSGRELGEYGLQAYTEVKTVTVKVPQKNS

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001204889

ORF Size: 1410 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204889.2](#)

RefSeq Size: 1935 bp


RefSeq ORF: 1413 bp

Locus ID: 217

UniProt ID: [P05091](#)

Cytogenetics:	12q24.12
Protein Families:	Druggable Genome
Protein Pathways:	Arginine and proline metabolism, Ascorbate and aldarate metabolism, beta-Alanine metabolism, Butanoate metabolism, Fatty acid metabolism, Glycerolipid metabolism, Glycolysis / Gluconeogenesis, Histidine metabolism, Limonene and pinene degradation, Lysine degradation, Metabolic pathways, Propanoate metabolism, Pyruvate metabolism, Tryptophan metabolism, Valine, leucine and isoleucine degradation
Gene Summary:	This protein belongs to the aldehyde dehydrogenase family of proteins. Aldehyde dehydrogenase is the second enzyme of the major oxidative pathway of alcohol metabolism. Two major liver isoforms of aldehyde dehydrogenase, cytosolic and mitochondrial, can be distinguished by their electrophoretic mobilities, kinetic properties, and subcellular localizations. Most Caucasians have two major isozymes, while approximately 50% of East Asians have the cytosolic isozyme but not the mitochondrial isozyme. A remarkably higher frequency of acute alcohol intoxication among East Asians than among Caucasians could be related to the absence of a catalytically active form of the mitochondrial isozyme. The increased exposure to acetaldehyde in individuals with the catalytically inactive form may also confer greater susceptibility to many types of cancer. This gene encodes a mitochondrial isoform, which has a low Km for acetaldehydes, and is localized in mitochondrial matrix. Alternative splicing results in multiple transcript variants encoding distinct isoforms.[provided by RefSeq, Nov 2016]

Product images:



Circular map for RG234197