

Product datasheet for **RG234252**

ZNF345 (NM_001242474) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZNF345
Synonyms:	HZF10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RG234252 representing NM_001242474
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGAAACCTTACAAAACACAGCATTGAGTGTTCAAGTTTCAGAGGTGATTGGGAATGTAAAAACCAGT
 TTGAGAGAAAACAGGGATCTCAGGAAGGACATTTCAAGTGAATGATATTTACTCCTGAAGACATGCCAC
 TTTCAGTATCCAGCATCAGAGAATTCATACTGATGAGAAACCTTGAATGTAAGGAATGTGGGAAGGAT
 TTTAGTTTTGTATCAGTCCTTGTTCGACATCAGCGAATTCATACTGGTGAGAAACCTTATGAATGCAAAG
 AATGTGGCAAGGCCTTTGGTAGTGGTGCAAACCTTGCTTACCATCAAAGAATTCATACTGGTGAGAAGCC
 TTTTGAATGTAAAGAATGTGGGAAGGCCTTTGGTAGTGGCTCAAACCTTACTCACCATCAGAGAATTCAT
 ACTGGTGAGAAACCTATGAGTGTAAGGAATGTGGGAAGCCTTTAGTTTTGGATCAGGCCTTATTCGAC
 ATCAGATCATTACAGTGGTGAGAAGCCTTATGAGTGTAAGGAATGTGGGAAGTCCTTTAGTTTTGAATC
 AGCCCTTATTCGGCATCACAGAATTCACACAGGTGAGAAACCTTATGAATGTATAGATTGTGGTAAAGCC
 TTTGGCAGTGGTTCAAACCTTACTCAACATCGGCGGATTTCATACTGGTGAGAAACCTTATGAATGCAAAG
 CATGTGGAATGGCCTTTAGCAGTGGTTCGGCTCTTACTCGGCATCAGAGAATTCATACCGGTGAGAAAC
 ATATATATGTAATGAATGTGGTAAGGCCTTTAGTTTTGGATCAGCCCTTACTCGACATCAAAGAATTCAT
 ACTGGTGAGAAACCTTATGATGTAAAGGAATGTGGGAAGGCCTTTAATAGTGGCTCAGATCTCACTCAGC
 ATCAGAGAATTCACACTGGTGAGAAACCTATGAGTGTAAGGAGTGTGAGAAAGCCTTTAGAAGTGGTTC
 AAAACTTATTCAGCATCAAAGAATGCATACTGGAGAGAAACCTTATGAATGTAAGGAATGTGGGAAGACC
 TTTAGTAGTGGTTCAGACCTTACTCAACATCACAGAATTCATACTGGTGAGAAACCTTATGAATGTAAGG
 AATGTGGGAAGGCCTTTGGTAGTGGCTCAAACCTTATCCAACACCAGCTAATCCATACTGGTGAAGACC
 CTATGAATGTAAAGAATGTGGAAAGTCCTTTAGTAGTGGTTCAGCTCTTAATCGGCACCAGAGAATACAC
 ACTGGTGAGAAACCTATGAATGTAAGGAGTGTGGGAAGGCTTTTATAGTGGCTCAAGCCTTACTCAGC
 ATCAGAGAATTCATACAGGTGAGAACTTTATGAATGTAAGAACTGTGGGAAGGCTTATGGGAGGATTTC
 AGAGTTTCAGCAACATAAGAAAAGTCATAATGGTAAGAACTCTGCGAATTGGAACCTATAAAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG234252 representing NM_001242474
 Red=Cloning site Green=Tags(s)

MENLTKHSIECSSFRGDWECKNQFERKQGSQEGHFSEMIPTEDMPTFSIQHQRIHTDEKLLCECKECKD
 FSFVSVLVRHQRIHTGEKPYECKECKGAFGSGANLAYHQRIHTGEKPFECKECKGAFGSGSNLTHHQRIH
 TGEKPYECKECKGAFSFGSLIRHQIHSGEKPYECKECKGSFSESALIRHHRIHTGEKPYECIDCGKA
 FGSGSNLTQHRIHTGEKPYECKACGMAFSSGSALTRHQRIHTGEKPYICNECKGAFSFGSALTRHQRIH
 TGEKPYVCKECKGAFNSGSDLTQHRIHTGEKPYECKEKEKAFRSGSKLIQHQMHTGEKPYECKECKGT
 FSSGSDLTQHRIHTGEKPYECKECKGAFGSGSKLIQHQLIHTGERPYECKECKGSFSSGSALNRHQRIH
 TGEKPYECKECKGAFYSGSSLTQHRIHTGEKLYECKNCGKAYGRDSEFQHKHSHNGKKLCELETIN

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001242474

ORF Size: 1464 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001242474.1](#), [NP_001229403.1](#)

RefSeq Size:	3110 bp
RefSeq ORF:	1467 bp
Locus ID:	25850
UniProt ID:	Q14585
Cytogenetics:	19q13.12
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG234252