

Product datasheet for **RG234261**

RGS6 (NM_001204424) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RGS6 (NM_001204424) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	RGS6
Synonyms:	GAP; HA117; S914
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG234261 representing NM_001204424
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCTCAAGGATCCGGGGATCAAAGAGCAGTGGGGTTGCTGACCCAGAGGAGAGTTCTCCAAACATGA
 TCGTTTACTGCAAAATTGAAGACATCATTACAAAGATGCAAGATGACAAGACAGGGGGTGTGCCATCAG
 AACAGTCAAGAGCTTTCTCTCCAAAATCCCCAGTGTCTGTCACAGGTAAGTACTGACATTGTGCAGTGGCTTATG
 AAGAACCCTTCCATTGAGGACCCAGTTGAAGCAATACACTTGGGGAGCCTTATCGCTGCCAGGGCTACA
 TCTTTCCAATCTCAGACCATGTTCTCACCATGAAGGATGATGGCACCTTTTATCGTTTCCAGGCTCCGTA
 CTTCTGGCCTTGAAGTCTGGGAACCTGAAAACACTGACTATGCCATCTATCTCTGTAAGAGGACAATG
 CAAAATAAAGCAAGGCTGGAAGTGGCAGATTATGAAGCAGAAAACCTTAGCAAGACTCCAGAGGGCCTTTG
 CGAGGAAGTGGGAATTCATCTTTATGCAAGCAGAAGCACAAGTAAAGATTGACCGGAAAAAGACAAGAC
 AGAAAAGAAAATTTGGATAGTCAAGAAGCAGCCTTTTGGGATGTCCACAGGCCTGTGCCAGGCTGTGTG
 AACACAACAGAAATGGATATCCGAAAATGTCGACGTTTGAAGAATCCACAAAAGGTTAAAAAGTCCGGTGT
 ATGGCGTGACTGAAGAGTCCCAGGCACAGAGCCCGGTGCATGTAAGTACTCAGCCAACCAATCAGGAAAAACA
 AAAAGAGGACATCCGAAAACAGATAACATTTTTGAACGCACAGATCGACAGACATTTGTTTAAAAATGTCC
 AAAGTGGCTGAAAGTTAATTGCCTACACGGAACAATAATGTGGAATATGACCCCTTGATAACACCAGCTG
 AGCCATCCAACCTTTGGATCAGCGATGACGTTGCTTTGTGGGACATAGAGATGAGCAAAGAGCCAGCCA
 ACAGCGAGTAAAAAGATGGGGCTTCTCTTTGATGAGATATTGAAGGACCAGTGGGGCGGGACCAAGTTT
 CTACGATTCCTGGAGTCCGAATTCAGTTTCAAAAACCTCAGGTTCTGGCTGGCTCCAAGATCTTAAGA
 AACAAACCCTACAGGATGTGCCAAGAGGGTAGAAGAAATCTGGCAAGAGTTTCTGGCTCCAGGGGCTCC
 AAGTGCATCAACTGGATTCTCACAGCTATGAGATAACCAAGTCAAAATGTCAAAGATGGAGGGAGATAT
 ACATTTGAAGACGCCAGGAGCAGTCTACAAGCTGATGAAGAGTACAGCTATGCCCGCTTCTCCGGT
 CAAATGCTTACCAGGATTTGCTGCTGGCCAAGAAGAAGCCAGAAAGTGAAGCAAGGTCGTAGAAGTCCCT
 AGAAAAGTCACTCGCAGTGTGGAAAGTCTGGCGGGCAAGCGCCTCACGGGCTGATGCAGTCTCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG234261 representing NM_001204424
 Red=Cloning site Green=Tags(s)

MAQSGDQRAVGADPEESSPNMIVYCKIEDIITKMDDKTGGVPIRTVKSFLSKIPSVVTGTDIVQWLM
 KNLSIEDPVEAHLGSLIAAQGYIFPISDHVLTMKDDGTFYRFQAPYFWSNCWEPENTDYAIYLCRMT
 QNKARLELADYEENLRLQRAFARKWEFIFMQAEQVKIDRKKDKTERKILDSQERAFWDVHRPVPGCV
 NTTEMDIRKCRRLKNPQKVKSVYGVTEESQAQSPVHVLSPQIRKTTKEDIRKQITFLNAQIDRHCLKMS
 KVAESLIAYTEQYVEYDPLITPAEPSNPWISDDVALWDIEMSKEPSQQRVCRWGFDFEILKQVGRDQF
 LRFLESEFSSENLRFWLAQDLKKQPLQDVAKRVEEIQWQFLAPGAPSAINLDSHSYEITSQNVKDGGRY
 TFEDAQEHYKLMKSDSYARFLRSNAYQDLLLAKKKPESEQRRTSLEKFTRSVSKSLAGKRLTGLMQSS

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001204424

ORF Size: 1470 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204424.2](#)

RefSeq Size: 6073 bp

RefSeq ORF: 1473 bp

Locus ID: 9628


UniProt ID: [P49758](#)

Cytogenetics: 14q24.2

Protein Families: Druggable Genome

Gene Summary: This gene encodes a member of the RGS (regulator of G protein signaling) family of proteins, which are defined by the presence of a RGS domain that confers the GTPase-activating activity of these proteins toward certain G alpha subunits. This protein also belongs to a subfamily of RGS proteins characterized by the presence of DEP and GGL domains, the latter a G beta 5-interacting domain. The RGS proteins negatively regulate G protein signaling, and may modulate neuronal, cardiovascular, lymphocytic activities, and cancer risk. Many alternatively spliced transcript variants encoding different isoforms with long or short N-terminal domains, complete or incomplete GGL domains, and distinct C-terminal domains, have been described for this gene, however, the full-length nature of some of these variants is not known.[provided by RefSeq, Mar 2011]

Product images:



Circular map for RG234261