

Product datasheet for **RG234273**

TMEM200A (NM_001258277) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	TMEM200A
Synonyms:	KIAA1913; TTMA; TTMC
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >RG234273 representing NM_001258277
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGATAGCAACTGGTGGAGTGATAACTGGCCTGGCCGCTTAAAAAGGCAAGACTCTGCCAGATCACAGC
 AGCATGTCAACCTCAGCCCGTCTCCTGCTACCCAAGAGAAGAAGCCCATCAGGCGCCGGCCCCGGGCAGA
 TGTTGTGGTTGTTCTGTTGCAAAATCCGGCTTTATCCCATCTGGTTTTTTCTATTTTAGGAGTGCTC
 ATCTCCATTATAGGAATTGCTATGGCCGTCTTGGATATTGGCCCCAAAAGAACATTTTATTGATGCTG
 AAACAACACTGTCAACAATGAAACTCAGGTCAATCGGAATGAAGGCGGTGTGGTGGTTCGCTTCTTTGA
 GCAGCATTTGCATTCTGATAAGATGAAAAATGCTTGGCCCATTCACCATGGGGATTGGCATTTCATTTTC
 ATTTGTGCTAATGCCATTCTTCATGAAAACCGTGACAAAGAGACCAAAATCATACATGAGGGATATCT
 ATTCCACAGTCATTGACATTACACGCTAAGAATCAAGGAGCAAAGGCAATGAACGGCATGTACACTGG
 TTTGATGGGAGAAACAGAAGTAAACAGAATGGGAGCTCCTGTGCCTCGAGATTGGCAGCAATACGATC
 GCCTCTTTCTCGGGTTTTTCGGAGCAGTTTTCGAATGGACAGCTCCGTGGAGGAGGATGAACCTATGTTAA
 ATGAAGGTAAGAGTTCTGGGCATCTTATGCCCCCTTTGCTCTCTGACAGCTCTGTGTCTGCTTTGGCCT
 CTATCCACCTCCTTCCAAGACAACGATGATAAGACCAGCGGCTCTAAGAAATGTGAAACCAAGTCAATT
 GTGTCATCGTCCATCAGTGCTTTTACATTGCCTGTGATCAACTTAATAACTGTGTTATTGATGAGCCCA
 GTATAGATAACATCACTGAAGATGCTGACAACCTCAAAAGTAGGTCAAGGAATTTGTCAATGGATTCCCT
 TGTGGTTCTTTGCCCAACACCAAGTGAATCCTTCCAGCCCGTCAGCACAGTGCTACCAAGGAATAATTCC
 ATTGGGGAGTCGTTGTGAGTCAGTACAAGTCATCTATGGCTCTCGGACCTGGGGCTGGACAGCTCTTGT
 CTCCTGGGGCTGCCAGAAGACAGTTTGGGTCCAATACATCCTTGCAATTTGCTCTCGTCACACTCAAAGTC
 CTTGGACTTAGACCGGGTCCCTCCACTCTAACTGTTCAAGCAGAACACGGAACATCCAAGTTGGCCT
 AGGTTGGATCGGAACAACAGCAAGGGATATATGAACTAGAGAACAAAGAACCCGATGGATAGGTTGC
 TTGTGCCCCAAGTTGCCATCAAAAGGACTTTACCAATAAGGAGAAGCTTCTTATGATTTCAGATCTCA
 CAATAATTTGAGTTTTGAACATGATGAGTTTTTGAGTAACAACCTAAAGAGGGGAACCTCTGAAACAAGG
 TTT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG234273 representing NM_001258277
 Red=Cloning site Green=Tags(s)

MIATGGVITGLAALKRQDSARSQQHVNLSPSPATQEKPIRRRPRADVVRGKIRLYSPSGFFLILGVL
 ISIIGIAMAVLGYWPQKEHFIDAETTLSTNETQVIRNEGGVVVRFFEQLHSDKMKMLGPFTMGIGIFIF
 ICANAILHENRDKETKIIHMRDIYSTVIDIHLRIKEQRQMNGMYTGLMGETEVKQNGSSCASRLAANTI
 ASFSGFRSSFRMDSSVEEDELMLNEGKSSGHLMPPLLSDSSVSFVGLYPPPSKTTDDKTSKSKKCTKSI
 VSSSISAFITLPVIKLNNCVIDEPSIDNITEDADNLKRSRNLSDSLVPLPNTSESFQPVSTVLPNNNS
 IGESLSSQYKSSMALGPGAGQLLSPGAARRQFGSNTSLHLLSSHKSLLDLDRGPSTLTVQAEQRKHPSWP
 RLDNRNSKGYMKLENKEDPMDRLLPVQVAIKKDFTNKEKLLMISRSHNNLSFEHDEFLSNNLKRGTSER
 F

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001258277

ORF Size: 1473 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001258277.2](#)

RefSeq Size: 2988 bp

RefSeq ORF: 1476 bp
Locus ID: 114801
UniProt ID: [Q86VY9](#)
Cytogenetics: 6q23.1
Protein Families: Transmembrane

Product images:



Circular map for RG234273