

Product datasheet for **RG234297**

KAT7 (NM_001199157) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KAT7 (NM_001199157) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	KAT7
Synonyms:	HBO1; HBOA; MYST2; ZC2HC7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG234297 representing NM_001199157
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCCCGAAGGAAGAGGAATGCAGGCAGTAGTTCAGATGGAACCGAAGATTCCGATTTTCTACAGATC
 TCGAGCACACAGACAGTTTCAGAAAGTGATGGCACATCCCGACGATCTGCTCGAGTACCCCGCTCCTCAGC
 CAGGCTAAGCCAGAGTTCTCAAGATTCCAGTCTGTTTCGAAATCTGCAGTCTTTTGGCACTGAGGAGCCT
 GCTTACTCTACCAGAAGAGTGACCCGTAGTCAGCAGCAGCCTACCCAGTGACACCGAAAAATACCCTC
 TTCGGCAGACTCGTTCATCTGGTTCAGAAACTGAGCAAGTGGTTGATTTTTCAGATAGAGGACACCTTAC
 AGGAAAACATGAGAGACATTTCTCCATCTCAGGATGCCCACTGTATCATAACCTCTCAGCTGACGAATGC
 AAGGCACCAACGGAGAGGCAGCTTCGATATAAGGAAAAAGTGGCTGAACCTCAGGAAGAAAAGAAATCTG
 GACTGAGCAAAGAACAGAAAGAGAAATATATGGAACACAGACAGACCTATGGGAACACACGGGAACCTCT
 TTTAGAAAACCTGACAAGCGAGTATGACTTGGATCTTTTCCGAAGAGCACAAAGCCGGGCTTCAGAGGAT
 TTGGAGAAGTTAAGGCTGCAAGGCCAAATCACAGAGGGAAGCAACATGATTAACAATTGCTTTTGGCC
 GCTATGAGCTTGATACCTGGTATCATTCTCCATATCCTGAAGAATATGCACGGCTGGGACGTCTCTATAT
 GTGTGAATTCTGTTTAAATATATGAAGAGCCAAACGATACTCCGCCGGCACATGGCCAAATGTGTGTGG
 AAACACCCACCTGGTGATGAGATATATCGCAAAGTTCATCTCTGTGTTGAAGTGGATGGCAAGAAAA
 ACAAGTCTACTGCCAAAACCTGTGCCCTGTGGCCAACTTTTCTGGACCACAAGACATTATATATGA
 TGTGGAGCCCTTCTGTCTATGTTATGACAGAGGCGGACAACACTGGCTGTACCTGATTGGATATTTT
 TCTAAGGAAAAGAATTCATTCTCAACTACAACGTCCTGTATCCTTACTATGCCTCAGTACATGAGAC
 AGGGCTATGGCAAGATGCTTATTGATTTTCAGTTATTTGCTTTCCAAAGTGAAGAAAAGTTGGTCCCC
 AGAACGTCCTACTCTCAGATCTGGGGCTTATAAGCTATCGCAGTTACTGGAAAGAAGTACTTCTCCGCTAC
 CTGCATAATTTTCAAGGCAAAGAGATTTCTATCAAAGAAATCAGTCAGGAGACGGCTGTGAATCCTGTGG
 ACATTGTCAGCACTCTGCAAGCCCTCAGATGCTCAAATACTGGAAGGGAAAAACCTAGTTTTAAAGAG
 ACAGGACCTGATTGATGAGTGGATAGCCAAAGAGGCCAAAGGTCCAACCTCAATAAAACCATGGATCCC
 AGCTGCTTAAAATGGACCCCTCCAAGGGCACT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG234297 representing NM_001199157
 Red=Cloning site Green=Tags(s)

MPRRKRNAGSSSDGTEDSDFSTDLEHTDSESDGTSRRSARVTRSSARLSQSSQDSSPVRNLQSFQTEEP
 AYSTRRVTRSQQPTPVPKYPKYLRLQTRSSGSETEQVVDVSDRGHLTGKHERHFSISGCPLYHNLSADEC
 KAPTERQLRYKEKVAELRKKRNSGLSKEQKEKYMHRQTYGNTREPLENLTSEYDLDFRRAQARASED
 LEKLRLLQGQITEGSMIKTIAFGRYELDTWYHSPYPEEYARLGRLYMCEFLKYMKSQTILRRHMAKCVW
 KHPPGDEIYRKSISVFEVDGKKNKIYQNLCLLAKLFLDHKTLYYDVEPFLFYVMTEADNTGCHLIGYF
 SKEKNSFLNYNVSCILTMPQYMRQYGKMLIDFSYLLSKVEEKVGSPPERPLSDLGLISYRSYWKVLLRY
 LHNFGQKEISIKEISQETAVNPVDIVSTLQALQMLKYWKGKHLVLRQDLIDEWIAKEAKRSNSNKTMDP
 SCLKWTPPKGT

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001199157

ORF Size: 1503 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001199157.1](#), [NP_001186086.1](#)

RefSeq Size: 3334 bp

RefSeq ORF: 1506 bp

Locus ID: 11143


UniProt ID: [O95251](#)

Cytogenetics: 17q21.33

Protein Families: Druggable Genome, Stem cell - Pluripotency, Transcription Factors

Gene Summary: The protein encoded by this gene is part of the multimeric HBO1 complex, which possesses histone H4-specific acetyltransferase activity. This activity is required for functional replication origins and is involved in transcriptional activation of some genes. In both cases, the acetylation of histone H4 helps unfold chromatin so that the DNA can be accessed and replicated or transcribed. [provided by RefSeq, Oct 2016]

Product images:



Circular map for RG234297