

Product datasheet for **RG234336**

PDE1A (NM_001258313) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PDE1A (NM_001258313) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PDE1A
Synonyms:	CAM-PDE-1A; CAM-PDE 1A; HCAM-1; HCAM1; HSPDE1A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG234336 representing NM_001258313
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGATGACCATGTCACAATCAGGAAGAAACATCTCCAAAGACCCATCTTTAGACTAAGATGCTTGGTGA
 AGCAGCTGAAAAGAGGTGATGTTAACGTCGTCGACTTAAAGAAGAATATTGAATATGCGGCATCTGTGCT
 GGAAGCAGTTTATATCGATGAAACAAGAAGACTTCTGGATACTGAAGATGAGCTCAGTGACATTCAGACT
 GACTCAGTCCCATCTGAAGTCCGGGACTGGTTGGCTTCTACCTTTACACGGAAAAATGGGGATGACAAAA
 AGAAACCTGAGGAAAAACAAAAATTCGGAGCATTGTGCATGCTGTTCAAGCTGGAATTTTTGTGGAAG
 AATGTACCGAAAAACATATCATATGGTTGGTTGGCATATCCAGCAGCTGTCATCGTAACATTAAGGAT
 GTTGATAAATGGTCTTTTCGATGATTTGCCCTAAATGAAGCAAGTGGAGAGCATAGTCTGAAGTTTATGA
 TTTATGAACTGTTTACCAGATATGATCTTATCAACCGTTTCAAGATTCCTGTTTCTTGCCTAATCACCTT
 TGCAGAAGCTTTAGAAGTTGGTTACAGCAAGTACAAAAATCCATATCACAATTTGATTCATGCAGCTGAT
 GCACTCAAACCTGTGCATTACATAATGCTTCATACAGGTATCATGCACTGGCTCACTGAACCTGGAATTT
 TAGCAATGGTCTTTGCTGCTGCCATTATGATTATGAGCATACAGGGACAACAAACAACTTTACATTCA
 GACAAGGTGAGATGTTGCCATTTTGATAATGATCGCTCTGCTTGGAGAATCACACCGTGAAGTGCAGCT
 TATCGACTTATGCAAGAAGAAGAAATGAATATCTTGATAAATTTATCCAAAGATGACTGGAGGGATCTTC
 GGAACCTAGTGATTGAAATGGTTTTATCTACAGACATGTCAGGTCCTCCAGCAAAATAAAAATATAAG
 AAACAGTTTGCAGCAGCCTGAAGGGATTGACAGAGCCAAAACCATGTCCCTGATTCTCCACGCAGCAGAC
 ATCAGCCACCCAGCCAAATCCTGGAAGCTGCATTATCGGTGGACCATGGCCCTAATGGAGGAGTTTTTCC
 TGCAGGGAGATAAAGAAGCTGAATTAGGGCTTCCATTTTCCCACTTTGTGATCGGAAGTCAACCATGGT
 GGCCAGTCACAAATAGGTTTCATCGATTTCATAGTAGAGCCAACATTTTCTCTTCTGACAGACTCAACA
 GAGAAAATTGTTATTCCTTATAGAGGAAGCCTCAAAAGCCGAACTTCTTCTATGTGCAAGCAGCT
 CAACCACCATGTGGGTTACACATTGCTGATGCACTAAGACGATCAAATACAAAAGGCTCCATGAGTGA
 TGGGTCTATTCCCAGACTACTCCCTTGACAGCAGTGGACCTGAAGAGTTTCAAGAACAACCTGGTGGAC
 ATCATTACAGCAGAACAAGAGAGGTGAAAGAGTTAGCTGCACAAGAAGCAAGAACAGTTCACAGAAGT
 GTGAGTTTATTCATCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG234336 representing NM_001258313
 Red=Cloning site Green=Tags(s)

MDDHVTIRKKHLQRPIFRLRCLVKQLERGDVNVVDLKKNIEYAASVLEAVYIDETRRLLDTEDELSDIQT
 DSVPEVRDWLASTFTRKMGMTKKKPEEKPKFRSIVHAVQAGIFVERMYRKYTHMVGLAYPAAVIVTLKD
 VDKWDFVDFALNEASGEHSLKFMIEYELFTRYDLINRFKIPVSCLITFAEALVGVSKYKNPYHNLHAAD
 VTQTVHYIMLHTGIMHWLLEILAMVFAAAIHDYEHTGTTNNFHIQTRSDVAILYNDRSVLENHHVSAA
 YRLMQEEEMNILINLSKDDWRDLRNLVIEMVLSTDMSGHFQQIKNIRNSLQQPEGIDRAKTMSLILHAAD
 ISHPAKSWKLYHRWTMALMEEFFLQGDKEAELGLPFSPLCDRKSTMVAQSQIGFIDFIVEPTFSLTLDST
 EKIVIPLIEEASKAETSSYVASSSTTIVGLHIADALRRSNTKGSMSDGSYSPDYSLAAVDLKSFKNNLVD
 IIQONKERWKELAAQEARTSSQKCEFIHQ

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:



ACCN: NM_001258313

ORF Size: 1557 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001258313.1](#), [NP_001245242.1](#)

RefSeq Size: 2338 bp

RefSeq ORF: 1560 bp

Locus ID: 5136

UniProt ID: [P54750](#)


Cytogenetics: 2q32.1

Protein Families: Druggable Genome

Protein Pathways: Calcium signaling pathway, Progesterone-mediated oocyte maturation, Purine metabolism, Taste transduction

Gene Summary: Cyclic nucleotide phosphodiesterases (PDEs) play a role in signal transduction by regulating intracellular cyclic nucleotide concentrations through hydrolysis of cAMP and/or cGMP to their respective nucleoside 5-prime monophosphates. Members of the PDE1 family, such as PDE1A, are Ca(2+)/calmodulin (see CALM1; MIM 114180)-dependent PDEs (CaM-PDEs) that are activated by calmodulin in the presence of Ca(2+) (Michibata et al., 2001 [PubMed 11342109]; Fidock et al., 2002 [PubMed 11747989]).[supplied by OMIM, Oct 2009]

Product images:



Circular map for RG234336