

Product datasheet for **RG234472**

RUNX1T1 (NM_001198625) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RUNX1T1 (NM_001198625) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	RUNX1T1
Synonyms:	AML1-MTG8; AML1T1; CBFA2T1; CDR; ETO; MTG8; ZMYND2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG234472 representing NM_001198625
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGCCTGATCGTACTGAGAAGCACTCCACAATGCCAGACTCACCTGTGGATGTGAAGACGCAATCTAGGC
 TGA CTCTCCAACAATGCCACCTCCCCAACTACTCAAGGAGCTCCAAGAACCAGTTCATTTACACCGAC
 AACGTTAACTAATGGCACGAGCCATTCTCTACAGCCTTGAATGGCGCCCCCTCACCACCAATGGCTTC
 AGCAATGGGCTTCTCTTCTCTCTCTCTCTGGCTAATCAACAGCTGCCCCAGCCTGTGGTGCCA
 GGCAACTCAGCAAGCTGAAAAGTTCTTACTACCCTGCAGCAGTTTGGCAATGACATTTACCCGAGAT
 AGGAGAAAAGATTGCGACCCTCGTTCTGGGACTAGTGAACCTCACTTTGACAATTGAAGAATTTCAATTC
 AAAGTGAAGAAGCTACTAACTCCCACTGAGACCTTTTGTATCCCATTTTGAAGGCCAACTGCCCC
 TGCTGCAGCGTGAGCTCCTCCACTGCGCAAGACTGGCCAAACAGAACCCTGCCAGTACCTCGCCAGCA
 TGAACAGCTGCTTCTGGATGCCAGCACCTCACCTGTTGACTCCTCAGAGCTGCTTCTCGATGTGAAC
 GAAAACGGGAAGAGGCGAACTCCAGACAGAACCAAGAAAATGGCTTTGACAGAGAGCCTTTGCACTCAG
 AACATCCAAGCAAGCGACCATGCACTATTAGCCAGGCCAGCGGTACAGTCCAAATAACGGCTTATCCTA
 CCAGCCCAATGGCCTGCCTCACCTACCCACCTCCACCTCAGCATTACCGTTTGGATGATATGGCCATT
 GCCCACCCTACAGGGACTCCTATCGACACCCAGCCACAGGGACCTCAGGGACAGAAACAGACCTATGG
 GGTTGCATGGCACACGTCAAGAAGAAATGATTGATCACAGACTAACAGACAGAGAATGGGAGAAGAGTG
 GAAACATCTTGACCATCTGTTAACTGCATAATGGACATGGTAGAAAAACAAGGCGATCTCTACCCTGA
 CTAAGGCGGTGTCAAGAAGCAGACCGGGAAGAAATGAATTACTGGATCCGGCGGTACAGTGACGCCGAGG
 ACTTAAAAAAGGTGGCGCAGTAGCAGCAGCCACTCTAGGCAGCAGAGTCCCCTCAACCCAGACCCAGT
 TGCCTAGACCGCATCGGGAATTCCTTACAGGCCTGCGTCTGGATACGTGCCAGAGGAGATCTGGAAG
 AAAGCTGAGGAGGCCGTAATGAGGTGAAGCGCCAGGCGATGACGGAGCTGCAGAAGGCCGTGTCTGAGG
 CGGAGCGGAAAGCCACGACATGATCACACAGAGAGGGCCAAGATGGAGCGCACGGTCGCCGAGGCCAA
 ACGGCAGGCGGCGGAGGACGCACTGGCAGTTATCAATCAGCAGGAGGATTCAAGCGAGAGTTGCTGGAAT
 TGTGGCCGTAAGCGAGTGAAACCTGCAGTGGCTGTAAACAGCCCGATACTGTGGCTCATTTTGCACG
 ACAAGACTGGGAGAAGCACCATCACATCTGTGGACAGACCCTGCAGGCCAGCAGCAGGGAGACACACC
 TGCAGTCAGCTCCTGTGACGCCAACAGCGGGGCTGGGAGCCCGATGGACACACCACCAGCAGCCACT
 CCGAGGTCAACCACCCCGGAACCCCTCCACCATAGAGACAACCCCTCGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG234472 representing NM_001198625
 Red=Cloning site Green=Tags(s)

MPDRTEKHSTMPDSPVDVKTQSRLTPPTMPPPTTQGAPRTSSFPTTTLNNGTSHSPTALNGAPSPPNGF
 SNGPSSSSSSSLANQQLPPACGARQLSKLKRFLTTLQQFGNDISPEIGERVRLVGLVNSTLTIEEFHS
 KLQEATNFPLRPFVIFPLKANLPLLQRELLHCARLAKQNPAYLAQHEQLLLDASTTSPVDSSELLLDVN
 ENGRRTPDRTKENGFDREPLHSEHPSKRPTISPGQRYSPNGLSYQPNGLPHPTPPPPQHYRLDDMAI
 AHHRDYSYRHPSHRDLRDRNRMGLHGTRQEEMIDHRLTDREWAEWKHLHLLNCIMDMVEKTRRSLTV
 LRRQEAADREELNYWIRRYSDAEDLKKGGGSSSSSRQQSPVNPDPVALDAHREFLHRPASGYVPEEIKW
 KAEAEAVNEVKRQAMTELQKAVSEAERKAHDMITTERAKMERTVAEAKRQAAEDALAVINQQEDSSSESCWN
 CGRKASETCSGCNTARYCGSFCQHKDWEKHHHICGQTLQAQQQGDTPAVSSSVTPNSGAGSPMDTPPAAT
 PRSTTPGTPSTIETTPR

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001198625

ORF Size: 1731 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001198625.2](#)

RefSeq Size: 7769 bp

RefSeq ORF: 1734 bp

Locus ID: 862

UniProt ID: [Q06455](#)


Cytogenetics: 8q21.3

Protein Families: Transcription Factors

Protein Pathways: Acute myeloid leukemia, Pathways in cancer

Gene Summary: This gene encodes a member of the myeloid translocation gene family which interact with DNA-bound transcription factors and recruit a range of corepressors to facilitate transcriptional repression. The t(8;21)(q22;q22) translocation is one of the most frequent karyotypic abnormalities in acute myeloid leukemia. The translocation produces a chimeric gene made up of the 5'-region of the runt-related transcription factor 1 gene fused to the 3'-region of this gene. The chimeric protein is thought to associate with the nuclear corepressor/histone deacetylase complex to block hematopoietic differentiation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2010]

Product images:



Circular map for RG234472