

Product datasheet for **RG234490**

TNIP1 (NM_001252386) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TNIP1 (NM_001252386) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TNIP1
Synonyms:	ABIN-1; NAF1; nip40-1; VAN
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG234490 representing NM_001252386
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAAGCGACCAGGCTCCGGCAGAAGGCAGAGGAGCTAGTGAAGGACAACGAGCTGCTCCACCACCTT
 CTCCTCCTTGGGCTCCTTCGACCCCTGGCTGAGCTCACAGAAAGGACTCAAATGTCACAGCATCTCC
 CACAGCCCTGCATGCCCAAGTGAAGCCAGCACCAGTCCAGAAGCCTCCATCCAGTGGACCTCTCT
 GAATTTGAAGTGGTCACTCCTGAGGAGCAGAATTCACCAGAGAGCAGCAGCCATGCCAATGCGATGGCGC
 TGGGCCCTGCCCGTGGGACGGCACTGATGCTGCACCTGCAGCGCTGGAGACCACGCTGAGTGT
 GTGTGCCGAGGAGCCGACCACGGCCAGCTCTTACCCACCTGGGCCGATGGCCCTGGAGTTCAACCGA
 CTGGCATCCAAGGTGCACAAGAATGAGCAGCGCACCTCCATTCTGCAGACCCTGTGTGAGCAGCTTCGGA
 AGGAGAACGAGGCTCTGAAGCCAAGTTGGATAAGGGCTGGAACAGCGGGATCAGGCTGCCGAGAGGCT
 GCGGGAGGAAAAATTTGGAGCTCAAGAAGTTGTTGATGAGCAATGGCAACAAGAGGGTGCCTCTGGCGG
 CCAAGGCTACCGAAGATGGAAGGACAGGCAAGAAGGCAGTGGCTGGACAGCAGCAGGCTAGTGTGACGG
 CAGGTAAGGTCCCAGAGGTGGTGGCCTTGGGCGCAGCCGAGAAGAAGGTGAAGATGCTGGAGCAGCAGCG
 CAGTGAGCTGCTGGAAGTGAACAAGCAGTGGGACCAGCATTTCCGGTCCATGAAGCAGCAGTATGAGCAG
 AAGTCACTGAGCTGCGTCAGAAGCTGGCTGATTTGCAGAAGCAGGTGACTGACCTGGAGGCCGAGCGGG
 AGCAGAAGCAGCGTGACTTTGACCGCAAGCTCCTCGCCAAGTCCAAGATTGAAATGGAGGAGACCGA
 CAAGGAGCAGCTGACAGCAGAGGCCAAGGAGCTGCGCCAAAAGGTCAAGTACCTGCAGGATCAGCTGAGC
 CCACTACCCGACAGCGTGAGTACCAGGAAAAGGAGATCCAGCGGCTCAACAAGGCCCTGGAGGAAGCAC
 TGAGCATCAAACCCCGCCATCATCTCCACCAACAGCATTGGGAGCCCAAGAAGGACAGGGGCCCTCCT
 AAGGAAAACAGGAGCTGGTCACGCAAGAATGAGTTGCTGAAACAGCAGGTGAAGATCTTCGAGGAGGACTTC
 CAGAGGGAGCGCAGTGATCGTGAGCGCATGAATGAGGAGAAGGAAGAGCTGAAGAAGCAAGTGGAGAAGC
 TGCAGGCCAGGTCAACCTGTCAAATGCCAGCTAAAAGCATTCAAAGATGAGGAGAAGGCAAGAGAAGC
 CCTCAGACAGCAGAAGAGGAAAGCAAAGGCCTCAGGAGAGCGTTACCATGTGGAGCCCCACCCAGAACAT
 CTCTGCGGGCCCTACCCCTACGCTACCCGCCATGCCAGCCATGGTGCCACACCATGGCTTCGAGGACT
 GGTCCCAGATCCGCTACCCCTCCCCCATGGCCATGGAGCACCAGCCCACTCCCAACTCGCGCT
 CTTCCATCTGCCGAATACACCTGGCGTCTACCCTGTGGAGGGTTTCAAATCAAATCAGAGCTCCCAA
 GTGATGGACCCTCCACAGCCAGGCCTACAGAACCAGAGTCTCAAAAAATGACCGTGGGGCCCTCAG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG234490 representing NM_001252386
 Red=Cloning site Green=Tags(s)

MEATRLRQKAELVKDNELLPPSPSLGSFDPLAELTGKDSNVTASPTAPACPSDKPAPVQKPPSSGTSS
 EFEVVTPEEQNSPESSSHANAMALGPLPREDGNLMLHLQRLETTLSVCAEEDPHGQLFTHLGRMALEFNR
 LASKVHKNEQRTSILQTLCEQLRKENEALKAKLDKLEQRDQAAERLREENLELKKLLMSNGNKEGASGR
 PGSPKMEGTGKKAVAGQQQASVTAGKVPEVVALGAAEKVKMLEQQRSELLEVNKQWDQHFRSMKQQYEQ
 KITELRQKLADLQKQVTDLEAEREQKQRFDRKLLAKSKIEMEETDKEQLTAEAKELRQKVYLQDQLS
 PLTRQREYQEKIQRLNKALEEALSIQTPPSSPPTAFGSPEGAGALLRKQELVTQNELLKQVVKIFEEDF
 QRERSDRERMNEEKEELKKQVEKLQAQVTL SNAQLKAFKDEEKAREALRQQKRKAKASGERYHVEPHPEH
 LCGAYPYAYPPMPAMVPHHGFEDWSQIRYPPPPMAMEHPPPLPNSRLLFHLPEYTWRLPCGGVRNPNQSSQ
 VMDPPTARPTPEPSPKNDREGPQ

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001252386

ORF Size: 1749 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001252386.2](#)

RefSeq Size: 2761 bp

RefSeq ORF: 1752 bp


Locus ID: 10318

UniProt ID: [Q15025](#)

Cytogenetics: 5q33.1

Gene Summary:

This gene encodes an A20-binding protein which plays a role in autoimmunity and tissue homeostasis through the regulation of nuclear factor kappa-B activation. Mutations in this gene have been associated with psoriatic arthritis, rheumatoid arthritis, and systemic lupus erythematosus. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2011]

Product images:


Circular map for RG234490