

Product datasheet for **RG234561**

SLC5A11 (NM_001258413) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC5A11 (NM_001258413) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SLC5A11
Synonyms:	KST1; RKST1; SGLT6; SMIT2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG234561 representing NM_001258413
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGGTGTGGTGGCCAGTGGGTGCATCCTTGTGTTGCCAGCAATGTTGGAAGTGGACATTTTCATTGGCCTGG
 CAGGGTCAGGTGCTGCTACGGGCATTTCTGTATCAGCTTATGAACCTAATGGCTTGTCTGTGCTGAT
 GTTGGCCTGGATCTTCTACCATCTACATTGCTGGTCAGGTCAACACGATGCCAGAATACCTACGGAAG
 CGCTTCGGTGGCATCAGAATCCCATCATCCTGGCTGACTCTACCTATTTATCTACATCTTCACCAAGA
 TCTCGGTAGACATGTATGCAGGTGCCATCTTCACAGCAGTCTTGCACCTGGATCTGTACCTGGCCAT
 AGTTGGGCTACTGGCCATCACTGCTGTATACACGGTGTGGTGGCCTGGCTGCTGTGATCTACACGGAT
 GCCCTGCAGACGCTGATCATGCTTATAGGAGCGCTCACCTTGATGGGCTACAGTTTCGCCCGGTTGGTG
 GGATGGAAGGACTGAAGGAGAAGTACTTCTGGCCCTGGTAGCAACCGAGTGAGAACAGCAGCTGCGG
 GCTGCCCGGGAAGATGCCTTCCATATTTCCGAGATCCGCTGACATCTGATCTCCGTCGCCGGGGGTC
 CTATTTGGAATGTCCATCCCATCCCTCTGGTACTGGTGCACGGATCAGGTGATTGTCCAGCGGACTCTGG
 CTGCCAAGAACCTGTCCCATGCCAAAGGAGGTGCTCTGATGGCTGCATACCTGAAGGTGCTGCCCTCTT
 CATAATGGTGTTCCTGGGATGGTCAGCCGCATCCTCTTCCAGATCAAGTGGCCTGTGCAGATCCAGAG
 ATCTGCCAGAAGATCTGCAGCAACCCCTCAGGCTGTTCCGACATCGCGTATCCCAAACCTCGTGCTGGAAC
 TCCTGCCACAGGGCTCCGTGGGCTGATGATGGCTGTGATGGTGGCGGCTCTCATGTCTCCCTCACCTC
 CATCTTTAACAGTGCCAGCACCATCTTCACCATGGACCTCTGGAATCACCTCCGGCTCGGGCATCTGAG
 AAGGAGCTCATGATTGTGGGCAGGGTGTGTTGCTGCTGCTGGTCTCCATCCTCTGGATCCCTG
 TGGTCCAGGCCAGCCAGGGCGGCCAGCTCTTCATCTATATCCAGTCCATCAGCTCCTACCTGCAGCCGCC
 TGTGGCGGTGGTCTTCATCATGGGATGTTTCTGGAAGAGGACCAATGAAAAGGGTGCCTTCTGGGGCTG
 ATCTCGGGCTGCTCCTGGGCTTGGTTAGGCTGGTCTGGACTTTATTTACGTGCAGCCTCGATGCGACC
 AGCCAGATGAGCGCCCGGTCTGGTGAAGAGCATTCACTACCTCTACTTCTCCATGATCTGTCCACGGT
 CACCTCATCACTGTCTCCACCGTGAGCTGGTTCACAGAGCCACCCTCCAAGGAGATGGTCAGCCACCTG
 ACCTGGTTTACTCGTCACGACCCCGTGGTCCAGAAGGAACAAGCACCACCAGCAGCTCCCTTGCTCTTA
 CCCTCTCTCAGAACGGGATGCCAGAGGCCAGCAGCAGCAGCGTCCAGTTCGAGATGGTTCAAGAAAA
 CACGTCTAAAACCCACAGCTGTGACATGACCCCAAAGCAGTCAAAGTGGTGAAGGCCATCTGTGGCTC
 TGTGGAATACAGGAGAAGGGCAAGGAAGAGCTCCCGGCCAGAGCAGAAGCCATCATAGTTTCCCTGGAAG
 AAAACCCCTTGGTGAAGACCCTCTGGACGTCAACCTCATTTTCTGCGTGAGCTGCGCCATCTTTATCTG
 GGGCTATTTTGCT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG234561 representing NM_001258413
 Red=Cloning site Green=Tags(s)


MVWVPVGASLFASNVGSGHFIGLAGSGAATGISVSAYELNGLFSVLMWIFLPIYIAGQVTTMPEYLRK
 RFGGIRIPIILAVLYLFYIFTKISVDMYAGAFIQQLHLDLYLAIVGLLAITAVYTVAGGLAAVIYTD
 ALQTLIMLIGALTLMGYSFAAVGGMEGLKEYFLALASNRSENSSCGLPREDAFHIFRDPLTSDLWPWPGV
 LFGMSIPSLWYWCTDQVIVQRTLAANKLSHAKGGALMAAYLKVLPLFIMVFPGMVSRIIFPDQVACADPE
 ICQKICSNPSGCSDIAYPKLVLELLPTGLRGLMMVMAALMSSLTSIFNSASTIFTMDLWNHLRPRASE
 KELMIVGRVFLVLLVLSILWIPVVQASQGGQLFIYIQSISSYLQPPVAVVFIMGCFWKRTNEKGAFWGL
 ISGLLLGLVRLVLDFIYVQPRCDQDERPVLVKSIIHYLYFSMILSTVTLITVSTVSWFTEPPSKEMVSHL
 TWFTRHDPVVQKEQAPPAAPLSLTLSQNGMPEASSSSSVQFEMVQENTSKTHSCDMTPKQSKVVKAILWL
 CGIQEKGEELPARAEAIIVSLEENPLVKTLLDVNLIFCVSCAIFIWGYFA

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001258413

ORF Size: 1833 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001258413.3](#)

RefSeq Size: 2353 bp

RefSeq ORF: 1836 bp


Locus ID: 115584

UniProt ID: [Q8WWX8](#)

Cytogenetics: 16p12.1

Gene Summary: Cotransporters, such as SLC5A11, represent a major class of proteins that make use of ion gradients to drive active transport for the cellular accumulation of nutrients, neurotransmitters, osmolytes, and ions Roll et al. (2002) [PubMed 12039040]. [supplied by OMIM, Mar 2008]

Product images:



Circular map for RG234561