

Product datasheet for **RG234633**

SLC5A11 (NM_001258411) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC5A11 (NM_001258411) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SLC5A11
Synonyms:	KST1; RKST1; SGLT6; SMIT2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG234633 representing NM_001258411
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGAGCGGCACCAGCAGCCCTCAGCCTCCACAGTTAGATCCCCTGGATGCGTTTCCCAGAAGGGCT
 TGGAGCCTGGGGACATCGCGGTGCTAGTTCTGTACTTCTCTTTGTCTGGCTGTTGGACTATGGTCCAC
 AGTGAAGACAAAAGAGACACAGTGAAAGGCTACTTCTGGCTGGAGGGGACATGGTGTGGTGGCCAGTG
 GGTGCATCCTTGTGGCCAGCAATGTTGGAAGTGGACATTTTCATTGGCCTGGCAGGGTCAAGTGCTGCTA
 CGGGCATTCTGTATCAGCTTATGAACCTAATGGCTTGTCTGTGCTGATGTTGGCCTGGATCTTCTCT
 ACCCATCTACATTGCTGGTCAGGTAGACATGTATGCAGGTGCCATCTTCATCCAGCAGTCTTTCACCTG
 GATCTGTACCTGGCCATAGTTGGGCTACTGGCCATCACTGCTGTATACAGGTTGCTGGTGGCCTGGCTG
 CTGTGATCTACACGGATGCCCTGCAGACGCTGATCATGCTTATAGGAGCGCTCACCTTGATGGGCTACAG
 TTTCCGCCGCGTTGGTGGGATGGAAGGACTGAAGGAGAAGTACTTCTGGCCCTGGCTAGCAACCGGAGT
 GAGAACAGCAGCTGCGGGCTGCCCGGGAAGATGCCTTCCATATTTCCGAGATCCGCTGACATCTGATC
 TCCCGTGGCCGGGGTCTATTTGGAATGTCCATCCCCTCCCTCTGGTACTGGTGACCGGATCAGGTGAT
 TGTCCAGCGGACTCTGGCTGCCAAGAACCTGTCCCATGCCAAGGAGGTGCTCTGATGGCTGCATACCTG
 AAGGTGCTGCCCTCTTCATAATGGTGTCCCTGGGATGGTCAGCCGCATCTCTTCCCAGATCAAGTGG
 CCTGTGCAGATCCAGAGATCTGCCAGAAGATCTGCAGCAACCCCTCAGGCTGTTCCGGACATCGCGTATCC
 CAAACTCGTGCTGGAACCTGCCACAGGGCTCCGTGGGCTGATGATGGCTGTGATGGTGGCGGCTCTC
 ATGTCTCCCTCACCTCCATCTTAAACAGTGCCAGCACCATCTTCACCATGGACCTCTGGAATCACCTCC
 GGCTCGGGCATCTGAGAAGGAGCTCATGATTGTGGCAGGGTGTGCTGCTGCTGGTCTGTCTGCTC
 CATCTCTGGATCCCTGTGGTCCAGGCCAGGCCAGGGCGGCCAGCTCTTCATCTATATCCAGTCCATCAGC
 TCTACCTGCAGCCGCTGTGGCGGTGGTCTTCATCATGGGATGTTTCTGGAAGAGGACCAATGAAAGG
 GTGCCCTCTGGGCGCTGATCTCGGCGCTGCTCTGGGCTTGGTTAGGCTGGTCTGGACTTTATTTACGT
 GCAGCCTCGATGCGACCAGCCAGATGAGCGCCCGTCTGGTGAAGAGCATTCACTACCTCTACTTCTCC
 ATGATCTGTCCACGGTCACCTCATCACTGTCTCCACCGTGAGCTGGTTCACAGAGCCACCCTCCAAGG
 AGATGGTCAGCCACCTGACCTGGTTACTCGTCACGACCCCGTGGTCCAGAAGGAACAAGCACCACCAGC
 AGCTCCCTTGTCTTACCCTCTCTCAGAACGGGATGCCAGAGGCCAGCAGCAGCAGCAGCGTCCAGTTC
 GAGATGGTTCAAGAAAACAGCTCTAAACCCACAGCTGTGACATGACCCAAAGCAGTCCAAAGTGGTGA
 AGGCCATCCTGTGGCTCTGTGGAATACAGGAGAAGGGCAAGGAAGAGCTCCCGGCCAGAGCAGAAGCCAT
 CATAGTTTCCCTGGAAGAAAACCCCTTGGTGAAGACCCTCTGGACGTCAACCTCATTTTCTGCGTGAGC
 TCGCCATCTTTATCTGGGGCTATTTTGCT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:


>RG234633 representing NM_001258411
 Red=Cloning site Green=Tags(s)

MESGTSPPQPPQLDPLDAFPQKLEPGDIAVLVLYFLFVLAVGLWSTVTKRDTVKGYFLAGGDMVWWPV
 GASLFASNVGSGHFIFLAGSGAATGISVSAYELNGLFVLMWLAWIFLPIYIAGQVDMYAGAIIFIQSLHL
 DLYLAIVGLLAITAVYTVAGGLAAVIYTDALQTLIMLIGALTLMGYSFAAVGGMEGLKEKYFLALASNRS
 ENSSCGLPREDAFHIFRDPLTSDLWPWGVLFMSIPSLWYWCTDQVIVQRTLAAKNLSHAKGGALMAAYL
 KVLPLFIMVFGMVSRIIFPDQVACADPEICQKICSNPSGCSDIAYPKLVLELLPTGLRGLMMAMVMAAL
 MSSLTSIFNSASTIFTMDLWNHLRPRASEKELMIVGRVFLLLVLVSILWIPVVQASQGGQLFIYIQSIS
 SYLQPPVAVVFIMGCFWKRTNEKGAFWGLISGLLLGLVRLVLDIFIYVQPRCDQPDERPVLVKSIIHYLYFS
 MILSTVTLITVSTVSWFTEPPSKEMVSHLTWFTRHDPVVQKEQAPPAAPLSLTLSONGMPAEASSSSSVQF
 EMVQENTSKTHSCDMTPKQSKVVKAILWLCGIQEKKEELPARAEAIIVSLEENPLVKTLTLDVNLIFCVS
 CAIFIWGYFA

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001258411

ORF Size: 1920 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001258411.2](#), [NP_001245340.1](#)

RefSeq Size: 2407 bp

RefSeq ORF: 1923 bp


Locus ID: 115584

UniProt ID: [Q8WWX8](#)

Cytogenetics: 16p12.1

Gene Summary: Cotransporters, such as SLC5A11, represent a major class of proteins that make use of ion gradients to drive active transport for the cellular accumulation of nutrients, neurotransmitters, osmolytes, and ions Roll et al. (2002) [PubMed 12039040]. [supplied by OMIM, Mar 2008]

Product images:



Circular map for RG234633