

Product datasheet for **RG234685**

DIP13B (APPL2) (NM_001251904) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DIP13B (APPL2) (NM_001251904) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	DIP13B
Synonyms:	DIP13B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG234685 representing NM_001251904
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCCCGCCGTGGACAAGCTCCTGCTAGAGGAGCGTTGCAGGACAGCCCCAGACTCGCTCTTTACTGA
 GCGTGTGTTGAAGAAGATGCTGGCACCCTCACAGACTATACCAACCAGCTGCTCCAGGCAATGCAGCGCGT
 CTATGGAGCCCAGAATGAGATGTGCTGGCCACACAACAGCTTTCTAAGCAACTGCTGGCATATGAAAA
 CAGAACTTTGCTCTTGGCAAAGGTGATGAAGAAGTAATTTCAACACTCCACTATTTTTCCAAAGTGGTGG
 ATGAGCTTAATCTTCTCCATACAGAGCTGGCTAAACAGTTGGCAGACACAATGGTTCTACCTATCATACA
 ATTCGAGAAAAGGATCTCACAGAAGTAAGCACTTTAAAGGATCTATTTGGACTCGCTAGCAATGATGTC
 TGTTTTATTTTAGAGCATGACCTCTCAATGGCAAAATACAGCAGGCTGCCTAAGAAAAAGGAGAATGAGA
 AGGTGAAGACCGAAGTCGAAAAGAGGTGGCCGCGCCGCGGAAGCAGCACCTCTCTCCCTTCAGTA
 CTACTGTGCCCTCAACGCGCTGCAGTACAGAAAGCAAATGGCCATGATGGAGCCCATGATAGGCTTTGCC
 CATGGACAGATTAACCTTTTTTAAGAAGGGAGCAGAGATGTTTTCCAAACGTATGGACAGCTTTTTATCCT
 CCGTTGCAGACATGGTTCAAAGCATTCAAGTAGAACTGGAAGCCGAGGCGGAAAAAGATGCGGGTGTCCCA
 GCAAGAATTACTTTCTGTTGATGAATCTGTTTACACTCCAGACTCTGATGTGGCCGACCCACAGATCAAC
 AGGAACCTCATCCAGAAGGCTGGTTACCTTAATCTTAGAAACAAAACAGGGCTGGTCACCACCACCTGGG
 AGAGGCTTTATTTCTTACCCAAAGGCGGGAATCTCATGTGTGAGCCAGGGGAGCCGTGGCTGGAGGTTT
 GATCCAGGACCTGGACAACCTGCTCAGTGATGGCCGTGGATTGCGAAGACCGGCGCTACTGCTTCAGATC
 ACCACGCCCAATGAAAAATCGGGAATAATCCTCCAGGCTGAGAGCAGAAAGGAAAAATGAAGAGTGGATAT
 GTGCAATAAACAACTCTCCAGACAGATCTACCTGACCGACAACCCTGAGGCAGTCGCGATCAAGTTGAA
 TCAGACCGCTCTGCAAGCAGTGACTCCCATTAACAAGTTTTGGAAAAAACAAGAAAGCTCATGCCCCAGC
 CAGAACCTGAAAAATTCAGAGATGGAAAAATGAAAATGACAAGATTGTTCCCAAAGCAACAGCCAGTCTAC
 CTGAAGCAGAGGAGCTGATCGCGCCTGGAACGCCGATTCAATTGATATTGTGCTTCCTGCTACAGAATT
 CCTTGATCAGAACAGAGGGAGCAGGCGTACCAACCCTTTTGGTGAAGTGAAGATGAATCATTTCCAGAA
 GCAGAAGATTCTTTTGCAGCAGATGTTTATAGTTTCGGTTTTTGGGATCAATGGCAGTTAAACAGACA
 GCACTACTGAAGTGATTTATGAAGCGATGAGACAAGTATTGGCTGCTCGGGCTATTCATAACATCTCCG
 CATGACAGAATCCCATCTGATGGTCACCAGTCAATCTTTGAGGTTGATAGATCCACAGACTCAAGTATCA
 AGGGCCAATTTTGAACCTACCAGTGTACACAATTTGCTGCTCATCAAGAAAACAAGAGACTGGTTGGTT
 TTGTCATCCGTGTTCTGAATCCACTGGAGAAGATCTCTGAGTACATACATTTTTGAAAGCAACTCAGA
 AGGCGAAAAGATATGTTATGCTATTAATTTGGGAAAAGAAATTATTGAGGTTTCAAGGATCCAGAAGCA
 CTGGCTCAATTAATGCTGTCCATACCACTAACCAATGACGGAAAAATGTACTGTTAAACGATCAACCAG
 ATGACGATGATGGAAATCCAAACGAACATAGAGGCGCAGAATCCGAAGCA

ACGCGTACGCGGGCGCTCGAG – GFP Tag – GTTTAA


Protein Sequence: >RG234685 representing NM_001251904
 Red=Cloning site Green=Tags(s)

MPAVDKLLLEEALQDSPQTRSLLSVFEEDAGTLTDYTNQLLQAMQRYGAQNEMCLATQQLSKQLLAYEK
 QNFALGKGDEEVISTLHYFSKVDELNLLHTELAKQLADTMVLPPIQFREKDLTEVSTLKDLFGLASNDV
 CLFLEHDLMAKYSRLPKKKENEKVKTEVGKEVAAARRKQHLSSLQYYCALNALQYRKQAMMEPMIGFA
 HGQINFFKKGAEMFSKRMDSFLSSVADMVQSIQVELEAEAEKMRVSQQELLSDVESVYTPDSDVAAPQIN
 RNLIQKAGYLNLRNKTGLVTTTWERLYFFTQGGNLMCQPRGAVAGGLIQDLDNCSVMAVDCEDRRYCFQI
 TTPNGKSGIILQAESRKENEWICAINNISRQIYLTDNPEAVAIAKLNTALQAVTPITSFGKKQESSCPS
 QNLKNSEMNENDKIVPKATASLPEAEELIAPGTPIQFDIVLPATEFLDQNRGSRRTNPFGETEDESFP
 AEDSLLQQMFIIVRFLGSAVKTDSTTEVIYEAMRQVLAARAIHNIFRMTESHLMTVSQSLRLIDPQTQVS
 RANFELTSVTQFAAHQENKRLVGFVIRVPESTGEESLSTYIFESNSEGEKICYAINLGKEIEVQKDPEA
 LAQLMLSIPLTNDGKYVLLNDQDDDDGNPNEHRGAEESEA

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001251904

ORF Size: 2010 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001251904.1](#), [NP_001238833.1](#)

RefSeq Size: 3307 bp

RefSeq ORF: 2013 bp


Locus ID: 55198

UniProt ID: [Q8NEU8](#)

Cytogenetics: 12q23.3

Gene Summary: The protein encoded by this gene is one of two effectors of the small GTPase RAB5A/Rab5, which are involved in a signal transduction pathway. Both effectors contain an N-terminal Bin/Amphiphysin/Rvs (BAR) domain, a central pleckstrin homology (PH) domain, and a C-terminal phosphotyrosine binding (PTB) domain, and they bind the Rab5 through the BAR domain. They are associated with endosomal membranes and can be translocated to the nucleus in response to the EGF stimulus. They interact with the NuRD/MeCP1 complex (nucleosome remodeling and deacetylase /methyl-CpG-binding protein 1 complex) and are required for efficient cell proliferation. A chromosomal aberration t(12;22)(q24.1;q13.3) involving this gene and the PSAP2 gene results in 22q13.3 deletion syndrome, also known as Phelan-McDermid syndrome. [provided by RefSeq, Oct 2011]

Product images:



Circular map for RG234685